Sonic Rubber Series Manual for Multipurpose phantom N-365



Targets layout



This phantom is designed to allow scanning from all 4 side walls.

Axtial resolution



185/1861 12Hz

Depth of the shallowest targets: 15,30,50mm depth

Target diameter: 0.05mm

Spaces between targets:

0.5,1,2,3mm

Linear Probe Image



Example

Axial resolution:1mm

2 targets with 0.5mm clearance are not recognized separately.

Angular resolution



Linear Probe Image



Use the same targets as axial resolution

Depths:

15,30,50mm depth

Target diameter: 0.05mm

Spaces between target:

0.5,1,2,3mm

Example

Angular resolution:1mm

3 targets with 0.5mm, 1mm clearance are not recognized separately.

Close range resolution





Depth range of 1-10 mm depth (10 targets are embedded) Target diameter 0.1mm Spaces between target :5 mm

Scan the targets with the minimum depth of the view setting and identify the target visible at the shallowest point.



Example

Graph showing relationship between the actual depth of targets and the results of automatic measurement by an instrument.

Cyst targets



Target diameters 1, 2, 3, 4mm Target depths:

10, 30, 50, 70, 100, 150mm



Convex Probe Image

Find the deepest point where the non-resonance cyst targets in various diameters are visible against the background of speckle pattern.

This phantom is designed to allow the scanning from 4 sides, which provides more variety of checkup options.

String targets



Target diameter: 0.1mm Target depths: from 10 to 200mm Spaces between targets: 10mm Horizontal spaces between targets: 1, 2, 3, 4, 5, 10mm



Convex Probe Image

Measure the intervals between the targets and compare the result with actual distances.

This phantom is designed to allow the scanning from 4 sides, which provides more variety of checkup options.

Gray scale



Target depth 20mm Target diameter 10mm Spaces between targets : 12.5 mm

Echogenicity: 7 steps

Convex Probe Image



Example

A graph created by histogram function of an instrument.

Quality Control with N-365

- Before you start your periodic checking, wait more than 15 minutes after turning on the main power to your ultrasound scanner.
- Ensure to scan the phantom always at the same phantom temperature shown on the thermo mentor on the phantom wall, since the speed of the sound may depend on the temperature; the higher the temperatures, the slower the speed of sound.

To stabilize the temperature, it is recommended to keep the phantom under the room temperature same as where your going to perform the QC procedures for more than 6 hours before scanning.

• Ensure to hold you prove completely straight and vertically against the targets you intend to scan.

Record the image when the targets shown the smallest, or you recognize the highest resolution).

- At the time of your first checking, find the optimal gain so that all gray scales targets can be shown clearly and record the setting. Then use this gain for all other targets and use the same setting for your second periodic checking and after.
- To monitor the change of device across the ages, compare the latest date with the first time checking data created with the same setting and procedures.

Ultrasound QA phantom N-365 Multipurpose phantom Specifications

- Speed of sound1440m/sec (25 degrees C)
- Attenuation rate 0.57dB/cmMHZ
- Phantom size 190 × 220 × 70mm

Axial resolution, angular resolution Depth at 15, 30, 50mm Target diameter 0.05mm Spaces between targets 0.5, 1, 2, 3mm

Close range resolution Depth range of 1–10 mm depth (10 targets are embedded) Target diameter 0.1mm Spaces between targets :5 mm

String targets Target diameter: 0.1mm Target depths: from 10 to 200mm Spaces between targets: 10mm Horizontal spaces between targets: 1, 2, 3, 4, 5, 10mm

Cyst targets Target diameters: 1, 2, 3, 4mm Axial spaces between targets :10, 30, 50, 70, 100, 150mm

Gray scale Target depth 20mm Target diameter 10mm Spaces between targets :12.5 mm

Echogenicity: 7 steps