

ISO 9001:2008
ISO 13485:2003
CERTIFIED

"The Clinical Advantage"™

CATALOG 117
SUMMER 2012

BIODEX



Nuclear Medicine

SUPPLIES AND ACCESSORIES

BIODEX

www.biodex.com

1-800-224-6339

Int'l call 631-924-9000

"The Clinical Advantage"™

DEPEND ON BIODEX...
FOR ALL YOUR **PET** ACCESSORIES...
why settle for less!



Selecting the appropriate PET supplies
and accessories can take time...
choosing the right company to get them from shouldn't!

BIODEX

www.biodex.com

1-800-224-6339

Int'l 631-924-9000

“The Clinical Advantage”™




Biodex Medical Systems, Inc., has been providing customers with innovative products and service excellence for 60 years. Our dedicated employees work as a team to bring the promise of functional and elegant design to life.

It all begins with our belief in science-based solutions. Once the Biodex development team isolates a problem or requirement in the field we begin exploring possible responses. After a review of the literature and feedback from field luminaries, our engineers use cutting-edge technology to create products that are exceptionally functional, durable and user-friendly.



Having passed a design review and rigorous testing, the products are manufactured in our state-of-the-art, 120,000 sq. ft facility. We are committed to the continuous education of our employees. Our manufacturing team follows strict quality control guidelines in making each item ready for use. Biodex is certified for ISO 9001:2008, ISO 13485:2003 with specific products certified to UL 60601-1, EN 60601-1 and EN 60601-1-2, or other appropriate standards.

At Biodex, customer satisfaction drives every decision. Over 200 employees strive to keep our customers at the forefront of the art and science of medicine. No wonder so many world-class facilities call Biodex first.



If you are an existing customer, we thank you for your support and promise to continue our tradition of excellence. If you are a new customer, we invite you to try our products and realize they are engineered and manufactured to the highest standards. Dependability and exceptional performance allow you to forget about the product and focus on the patient – that’s your “Clinical Advantage.”™



• Syringe and Vial Shields	SECTION 1 4-18
• Shielding and Storage	SECTION 2 19-30
• Personal Radiation Protection	SECTION 3 31-32
• Lead-Lined Laboratory Furniture	SECTION 4 33-41
• Radiopharmacy	SECTION 5 42-63
• Dose Calibration	SECTION 6 64-75
• Thyroid Uptake System	SECTION 7 76-79
• Lung Ventilation Systems and Accessories	SECTION 8 80-92
• Sources	SECTION 9 93-95
• Phantoms	SECTION 10 96-103
• Radiation Detection	SECTION 11 104-113
• Decontamination and Sampling	SECTION 12 114-116
• Imaging Accessories	SECTION 13 117-120
• Signs, Labels and Tags	SECTION 14 121-123
• Numeric Index	SECTION 15 124
• Alphabetical Index	SECTION 16 125-126
• Conditions of Purchase	125

PRO-TEC® II SYRINGE SHIELD



- Tungsten shielding
- Lead glass window, easy to replace
- Lightweight
- Fits most disposable syringes

The Pro-Tec® II Syringe Shield reduces hand exposure. It is designed for the clinician who has previously found it difficult to use conventional syringe shields. The Pro-Tec® II is thin, lightweight and easy to use.

The barrel of the shield is constructed of .06" (1.5 mm) thick tungsten that will reduce radiation exposure from Tc-99m by more than 99% attenuation for Tc-99 tested with TLD chips. A 5.05 density (1.5 mm lead equivalency) lead glass window provides protection and visibility. A white reflective surface on the shield interior improves viewing of the syringe's markings and fluid content. A clip secures the lead glass to the tungsten barrel and makes it easy to replace scratched or broken glass. A thumbscrew holds syringes firmly in place.

Pro-Tec® II Syringe Shields accommodate the standard sized 1 cc to 12 cc syringes.



Replacing scratched or broken glass is simple.



White coated interior makes it easy to read syringe contents.

Pro-Tec® II Syringe Shields with lead glass window:

- 007-700** Syringe Shield, 1 cc
- 007-800** Syringe Shield, 3 cc
- 007-900** Syringe Shield, 5 cc
- 007-909** Syringe Shield, 10 cc
- 007-913** Syringe Shield, 12 cc

U.S. Patent #4,062,353

Pro-Tec® II Syringe Shields without lead glass window:

- 007-701** Syringe Shield, 1 cc
- 007-801** Syringe Shield, 3 cc
- 007-901** Syringe Shield, 5 cc
- 007-911** Syringe Shield, 10 cc

Replacement Glass for Pro-Tec® II Syringe Shields:

- 127-700** Syringe Shield Replacement Glass, 1 cc

For 007-700

- 127-789** Syringe Shield Replacement Glass,
3 cc and 5 cc

For 007-800, -900, -835, -837, -903, -915

- 127-999** Syringe Shield Replacement Glass, 10 cc

For 007-909, -916


Note: Syringe Shields available for a selection of international syringes. Contact Biodex at 631-924-9000 or e-mail sales@biodex.com.

Sharps Container Shields



View our full line of Sharps Container Shields
www.biodex.com/syringeshields

PRO-TEC® III SYRINGE SHIELD



Design accommodates most 1 cc, 2.5 cc, 3 cc, 5 cc, 10 cc and 20 cc syringes

- Tungsten shielding
- Lead glass window
- Lightweight
- Fits most disposable syringes
- Safe-T-Lock design reduces exposure with faster handling

The Pro-Tec® III Syringe Shield reduces hand exposure and is convenient to use.

The barrel of the shield is constructed of .08" (2 mm) thick tungsten that will reduce radiation exposure from Tc-99m by more than 99% attenuation for Tc-99 tested with TLD chips. A 5.05 density (1.5 mm lead equivalency) lead glass window provides protection and visibility. A white reflective surface on the shield interior improves viewing of the syringe's markings and fluid content. A bevel around the lead glass helps protect it from scratching or breaking.

The Safe-T-Lock design provides increased protection by allowing minimal handling. The Safe-T-Lock design grips and secures the syringe upon insertion and will release the syringe at the press of a button. Disposing used syringes is easy; invert the syringe shield over a sharps container, press the release button and the syringe freely disengages.

Pro-Tec® III Syringe Shields accommodate the standard sized 1 cc, 3 cc, 5 cc, 10 cc and 20 cc syringes.



Safe-T-Lock design quickly engages or releases syringe


Pro-Tec® III Syringe Shields:

- 007-723** Syringe Shield, 1 cc (BD luer lock)
- 007-734** Syringe Shield, 1 cc (press fit)
- 007-755** Syringe Shield, 2.5 cc
- 007-735** Syringe Shield, 3 cc
- 007-736** Syringe Shield, 5 cc
- 007-738** Syringe Shield, 10 cc
- 007-747** Syringe Shield, 20 cc

Note: Syringe Shields available for a selection of international syringes. Contact Biodex at 631-924-9000 or e-mail sales@biodex.com.

Shielded Syringe Carriers

Two sizes, each offered in .125" and .25" lead shielding for added protection



See page 12 for details.

PRO-TEC® IV SYRINGE SHIELD

*High density lead glass with a full 360° view -
all the features you want combined in a single shield*



- *High density (5.6) lead glass shielding, 360-degree barrel view*
- *Lightweight*
- *Fits most disposable syringes*
- *Safe-T-Lock design reduces exposure with faster handling*

The Pro-Tec® IV Syringe Shield reduces hand exposure, maximizes the viewing area and is convenient to use. The barrel of the shield is constructed of optically clear high density (5.6) lead glass. The end of the lead glass barrel is tapered so that the contents of the syringe can be viewed. The 360° view and taper make the venipuncture faster and easier. The high density lead glass significantly reduces radiation exposure from Tc-99m by more than 99% attenuation for Tc-99 tested with TLD chips.

The Safe-T-Lock design provides increased protection by allowing minimal handling. The Safe-T-Lock design grips and secures the syringe upon insertion and will release the syringe at the press of a button. Disposing used syringes is easy; invert the syringe shield over a sharps container, press the release button and the syringe freely disengages.

Pro-Tec® IV Syringe Shields accommodate standard sized 1 cc to 10 cc syringes.

Pro-Tec® IV Syringe Shields:

- 007-670** Syringe Shield, 1 cc
- 007-671** Syringe Shield, 2 cc
- 007-675** Syringe Shield, 3 cc
- 007-680** Syringe Shield, 5 cc
- 007-685** Syringe Shield, 10 cc

Note: Syringe Shields available for a selection of international syringes. Contact Biodex at 631-924-9000 or e-mail sales@biodex.com.

HIGH DENSITY LEAD GLASS SYRINGE SHIELD



Offers 360° syringe visibility; model 007-635 shown

The High Density Lead Glass Syringe Shield reduces hand exposure from Tc-99m by more than 99%, allows a large viewing area and is easy to use. The barrel of the syringe shield is constructed of optically clear high density (5.6) lead glass and offers complete 360° visibility. The end of the barrel is tipped with steel to protect it from breaking or scratching. This lightweight shield features quick and smooth syringe insertion with an O-ring seal and an anti-roll cap.

High Density Lead Glass Syringe Shields:

Accommodates luer and non-luer lock syringes


- 007-620** Syringe Shield, 1 cc
- 007-635** Syringe Shield, 3 cc
- 007-652** Syringe Shield, 5 cc and 6 cc
- 007-612** Syringe Shield, 10 cc and 12 cc

Quality Shielding & Storage

There is a difference between Safe and Biodex Safe.™

www.biodex.com/shielding

MYOVIEW SYRINGE SHIELD



Pro-Tec III shown



Pro-Tec IV shown

These Pro-Tec® Syringe Shields accommodate the HSW brand of syringes that are recommended for drawing and administering Myoview. Lightweight and convenient to use, these syringe shields will reduce hand exposure.

The Safe-T-Lock design provides increased protection by allowing minimal handling.

007-755 Syringe Shield, Pro-Tec III, 2.5 cc

007-671 Syringe Shield, Pro-Tec IV, 2 cc
Suitable for Myoview

QUADRAMET SYRINGE SHIELD

Pro-Tec IV Syringe Shield used for Quadramet



- *High density (5.6) lead glass shielding reduces exposure from Samarium 153 by 99%*
- *360° barrel view*
- *Lightweight*
- *Fits most disposable syringes*
- *Safe-T-Lock design reduces exposure with faster handling*

The Pro-Tec® IV 5 cc Syringe Shield is an ideal shield for the administration of Quadramet®. It reduces hand exposure, maximizes the viewing area and is convenient to use. The barrel of the shield is constructed of optically clear high-density (5.6) lead glass. The end of the lead glass barrel is tapered so that the contents of the syringe can be viewed. The 360° view and taper make venipuncture faster and easier. The high-density lead glass significantly reduces radiation exposure from Sm-153 by greater than 99%.

The Safe-T-Lock design provides increased protection by allowing minimal handling. It grips and secures the syringe upon insertion and releases the syringe at the press of a button. Disposing used syringes is easy; invert the syringe shield over a sharps container, press the release button and the syringe freely disengages.

007-680 Syringe Shield, Pro-Tec IV, 5 cc
Suitable for Quadramet

COLOR-CODED SYRINGE SHIELD

Use in combination with color-coded vial shields



- *High density (5.6) lead glass shielding, 360-degree barrel view*
- *Lightweight*
- *Fits most disposable syringes*
- *Color-coded for easy identification, coordinated with color-coded vial shields*

This 5.6 density lead glass syringe shield features an optically clear barrel that permits an unobstructed, 360° view of syringe contents. The tapered tip is designed for maneuverability making venipuncture faster and easier.

A color-coded acrylic collar with a thumbscrew holds the syringe firmly in place. The syringe shield weighs only three ounces.

Color-Coded, 3 cc, Syringe Shields:

007-241 Syringe Shield, Red

007-242 Syringe Shield, Purple

007-243 Syringe Shield, Orange

007-244 Syringe Shield, Blue

007-245 Syringe Shield, Brown

007-246 Syringe Shield, Yellow



For our full line of Color-Coded Vial Shields, click into our website:

www.biodex.com/syringeshields

PRO-TEC® PET SYRINGE SHIELD



- Constructed of .34" thick (9 mm) tungsten, attenuates FDG F-18 by 88%
- Available with or without a high density lead glass window
- Fits most disposable syringes

The Pro-Tec® PET Syringe Shield reduces hand exposure from syringes containing 511 keV radionuclides. The barrel of the shield is constructed of .34" thick (9 mm) tungsten that attenuates FDG F-18 by 88%.

The syringe shield is offered with or without a high density (5.6) flush mounted lead glass window that provides protection and visibility. A white reflective surface on the shield interior improves viewing of the syringe's markings and fluid content. A thumbscrew holds syringes firmly in place.

Pro-Tec® PET Syringe Shields accommodate the standard sized 1 cc, 3 cc, 5 cc and 10 cc syringes. The Manual Dose Injector is an ideal companion, providing both additional shielding and distance.

SPECIFICATIONS:

Shielding: .34" thick (9 mm) tungsten

Lead Glass: 5.6 density

Weight:

007-973 & 007-985: 1.4 lb (.64 kg)

007-975 & 007-990: 1.7 lb (.83 kg)

007-980 & 007-995: 2.3 lb (1.05 kg)

Pro-Tec® PET Syringe Shields with lead glass window:

007-973 Syringe Shield, 3 cc

007-975 Syringe Shield, 5 cc

007-980 Syringe Shield, 10 cc

Pro-Tec® PET Syringe Shields without lead glass window:

007-983 Syringe Shield, 1 cc

007-985 Syringe Shield, 3 cc

007-990 Syringe Shield, 5 cc

007-995 Syringe Shield, 10 cc

Replacement Glass:

007-974 Glass, Replacement

For Pro-Tec PET and Gaard Lock

PET Syringe Shields 007-973, 007-975,

007-980, 007-716, 007-717 and 007-718

Note: Syringe Shields available for a selection of international syringes. Contact Biodex at 631-924-9000 or e-mail sales@biodex.com.

Z-PET SYRINGE SHIELD



Extra thick wall for extra protection

- Constructed of .55" thick (14 mm) tungsten, attenuates FDG F-18 by 97%

The Z-PET Syringe Shield greatly reduces hand exposure from syringes containing 511 keV radionuclides. The barrel of the shield is constructed of .55" thick (14 mm) tungsten that attenuates FDG F-18 by 97%. The shield accommodates standard 5 cc syringes. The Manual Dose Injector is an ideal companion, providing both additional shielding and distance.

SPECIFICATIONS:

Dimensions: 2.75" l x 1.7" dia (7 x 4.3 cm)

Shielding: .55" thick (14 mm) tungsten

Weight: 3.7 lb (1.7 kg)

007-945 Syringe Shield, Z-PET, 5 cc*

*Z-PET Syringe Shield was conceived by Michael Zimmer, Ph.D.

Lead-Lined Unit Dose Cabinet




Designed for PET Hot Labs with limited space.

See page 40 for details.

GAARD LOCK™ PET SYRINGE SHIELD



Available with lead glass window



Twist to lock syringe in place with one hand

- *Unique flange locking design reduces exposure with faster handling*
- *Constructed of .34" thick (9 mm) tungsten, attenuates FDG F-18 by 88%*
- *Available with or without a high density lead glass window*

The Gaard Lock™ PET Syringe Shield reduces hand exposure from syringes containing 511 keV radionuclides FDG F-18. The barrel of the shield is constructed of .34" thick (9 mm) tungsten that attenuates FDG F-18 by 88%.

The syringe shield is offered with or without a high density (5.6) flush mounted lead glass window that provides protection and visibility. A white reflective surface on the shield interior improves viewing of the syringe's markings and fluid content.

The syringe shield features a unique flanged locking mechanism that speeds the loading and unloading of the syringe, further reducing hand exposure.

Gaard Lock PET Syringe Shields accommodate the standard sized 3 cc, 5 cc and 10 cc syringes. The Gaard Lock PET Syringe Shield fits the Manual Dose Injector.

SPECIFICATIONS:

Shielding: .34" tungsten (9mm)

Lead glass: 5.6 density

007-711 & 007-716 Gaard Lock Syringe Shield, 3 cc

Dimensions: 2.9" l (74 mm)

Weight: 1.7 lb (.77 kg)

007-712 & 007-717 Gaard Lock Syringe Shield, 5 cc

Dimensions: 3" l (76 mm)

Weight: 2 lb (.91 kg)

007-713 & 007-718 Gaard Lock Syringe Shield, 10 cc

Dimensions: 3.6" l (91 mm)

Weight: 3 lb (1.4 kg)

Gaard Lock™ PET Syringe Shields with lead glass window:

007-716 Syringe Shield, 3 cc

007-717 Syringe Shield, 5 cc

007-718 Syringe Shield, 10 cc

Gaard Lock™ PET Syringe Shields without lead glass window:

007-711 Syringe Shield, 3 cc

007-712 Syringe Shield, 5 cc

007-713 Syringe Shield, 10 cc

Replacement Glass:

007-974 Glass, Replacement

For Pro-Tec PET and Gaard Lock
PET Syringe Shields 007-973, 007-975,
007-980, 007-716, 007-717 and 007-718

BETA SYRINGE SHIELD



- *Lead and plastic shielding, designed for Beta*
- *Clear plastic window*
- *Fits most disposable syringes*
- *Ideal for Zevalin Y-90*

This Beta Syringe Shield reduces hand exposure from syringes containing Strontium 89, P-32, and other beta-emitting radiopharmaceuticals. The barrel of the syringe shield is constructed of .375" clear plastic with a .062" thick embedded lead lining, which attenuates Beta emission and errant bremsstrahlung. The lead lining is interrupted at the viewing window for clear visibility when drawing and administering the dose. A thumbscrew holds syringes firmly in place.

For pharmacies preparing a dose of Zevalin, both a 1 cc and 10 cc syringe shield are offered for the two stages of preparation. Use the 1 cc syringe shield to draw Zevalin Y-90 directly from the pharmacy's supplied vial shield, count in a dose calibrator and inject to the reaction vial. Use the 10 cc syringe shield to draw the prepared dose from the reaction vial for patient administration.

The Manual Dose Injector is an ideal companion, providing both additional shielding and distance.

SPECIFICATIONS:

Dimensions: 1.625" dia x 3.25" l (4.1 x 8.3 cm)

Shielding:

- 007-956: .422" plastic, .062" lead, .125" plastic
(1.07 x .16 x .32 cm) with .86" (2.2 cm) plastic window
- 007-957: .25" plastic, .062" lead, .125" plastic
(.64 x .16 x .32 cm) with .70" (1.8 cm) plastic window

Weight:

- 007-956: 10.4 oz (297 g)
- 007-957: 9.6 oz (274 g)

- 007-956** Syringe Shield, Beta, 1 cc
- 007-957** Syringe Shield, Beta, 10 cc

Related:

- 007-997** Dose Injector, Manual, 1" lead

How the use of syringe shields can reduce hand exposure

*Health Physics, Department of Clinical Physics and Bio-Engineering,
Western Infirmary, Glasgow G11 6NY, UK.*

Mark.Whitby@NorthGlasgow.Scot.NHS.UK

Staff preparing and injecting radiopharmaceuticals in hospitals may receive significant radiation doses to their hands. These doses may be high enough to warrant that they be classified as radiation workers. The influence of local shielding on finger doses has been investigated. Staff preparing radioactive liquids in a radionuclide dispensary and drawing up and injecting radiopharmaceuticals in a nuclear medicine department have been studied. Measurements have been recorded with an electronic extremity dose monitor, an advanced extremity gamma instrumentation system (AEGIS), worn near to the finger tip. The electronic dosimeter allows the pattern of doses received during different procedures to be determined. Doses received for individual manipulations during many routine sessions have been recorded for different staff members. Dose distributions around shielded vials and syringes have also been measured using AEGIS. In the radionuclide dispensary the vials from which radioactive liquids are dispensed are held in tungsten shields, whereas in nuclear medicine simple lead pots are used. Syringe shields are employed for some parts of dispensing and patient injections. Data on dose distributions have been used in interpretation of results from monitoring. Use of syringe shields during dispensing reduced the finger dose by 75-85%. The peaks in dose rate were 60% lower, and periods of exposure to high dose rates were reduced in length by a third because of the restriction in the region of high dose rate. The extremity doses to staff dispensing and injecting radiopharmaceuticals in nuclear medicine were of similar magnitude. Doses received during dispensing varied from 10 to 555 microGy depending upon whether the vial containing the radiopharmaceutical was directly handled or not. Doses received from individual injections varied from 1 to 150 microGy depending on the degree of difficulty experienced during the injection.

PMID: 12729420 [PubMed - indexed for MEDLINE]

For complete article, click into our website
www.biodex.com/syringeshields

BIODEX

SHIELDED SYRINGE HOLDER

Accommodates shielded and unshielded syringes



For added protection, store loaded syringes in the Shielded Syringe Holder

SPECIFICATIONS:

Dimensions: 6.5" h (16.5 cm)

I.D.: .84" dia. x 5.6 h (2.1 x 14.3 cm)

Lead Shielding: .25" - .5" thick (.64 cm - 1.3 cm)

Accommodates Syringe Shields:

Pro-Tec II: 1 cc, 3 cc, 5 cc

Pro-Tec III: 1 cc, 3 cc, 5 cc

Pro-Tec IV: 1 cc

Weight: 6 lb (2.7 kg)

009-205 Syringe Holder, Shielded

UNIT DOSE PIG WALL RACK



Holds pigs up to 2.25" diameter.

SPECIFICATIONS:

Dimensions: 17" w x 18" h (43.2 x 45.7 cm)

008-400 Wall Rack, Unit Dose Pig

Improve lab safety, efficiency and organization with the new Unit Dose Pig Wall Rack.

The wall rack improves work space with its pigeonhole design and reduces unnecessary handling. Unit doses can be identified at a glance. The rack's sturdy construction will hold up to 25 unit dose pigs.

SYRINGE SHIELD HOLDER

Protect your investment



Don't let your syringe shields roll around on the counter. The Syringe Shield Holder offers a means of protecting syringe shields from scratches or misplacement while freeing up extra work space. The Syringe Shield Holder will support up to eight

shields and is counterbalanced to prevent tipping. Know exactly where syringe shields are when you need them.

SPECIFICATIONS:

Dimensions: 7.5" w x 4" depth x 6" h (19 x 10.2 x 15.2 cm)

Weight: 5 lb (2.3 kg)

007-999 Syringe Shield Holder

SYRINGE RECAPPER

Syringe recapping device



Don't put yourself at risk with an accidental needle stick! The Syringe Recapper is a safe and inexpensive way to protect yourself when recapping a used syringe. Used either hand-held or placed on a flat surface, such as a procedure tray, the Recapper is made of a lightweight plastic that is easily carried anywhere.

008-300 Syringe Recapper

3.25" l x 2.75" w (8.3 x 7 cm)

SHIELDED SYRINGE CARRIERS

Two sizes, each offered in .125" and .25" lead shielding for added protection



Shielded Syringe Carriers, small and large

Shielded Syringe Carriers reduce exposure while storing or transporting radioactive material. The overlapping lid design with snap-latch closure prevents streaming. There are two sizes and thicknesses of lead to choose from. The ends of the carriers are double thick to reduce the exposure from the ends of syringes.

SPECIFICATIONS:

001-179 Shielded Syringe Carrier, Large

Dimensions:

I.D.: 8.25" l x 3" w x 2.9" h (21 x 7.6 x 7.4 cm)

O.D.: 9.5" l x 4.4" w x 3.5" h (24 x 11.2 x 8.9 cm)

Lead Shielding:

Sides, top and bottom: .125" thick (.32 cm)

Ends: .25" thick (.64 cm)

Weight: 11.3 lb (5.1 kg)

001-181 Shielded Syringe Carrier, Small

Dimensions:

I.D.: 8" l x 1.9" w x 1.97" h (20.3 x 4.8 x 5 cm)

O.D.: 9.25" l x 3.4" w x 2.6" h (23.5 x 8.6 x 6.6 cm)

Lead Shielding:

Sides, top and bottom: .125" thick (.32 cm)

Ends: .25" thick (.64 cm)

Weight: 7.5 lb (3.4 kg)

001-182 Shielded Syringe Carrier, Small

Dimensions:

I.D.: 7.5" l x 1.7" w x 1.7" h (19 x 4.3 x 4.3 cm)

O.D.: 9.25" l x 3.4" w x 2.6" h (23.5 x 8.6 x 6.6 cm)

Lead Shielding:

Sides, top and bottom: .25" thick (.64 cm)

Ends: .5" thick (1.3 cm)

Weight: 11 lb (4.9 kg)

001-180 Shielded Syringe Carrier, Large

Dimensions:

I.D.: 7.7" l x 2.6" w x 2.7" h (19.6 x 6.6 x 6.9 cm)

O.D.: 9.5" l x 4.4" w x 3.5" h (24 x 11.2 x 8.9 cm)

Lead Shielding:

Sides, top and bottom: .25" thick (.64 cm)

Ends: .5" thick (1.3 cm)

Weight: 17 lb (7.7 kg)

Shielded Syringe Carrier, .125" lead:

001-181 Syringe Carrier, Small

001-179 Syringe Carrier, Large


Shielded Syringe Carrier, .25" lead:

001-182 Syringe Carrier, Small

001-180 Syringe Carrier, Large

MANUAL DOSE INJECTOR

Reduces hand exposure



Manual Dose Injector further reduces hand exposure by adding shielding and distance during injection (shown with PET Syringe Shield).

The Manual Dose Injector is designed to reduce exposure when capping or uncapping syringes and injecting radionuclides. It is especially useful to protect fingers when administering high-energy nuclides. The injector is also useful when the injection protocol is a slow infusion as with Zevalin™.

Here's how it works:

1. **Remove the needle** - Place a shielded syringe into the Manual Dose Injector. Use forceps to hold the needle hub. Turn the lead base and the needle detaches.
2. **Attach a butterfly, etc.** - Use forceps to hold the butterfly hub. Place the hub onto the syringe luer lock. Turn the lead base and the butterfly is attached.
3. **Inject** - Use the plunger on the Dose Injector to push the syringe plunger. The plunger can also be drawn back to verify that the needle is still in the vein. This is important for Zevalin™.

The Manual Dose Injector features a spring-loaded clamp to hold the syringe shield barrel firmly in place. There is a 1" thick lead shield to protect hands. The syringe plunger is held by the injector shield so that the syringe plunger can move forward or back. The Dose Injector can be easily transported or placed wherever administrations are performed. It can also be used in pharmacies to remove needles or to recap.

SPECIFICATIONS:

Dimensions: 11.5" w x 5" depth x 10.75" h (29 x 12.7 x 27 cm)

Shielding:

Rotation Shield: 1" thick lead (2.5 cm)

Plunger: 0.5" thick tungsten (1.25 cm)

Finish: Powder coat

Weight: 8.5 lb (3.8 kg)

007-997 Dose Injector, Manual, 1" lead

Note: Syringe Shield sold separately

GRASEBY INJECTION PUMP SHIELD



Lead glass window offers 5.6 density with clear visibility of syringe.



The shield cover, when open, acts as a protective barrier; stand behind the pump to insert the syringe.

- 0.50" lead shielding
- Lead-glass viewing window
- Accommodates Graseby Syringe Pumps

The Injection Pump Shield is designed to reduce technologist exposure during the injection of radionuclides that require a very slow infusion protocol using Graseby #3100 or #3400 Syringe Pumps. The shield is constructed with 0.5"-thick lead to attenuate I-131. A lead glass viewing window allows the technologist to see the syringe during the infusion. The cover can slide lengthwise to reveal the syringe stopcock for adjustments, without exposing the syringe.

SPECIFICATIONS:

Dimensions: 20.3" l x 11.4" w x 8" h (51.6 x 29 x 20.3 cm)

Lead Shielding: 0.50" thick (1.3 cm)

Lead Glass Window:

Dimensions: 4.25" w x .50" h x .38" thick (10.8 x 1.27 x .97 cm)

Density: 5.6

Construction:

Exterior: Steel

Interior: Lead

Finish: Powder coat

Weight: 90 lb (41 kg)

Attenuation:

for Y-90 greater than 99%

for I-131 greater than 99%

for F-18 greater than 99%

For complete test results, visit www.biodex.com/syringeshields

007-040 Shield, Injection Pump

Note: Injection Pump sold separately

INJECTION STAND



Convenient, adjustable and portable for safe injections; includes a Lead-Lined Multi-Syringe Holder for convenient radionuclide injections



The Injection Stand allows fast, comfortable arm positioning for radionuclide injections. The clear plastic armrest rotates 180° facilitating convenient angling while the cradle design holds the patient's arm firmly in place. A utility tray sits adjacent to the armrest to place various supplies and includes a shielded Multi-Syringe Holder, which accommodates up to four syringes.

The stand is height adjustable to accommodate patient comfort. Smooth rolling casters allow the

stand to roll easily into position or out of the way for storage when not in use. The stand is structurally balanced to help prevent tipping.

SPECIFICATIONS:

Height Adjustable: 29" h to 44" h (74 to 112 cm)
 Stand Base: 17.5 x 17.5 (44.5 x 44.5 cm)
 Construction: Stable, chrome plated tubular steel
 Shipping Weight: 21 lb (10 kg)

135-022 Stand, Injection

INJECTION CHAIR



214-200 Adjustable pivot swing arm for ultimate comfort and positioning can be mounted on either side of the chair.



214-220 Injection Chair is also available with an easy to clean, one piece plastic drawer. Drawer and arm can be mounted on either side of the chair.

Designed specifically for injections and blood drawing, this chair is a comfortable solution for patient seating and positioning. The flat surface of the armrest prevents backbending of the elbow and subsequent flattening of the vein. The armrest is height adjustable to accommodate all patient sizes with a knob adjustable double pivot swing arm. The seat itself is one piece of plastic which makes it easy to clean. A stiffener bias limits backrest flexibility to ensure stable seating. Choose with or without accessory drawer.

SPECIFICATIONS:

Dimensions:

214-200: 23" l x 33" w (58.5 x 79 cm)

214-220: 23" l x 43" w (58.5 x 109 cm)

Seat Dimensions: 16" l x 17" w x 20" h (40.6 x 43 x 51 cm)

Arm Height (adjustable): 27" h to 33" h (69 to 84 cm)

Width between arms: 20" (51 cm)

Construction (frame): Heavy steel tubing 1.25" (3 cm) square, with chromed steel uprights for stability

Weight:

214-200: 32 lb (14.5 kg)

214-220: 35 lb (15.9 kg)

Patient Weight Capacity: 250 lb (113.4 kg)

214-200 Chair, Injection

214-220 Chair, Injection with Storage Drawer

VEINLITE LED

No more guesswork – this vein transilluminator accurately locates veins



- *Reduces error when injecting radioactivity in patients*
- *Reduces multiple needle sticks*
- *Reduces time required for vein access*

Veinlite LED helps locate the vein, and then secures it in place by stretching the skin for easy needle insertion. There are 24 diodes in two LED colors for vein imaging: 12 orange diodes are used for superficial veins and 12 red diodes are used for deeper veins or for dark-skinned patients. When time is of the essence, Veinlite LED eliminates the guesswork in locating the vein when injecting radioactivity into patients.

Disposable plastic covers are available to eliminate cross contamination.

SPECIFICATIONS:

Dimensions: 3.75" l x 2" w x 1" thick (9.5 x 5.5 x 2.4 cm)

Power: Rechargeable lithium-ion battery

Weight: 3 oz (70 gm)

Warranty: One year

007-275 Veinlite LED

Includes: 50 disposable plastic covers

Replacement:

007-276 Covers, Disposable Plastic, 50/pkg

INJECTION / RESTING CHAIR

A comfortable rest between injection and imaging



Most protocols require that prior to imaging, a patient rest for up to one hour after an injection of FDG F-18. Typical injection chairs are not designed for comfort, making rest difficult.

The Injection/Resting Chair is a comfortable resting spot for patients to wait between injection and imaging. The chair achieves three recline positions and features an infinitely adjustable back. For patient safety, deep recline and Trendelenburg positions can only be achieved by the attendant. A quick release handle that spans the entire width of chair back allows effortless position change from either side of the chair.

SPECIFICATIONS:

Dimensions: 35" w x 49" h (89 x 125 cm)

Seat: 20" depth x 25" w x 21" h (51 x 64 x 53 cm)

Back: 25" w x 34" h, above seat (64 x 86 cm)

Upholstery: Blueridge, fire retardant

Frame: Welded powder-coated steel

Casters: Four 5" heavy-duty with positive locking swivel and wheel brakes

Warranty: Five-year on frame

Patient Capacity: 450 lb (204 kg)

Weight: 116 lb (52.6 kg)

214-210 Chair, Injection / Resting

SHIELDED STORAGE CONTAINERS

For beta and gamma radiation



Eliminate clutter and reduce safety hazards with convenient Shielded Storage Containers



Molded plastic liner with threaded lid fits Storage Container 050-250.

SPECIFICATIONS:

050-200 Shielded Storage Container, Gamma, Small

Dimensions: 6.5" h x 5" dia (16.5 x 12.7 cm)

Lead Shielding: .125" thick (.32 cm)

Weight: 7 lb (3.2 kg)

050-250 Shielded Storage Container, Gamma, Large

Dimensions: 7" h x 6" dia (17.8 x 15.2 cm)

Lead Shielding: .125" thick (.32 cm)

Weight: 9 lb (4.1 kg)

050-205 Shielded Storage Container, Beta/Gamma, Small,

Dimensions: 6.5" h x 5" dia (16.5 x 12.7 cm)

Lead Shielding: .25" thick (.64 cm)

Aluminum Shielding: .0625" thick (1.6 mm)

Weight: 12 lb (5.4 kg)

050-200 Storage Container, Gamma, Small,
.125" lead

050-250 Storage Container, Gamma, Large,
.125" lead

050-205 Storage Container, Beta/Gamma, Small,
.25" lead/.0625" al

Related:

007-007 Liner, Molded Plastic, 12/pkg


040-315 Liner, Poly Bag, 50/pkg

Measures 6" x 3" x 15"

Fits 050-200, 050-250 and 050-205

SHARPS CONTAINER SHIELD

For beta and gamma radiation



This Sharps Container Shield is a simple way to dispose of used syringes that may be contaminated with low-energy gamma and beta radiation residue. The shield is constructed of steel lined with .125" (.32 cm) lead and .0625" (1.6 cm) aluminum. The shielding combination attenuates gamma radiation, beta radiation and errant bremsstrahlung.

The shield features a top with sliding port and a swing hand for easy transport. It is designed to accommodate small and medium Monoject Sharps containers.

SPECIFICATIONS:

Shield accommodates small and medium Monoject Sharps containers

Dimensions: 7.5" l x 11.75" w x 10" h (19 x 30 x 25 cm)

I.D.: 7" l x 11.25" w x 9.5" h (18 x 29 x 24 cm)

Lead Shielding: .125" thick (.32 cm)

Aluminum Shielding: .0625" thick (1.6 mm)

Security: Key-locked

Finish: Powder coat

Weight: 48 lb (22 kg)

039-330 Sharps Container Shield, Beta/Gamma

For Monoject Containers 039-338 and 039-387

Related:

039-338 Monoject Container, 4 qt., 10/pkg

Fits 039-325, 039-326, 039-350 and 039-330

039-387 Monoject Container, 8 qt., 10/pkg

Fits 039-325, 039-326 and 039-330

SHARPS CONTAINER SHIELDS

Key-lock design meets OSHA standards



Transfer the full Sharps Container to the Decay Module for syringe decay storage.

Sharps Container Shields are a simple, safe and convenient way to dispose of used syringes that may contain low-energy gamma radiation residue. The shields are constructed of steel lined with .125" (.32 cm) lead. They feature a top with sliding port and a swing handle that allows easy transport.


Two models are available for one sharps container:

- Model 039-325 accommodates small and medium Monoject Sharps Containers.
- Model 039-335 accommodates Sharpstainer* and Medium Monoject Sharps containers.

One model is available for two sharps containers:

- Model 039-350 accommodates two small Monoject Sharps Containers simultaneously for decay rotation.

The sharps containers are stacked inside the shield. When the top container is filled, the bottom decayed container is removed for proper waste disposal. The top container is then shifted to the lower position for decay while a fresh container is placed on top for immediate use. Spring clamps hold containers securely in place.



Dual Container Sharps Shield houses two containers for decay rotation.

SPECIFICATIONS:

039-325 Sharps Container Shield
Shield uses 039-338 and 039-387

Monoject Sharps containers

Dimensions: 7.5" l x 11.75" w x 10" h
(19 x 30 x 25 cm)

I.D.: 7" l x 11.25" w x 9.5" h
(18 x 29 x 24 cm)

Lead Shielding: .125" thick (.32 cm)

Security: Key-locked

Finish: Powder coat

Weight: 48 lb (22 kg)



039-335 Sharps Container Shield
Shield uses 039-341 and

039-387 Monoject Sharps containers

Dimensions: 7.5" l x 12.5" w x 13" h
(19 x 32 x 33 cm)

I.D.: 6.875" l x 10.75" w x 12.5" h
(17.5 x 27 x 32 cm)

Lead Shielding: .125" thick (.32 cm)

Security: Key-locked

Finish: Powder coat

Weight: 53 lb (24 kg)



039-350 Sharps Container Shield, Dual

Shield uses two 039-338 Monoject Sharps containers

Dimensions: 7.375" l x 12" w x 16.75" h
(19 x 30.5 x 42.5 cm)

I.D.: 7" l x 11.5" w x 14.875" h
(18 x 29 x 38 cm)

Lead Shielding: .125" thick (.32 cm)

Finish: Powder coat

Weight: 60 lb (27 kg)



Sharps Container Shields:

039-325 Sharps Container Shield, .125" lead

For Monoject containers 039-338 and 039-387

039-335 Sharps Container Shield, .125" lead

For Sharpstainer container 039-341 and
B-D Guardian container 039-321*

039-350 Sharps Container Shield, Dual, .125" lead

For Monoject containers 039-338

Related:

039-338 Monoject Container, 4 qt., 10/pkg

Fits 039-325, 039-326, 039-350 and 039-330

039-387 Monoject Container, 8 qt., 10/pkg

Fits 039-325, 039-326, 039-335 and 039-330

039-341 Sharpstainer* Container, 6.2 qt., (#182),
12/pkg

Fits 039-335

**Note: Sharpstainer formerly Winfield*

HIGH-ENERGY SHARPS CONTAINER SHIELD



This Sharps Container Shield is a simple, safe and convenient solution for disposal of used syringes that may be contaminated with high-energy isotopes such as I-131. The shield is constructed of steel lined with .5" of lead (1.3 cm).

The shield features a hinged top with a sliding port and side handles which allow for easy transport. It will accommodate both small and medium Monoject Sharps containers.

SPECIFICATIONS:

Dimensions:

I.D.: 7" l x 11.38 w x 11.25 h

O.D.: 9.19" l x 15.25" w x 13.31" h

Lead Shielding: .5" thick (1.3 cm)

Security: Key-locked

Finish: Powder coat

Weight: 160 lb (72.3 kg)

039-326 Sharps Container Shield, .5" lead
For Monoject Containers 039-338 and 039-387

Related:


039-338 Monoject Container, 4 qt., 10/pkg
Fits 039-325, 039-326, 039-350 and 039-330

039-387 Monoject Container, 8 qt., 10/pkg
Fits 039-325, 039-326 and 039-330

HIGH-ENERGY PET SHARPS CONTAINER SHIELD



Designed to stand on, or recess into, a countertop



Sharps Container Shield shown recess mounted in countertop.

The Sharps Container Shield for PET is a simple, safe and convenient solution for disposal of syringes that have been contaminated with high-energy radionuclides. The shield is constructed of steel and lined with 1" thick (2.5 cm) lead.

The shield is designed to be used with 039-413 sharps containers. It features a lockable sliding cover for container removal and a hinged top door for syringe disposal.

The shield can stand independently or can be recessed into a cabinet or countertop.

SPECIFICATIONS:

Dimensions: 12" h x 8.75" dia (30.5 x 22.2 cm)

Lead Shielding:

Sides and Bottom: 1" thick (2.5 cm)

Rotating Cover: .875" thick (2.2 cm)

Hinged Door: .625" thick (1.5 cm)

Security: Key-locked

Weight: 175 lb (79.4 kg)

Shipping Weight: 222 lb (100.6 kg)



039-412 Sharps Container Shield, PET, 1" lead
Uses one 039-413 Sharps Container

Related:

039-413 Sharps Container, 3.2 qt., 30/pkg
Fits 039-412


DELUXE L-BLOCK SHIELD

Ideal for PET pharmacies



- 2.4" lead shielding in front and in base
- 8" x 8" x 4" easily adjustable lead glass window
- Optional Lead Brick Cave for complete lateral shielding

The Deluxe L-Block Shield is designed for managing large quantities of high-energy radionuclides. The shield is constructed of 2.4" thick lead encased in steel, and features a large 8" x 8" x 4" lead glass window. A convenient lever allows quick adjustment of window to optimal angle for any user and procedures. The optional 042-417 Lead Brick Cave fits neatly into the sides of the vertical section to provide lateral shielding around the full perimeter of the L-Block's base. A special plate with a hex-shaped recess is mounted on the base to facilitate one-handed loading and unloading of dose pigs incorporating hex-shaped bottoms.



SPECIFICATIONS:

042-413 Deluxe L-Block Shield

Dimensions: 20" w x 17.88" depth x 28.75" h (50.8 x 45.4 x 73 cm)

Lead Shielding: 2.4" thick (6 cm)

Lead Glass Window:

Dimensions: 8" w x 8" h x 4" thick (20.3 x 20.3 x 10.2 cm)

Density: 5.2 g/cm³

Finish: Powder coat

Weight: 550 lb (250 kg)

Shipping Weight: 702 lb (318.4 kg)

042-417 Interlocking Lead Brick Cave

Dimensions:

I.D.: 15.3" w x 13.6" depth x 16" h (38.9 x 34.5 x 40.6 cm)

Lead Shielding: 2.4" thick (6 cm)

Finish: Paint

Weight: 611 lb (278 kg)

042-407 Steel Table

Dimensions: 36.75" w x 24" depth x 36" h (93.5 x 61 x 91.5 cm)

Shipping Weight: 195 lb (88 kg)

Detailed specifications on page 25

042-413 L-Block Shield, Deluxe, 2.4" lead

Related:

042-417 Lead Brick Cave, 3-wall, 2.4" lead
Fits 042-413 L-Block Shield

042-407 Table, Steel


Biodex L-Block Shields incorporate a hex-shaped plate to facilitate one-handed loading and unloading of Biodex PET Pigs.

COMPACT L-BLOCK WITH DOSE CALIBRATOR SHIELD

Space-saving design – Ideal for mobile units



Leaded frame prevents streaming around lead glass and provides shielding equal to the width of the L-Block Shield



- 1.5" thick lead shielding in front, 1" in base
- 8" x 8" x 4" adjustable lead glass window
- 1" thick lead shield surrounds calibration chamber
- Optional Lead Brick Cave for complete lateral shielding

The unique Compact L-Block with Dose Calibrator Shield is designed to maximize space in facilities receiving and preparing doses of high-energy nuclides such as FDG F-18. This unit provides convenient access and viewing of the work area and incorporates a built-in calibration chamber shield. The special shield is designed to accommodate a chamber that is through-mounted in a countertop (customer responsible for installation). The chamber shield accommodates all Atomlab chambers and many others (check chamber shield specifications to determine fit). This combination of L-Block and dose calibrator shield eliminates the need to purchase interlocking shielding rings. This unit is constructed of lead encased in steel. It features a large 8" x 8" x 4" lead glass window with adjustable window angle, 1.5" thickness lead shielding in front, and 1" thick lead in the base and in the chamber shield. A special plate with a hex-shaped recess is mounted on the L-Block base to facilitate one-handed loading and unloading of dose pigs incorporating hex-shaped bottoms. The optional 042-434 Lead Brick Cave fits neatly into the sides of the vertical section to provide lateral shielding around the full perimeter of the L-Block's base. For hot labs in mobile vans, the optional Brick Cave Cover will prevent the cave from shifting when the vehicle is in motion.

SPECIFICATIONS:

042-433 Compact L-Block with Dose Calibrator Shield

Dimensions: 18" w x 21.5" depth x 26" h (45.7 x 54.6 x 66 cm)

Lead Shielding:

Front: 1.5" thick (3.8 cm)

Base: 1" thick (2.5 cm)

Calibrator Shield: 1" thick (2.5 cm)

Calibrator Shield Inside Dimensions: 6.85" I.D. x 10.25" h (17.4 x 26 cm)

Lead Glass Window:

Dimensions: 8" w x 8" h x 4" thick (20.3 x 20.3 x 10.2 cm)

Density: 5.2 g/cm³

Finish: Powder coat

Weight: 570 lb (259 kg)

Shipping Weight: 590 lb (267.6 kg)

042-434 Interlocking Lead Brick Cave

Dimensions:

I.D.: 14" w x 20.5" depth x 16" h (35 x 52.1 x 40.6 cm)

Lead Shielding: 2" thick (5 cm)

Finish: Paint

Weight: 597 lb (271 kg)

042-433 Compact L-Block with
Dose Calibrator Shield, 1.5" lead
With built-in Dose Calibrator Shield

Related:

042-434 Lead Brick Cave, 3-wall, 2" lead
Fits 042-433 L-Block Shield

042-435 Lead Brick Cave Cover
Fits 042-434 Lead Brick Cave

Biodex L-Block Shields incorporate a hex-shaped plate to facilitate one-handed loading and unloading of Biodex PET Pigs.


L-BLOCK SHIELDS

For handling unit doses of high-energy radionuclides




042-428 L-Block Shield

Easy Assembly – no component weighs more than 50 pounds.



042-419 L-Block Shield



- **Lead shielding choices:**
1.5" – 2" thickness in front
1" thickness in base
- **Lead glass window choices:**
8" x 4" x 4"
8" x 8" x 4"
- **Installs easily**
- **Optional Lead Brick Cave for complete lateral shielding**

Designed for receiving and preparing unit doses of high-energy radionuclides, these L-Blocks provide a choice of 1.5" or 2" thick lead shielding in front, and 1" thick lead in the base. The L-Block may be ordered with either an 8" x 4" x 4" or 8" x 8" x 4" lead glass window. A special plate with a hex-shaped recess is mounted on the base to facilitate one-handed loading and unloading of dose pigs incorporating hex-shaped bottoms. An optional Lead Brick Cave may be added to provide lateral shielding around the full perimeter of the L-Block's base. These L-Blocks are shipped in modular form for easy installation without lifting equipment. No component weighs more than 50 pounds. After placing the base frame in its location, pre-cut lead sheets are loaded into the horizontal and vertical portions of the steel frame. The window module is mounted, and assembly is completed by securing window unit and steel end cap with four Phillips head screws. A Phillips screwdriver is the only tool needed. Assembly instructions are provided.

SPECIFICATIONS:

042-428 L-Block Shield (8" x 4" x 4" window)

Dimensions: 14" w x 15" d x 21.6" h (36 x 38 x 55 cm)

Lead Shielding:

Front: 1.5" (3.8 cm) thick

Base: 1" (2.5 cm) thick

Lead Glass Window:

Dimensions: 8" w x 4" h x 4" thick (20 x 10 x 10 cm)

Density: 5.2 g/cm³

Finish: Powder coat

Weight: 270 lb (122 kg)

Shipping Weight: 330 lb (149.6 kg)

042-419 L-Block Shield (8" x 8" x 4" window)

Dimensions: 14" w x 15" d x 24.7" h (36 x 38 x 62 cm)

Lead Shielding:

Front: 1.5" (3.8 cm) thick

Base: 1" (2.5 cm) thick

Lead Glass Window:

Dimensions: 8" w x 8" h x 4" thick (20 x 20 x 10 cm)

Density: 5.2 g/cm³

Finish: Powder coat

Weight: 290 lb (131 kg)

Shipping Weight: 355 lb (161 kg)

042-449 2" L-Block Shield

Dimensions: 14" w x 18" depth x 24.7" h (36 x 46 x 62 cm)

Lead Shielding:

Front: 2" thick (5.08 cm)

Base: 1" thick (2.5 cm)

Lead Glass Window:

Dimensions: 8" w x 8" h x 4" thick (20 x 20 x 10 cm)

Density: 5.2 g/cm³

Finish: Powder Coat

Weight: 345 lb (156 kg)

Shipping Weight: 420 lb (191 kg)

042-428 L-Block Shield, 1.5" lead

With 8" x 4" x 4" lead glass window

042-419 L-Block Shield, 1.5" lead

With 8" x 8" x 4" lead glass window

042-449 L-Block Shield, 2" lead

With 8" x 8" x 4" lead glass window

Related:

042-425 Lead Brick Cave, 3-wall, 2" lead

Fits 042-428 and 042-419 L-Block Shields

042-426 Lead Brick Cave, 3-wall, 2" lead

Fits 042-428, 042-449 and 042-419 L-Block Shields

Accommodates 042-466 PET Dose Drawing System

042-407 Table, Steel

Note: For detailed specifications on Interlocking

Lead Brick Caves see page 27.

TABLETOP SHIELDS

For Beta and Gamma Radiation




Protects face and body when working with radioactive material. Continuous stainless steel lining is easy to clean/decontaminate.

- Shields both Beta and Gamma Radiation
- Seamless Stainless Steel work area for easy cleaning and decontamination
- Work surface is designed with a spill proof lip for fast, contained clean up


This Tabletop Shield provides protection while working with mid to low-energy beta and gamma radioactive materials. The shield is designed to resist tipping and incorporate shielding overlap to eliminate radiation streaming. The entire interior work surface is one piece stainless steel for easy clean up and decontamination. The base of the work surface is designed with a spill proof lip for fast, contained clean up.

The front wall of the shield is constructed of .5" (1.2 cm) thick lead and completely encased with stainless steel on the work surface side and powder-coated cast aluminum on the technologist side. A lead glass window (2 mm LE) and a Plexiglas® panel are included. This shielding combination attenuates gamma and beta radiation, and errant bremsstrahlung, with unimpaired viewing of the work area. An optional 2 mm LE lead glass window is available to achieve 4.0 mm lead equivalency.

Optional side and bottom shields offer additional protection. Manufactured to be durable, the Beta/Gamma Tabletop Shield is framed in cast aluminum with a powder-coat finish.



Dimensions of Tabletop Shield



Tabletop Shield shown with optional side shield, 042-227

SPECIFICATIONS:

Tabletop Shields

Lead Shielding: .5" thick (1.2 cm)

Dimensions: See drawings featured

Viewing Panel:

Lead Equivalency: 2 mm (an optional lead glass panel can be paired to achieve 4.0 mm lead equivalency)

Clear Plexiglas®: For Beta Shielding

Viewing area: 11.1" x 11.1"

(28.2 x 28.2 cm)

Finish: Silver powder coat

Shipping Weight:

042-224: 79 lb (35.8 kg)

042-226: 10 lb (4.5 kg)

042-227: 42 lb (19 kg)

042-228: 10 lb (4.5 kg)

Side Shields

Lead Shielding: .25" thick (.64 cm)

Finish: Stainless Steel

Dimensions: See drawings featured

Shipping Weight: 42 lb (19 kg)

Tabletop Shields:

042-224 Shield, Tabletop, .5" lead

With 2.0 mm lead glass and Plexiglas panel

Optional:

042-228 Lead Glass Plate, 2.0 mm, 12" x 12"


042-227 Side Shield, .25" lead

042-226 Bottom Shield, .25" lead



Tabletop Shield shown with optional side shield.

RECTANGULAR LEAD BRICKS



Lead Bricks are a convenient solution to temporary or permanent shielding/storage situations. Easy to stack, unstack and relocate, lead bricks provide maximum protection wherever needed.

Produced from the highest quality lead, each brick is standard hardness with a flat, smooth surface, allowing a perfectly flush fit, even at sharp right angles.

SPECIFICATIONS:

002-246 Rectangular Lead Brick

Dimensions: 6" l x 4" w x 2" h (15 x 10 x 5 cm)

Weight: 21 lb (9.5 kg)

002-248 Rectangular Lead Brick


Dimensions: 8" l x 4" w x 2" h (20 x 10 x 5 cm)

Weight: 27 lb (12.5 kg)

002-246 Lead Brick, Rectangular, 6" long

002-248 Lead Brick, Rectangular, 8" long

GAMMA BENCH TOP SHIELD



The Bench Top Shield for gamma radiation is constructed of 1.375" thick (35 mm) lead acrylic, providing the equivalent of 0.060" (1.5 mm) lead for shielding gamma radiation. Lead acrylic can be used with low-energy gamma emitters. Lead acrylic is not recommended for beta emitters. This lightweight, sturdy shield is ideal for use in hoods, on laboratory benches or any "hot" area where local shielding is needed. Side legs give the shield exceptional stability.

SPECIFICATIONS:


Dimensions: 9" w x 12" h x 1.375" thick (230 x 305 x 35 mm)

Weight: 9 lb (4 kg)

Radionuclide	Approximate Attenuation Factor
I-125	10^{25}
Xe-133	10^{10}
Co-57	10^3
Tc-99m	75
I-123	6
Ga-67	2

042-710 Shield, Bench Top, Gamma

LEAD FOIL



Easily cut and shaped

Lead Foil has a myriad of uses around the laboratory. It is easily cut and formed into almost any shape for shielding areas or objects.

SPECIFICATIONS:

Dimensions: 6.7' l x 14" w x .03" thick (2 m x 36 cm x .76 mm)

Weight: 13.7 lb (6.2 kg)

104-030 Lead Foil

6.7' l x .030" thick (2 m x .76 mm)

LEAD VINYL SHEETS



Lead Vinyl Sheets are 10% lighter in weight than lead rubber with the same lead equivalent. The sheets feature uniform density throughout with an abrasion resistant, non-absorbing, smooth surface on both sides. Easily cut to any shape for protective applications, the Lead Vinyl Sheets are acid and alkali resistant, odorless and have an indefinite shelf life.

SPECIFICATIONS:

Dimensions: 2' l x 3' w (.61 x .92 m)

Nominal Thickness: .125" (0.3 cm)

Lead Equivalency: 1.00 mm

Color: Beige


Weight: 22 lb (10 kg)

055-903 Lead Vinyl Sheet



INTERLOCKING LEAD BRICKS

**DURABLE
PAINT
FINISH**

Interlocking design prevents leakage and assures stability



Interlocking Lead Bricks make it easy to erect, modify and relocate protective walls and caves of any size. Their V-shaped edges (A) eliminate the danger of leakage, common to all straight-edged bricks. Interlocking Lead Bricks create sturdy walls for a safe, protective enclosure.



This wall demonstrates the dimensions and alternative positions of all variations of the standard brick wall. The walls should be built with male chevrons directed upwards and to the right, as viewed from the outside of the cell.

Standard Bricks:		A x B x C
1	001-100 Brick, Standard, Wall <i>Weight: 13.2 lb (6 kg)</i>	2" x 4" x 4"
2	001-200 Brick, Standard, Base <i>Weight: 14.9 lb (6.9 kg)</i>	2" x 4" x 4"
3	001-300 Brick, Standard, Top <i>Weight: 5 lb (2.3 kg)</i>	2" x 4" x 2"

Corner Bricks:		A x B x C
4	001-400 Brick, Corner, Wall <i>Weight: 13.2 lb (6 kg)</i>	2" x 4" x 4"
5	001-500 Brick, Corner, Base <i>Weight: 14.9 lb (6.9 kg)</i>	2" x 4" x 4"
6	001-600 Brick, Corner, Top <i>Weight: 5 lb (2.3 kg)</i>	2" x 4" x 2"

Reverse Corner Bricks:		A x B x C
11	001-132 Brick, Reverse Corner, Wall <i>Weight: 13.2 lb (6 kg)</i>	2" x 4" x 4"
10	001-130 Brick, Reverse Corner, Base <i>Weight: 14.9 lb (6.8 kg)</i>	2" x 4" x 4"
12	001-135 Brick, Reverse Corner, Top <i>Weight: 5 lb (2.3 kg)</i>	2" x 4" x 2"


Left End Cap Bricks:		A x B x C
26	001-148 Brick, Left End Cap, Wall <i>Weight: 14.9 lb (6.9 kg)</i>	2" x 4" x 4"
16	001-147 Brick, Left End Cap, Base <i>Weight: 13.2 lb (6 kg)</i>	2" x 4" x 4"
14	001-142 Brick, Left End Cap, Top <i>Weight: 5.2 lb (2.4 kg)</i>	2" x 4" x 2"

Right End Cap Bricks:		A x B x C
27	001-149 Brick, Right End Cap, Wall <i>Weight: 14.9 lb (6.9 kg)</i>	2" x 4" x 4"
15	001-145 Brick, Right End Cap, Base <i>Weight: 13.2 lb (6 kg)</i>	2" x 4" x 4"
13	001-140 Brick, Right End Cap, Top <i>Weight: 5.2 lb (2.4 kg)</i>	2" x 4" x 2"

INTERLOCKING LEAD BRICK CAVES

**DURABLE
PAINT
FINISH**

Design your own enclosure or select from these most popular 3-walled caves.



042-434 Interlocking Lead Brick Cave,
2" lead (fits L-Block Shield 042-433)

The optional Lead Brick Caves fit neatly in the sides of the vertical section to provide lateral shielding around the perimeter of the L-Block's base. For hot labs in mobile vans, the optional Brick Cave Cover will prevent the cave from shifting when the vehicle is in motion.

SPECIFICATIONS:

042-434 Interlocking Lead Brick Cave

Dimensions:

I.D.: 14" w x 20.5" depth x 16" h (35 x 52.1 x 40.6 cm)

Lead Shielding: 2" thick (5 cm)

Finish: Paint

Weight: 597 lb (271 kg)

042-425 Interlocking Lead Brick Cave

Dimensions:

I.D.: 14" w x 15" depth x 16" h (35 x 38.1 x 40.6 cm)

Lead Shielding: 2" thick (6 cm)

Finish: Paint

Weight: 492 lb (223 kg)

042-417 Interlocking Lead Brick Cave

Dimensions:

I.D.: 15.3" w x 13.6" depth x 16" h (38.9 x 34.5 x 40.6 cm)

Lead Shielding: 2.4" thick (6 cm)

Finish: Paint

Weight: 611 lb (278 kg)

042-426 Interlocking Lead Brick Cave

Dimensions:

I.D.: 14" w x 17.8" depth x 13.8" h (35.5 x 45.3 x 34.6 cm)

Lead Shielding: 2" thick (5 cm)

Finish: Paint

Weight: 532 lb (241 kg)

042-434 Lead Brick Cave, 3-wall, 2" lead

Fits 042-433 L-Block Shield

042-425 Lead Brick Cave, 3-wall, 2" lead

Fits 042-428 and 042-419 L-Block Shields

042-417 Lead Brick Cave, 3-wall, 2.4" lead

Fits 042-413 L-Block Shield

042-426 Lead Brick Cave, 3-wall, 2" lead

*Fits 042-428, 042-449 and 042-419 L-Block Shields
Accommodates 042-466 PET Dose Drawing System*

Related:

042-435 Lead Brick Cave Cover

Fits 042-434 Lead Brick Cave

STEEL TABLE



- *Strong and sturdy*
- *Ideal for L-Block Shields and Lead Brick Caves*
- *Use in PET, nuclear medicine or radiation therapy departments*
- *Upper and lower shelf to accommodate small and large items, including PET shipping containers*

Sturdy as they come, this steel table can be used for just about any application requiring a strong, level platform. Ideal for holding heavy L-Block shields and caves, the surface is powder coated and the front legs feature adjustable levelers. Use the middle shelf to hold small items and the bottom shelf to support shipping containers or other large objects.

SPECIFICATIONS:

Dimensions: 36.75" w x 24" depth x 36" h (93.5 x 61 x 91.5 cm)

Front legs incorporate adjustable levelers

Finish: Powder coat

Weight Capacity: 1550 lb (730 kg)

Weight: 200 lb (90.9 kg)

Shipping Weight: 254 lb (116 kg)

042-407 Table, Steel

LEAD-LINED STORAGE SAFE



Conveniently loaded from the front, this Storage Safe is ideal for storing large quantities of high-activity radioisotopes. Shielded with a thickness of 2" of lead, the safe is encased in a powder-coated steel jacket and features an adjustable shelf. The lead-lined door is hung with heavy duty non-sagging hinges and is key-locked to prevent unauthorized access.

Transporting this half-ton safe is made easier with the built-in lifting handles for use with a hoist or other means.

SPECIFICATIONS:

Dimensions: 17.4" w x 17" depth x 19" h (44.2 x 43.2 x 48.3 cm)

I.D.: 12" w x 12" depth x 12" h (30.5 x 30.5 x 30.5 cm)

Lead Shielding: 2" thick (5 cm)

Finish: Powder coat

Door: Key-locked

Weight: 1050 lb (476 kg)

244-006 Safe, Storage, 2" lead

RADIOIODINE FUME HOOD

Ample work space for comfortable use



The Radioiodine Fume Hood meets the demands of iodination procedures. Constructed of 3/8" clear, rugged Plexiglas, the fume hood provides a large internal work area with spacious arm ports to allow maximum uninhibited manipulation of material within the unit. A swing-away front door permits easy placement and retrieval of items.

The air baffle assures even flow speed of air out of the box while negative air flow speed can be adjusted from 0 to a maximum of 50 CFM. The Fume Hood includes a disposable 12" x 12" x 1" metal frame filter that contains TEDA charcoal, trapping 90% of the radioiodine produced. The Fume Hood can accommodate two filters if needed.

SPECIFICATIONS:

Dimensions: 24" w x 21" depth x 36" h (61 x 53 x 91 cm)

Front Door (Swing-Away): 24" l x 13" h (61 x 33 cm)

Motor: 115 VAC, 50/60 Hz, 1 amp fuse; 230 VAC uses a separate external step down transformer

Weight: 90 lb (40.8 kg)

Shipping Weight: 114 lb (51.7 kg)

190-210 Radioiodine Fume Hood, 115V

Includes: Charcoal filter, TEDA

Related:

116-010 Gloves, Leaded Neoprene, pair

087-112 Converter, 230V

Replacement:

112-036 Filter, Charcoal (TEDA)

12" x 12" x 1" metal frame

SHIELDED WASTE CONTAINER

For low-energy gamma waste



The Shielded Waste Container is a must-have item for any facility that generates low-energy gamma radiation waste. Constructed entirely of 18-gauge stainless steel and lined with .125" (3 mm) lead, this 20-quart (18.9 liter) container can be placed on the floor or counter. Simply lift off the shielded cover for quick disposal of waste. Plastic liners make it easy to transfer waste to a decay or disposal site once the container is filled.

SPECIFICATIONS:

Dimensions: 11.9" w x 9.9" depth x 15.25" h (30.2 x 25.1 x 38.7)
Lead Shielding: .125" thick (3 mm)
Capacity: 20 qt (18.9 L)
Weight: 51 lb (22.6 kg)

039-106 Waste Container, Shielded, .125" lead
Includes 20 Poly Liners

Replacement:

040-108 Liner, Poly, 20/pkg

SHIELDED WASTE CONTAINER

For low-energy beta and gamma waste

KEY-LOCKED



This Shielded Waste Container is used in facilities that generate low-energy beta and gamma radiation waste.

The interior of the container is constructed of .063" aluminum (1.6 mm) and .25" lead (6 mm). The shielding combination attenuates gamma radiation, beta radiation and errant bremsstrahlung.

Extra protection is provided with a specially designed hatch door that protects the user from container contents even while open. Convenient side handles let you easily lift the container top to empty decayed contents. Sleek and sturdy, the exterior is constructed of steel, with a powder-coat finish.

SPECIFICATIONS:

Dimensions: 12" w x 22" h x 9" d (30 x 56 x 23 cm)
I.D.: 10.5" w x 14.5" h x 8.5" d (27 x 37 x 22 cm)
Lead Shielding: .25" thick (6 mm)
Aluminum Shielding: .063" (1.6 mm)
Finish: Powder coat
Shipping Weight: 120 lb (54 kg)

039-110 Waste Container, Shielded, .25" lead
Includes 20 Poly Liners

Replacement:

040-108 Liner, Poly, 20/pkg

MOBILE RADIATION SHIELD

Adjustable height, compact shield



The Compact Adjustable Height Radiation Shield is a versatile addition to any PET site. The 1" thick lead panel is 22" w x 19" h. The panel can be height adjusted within a 10" range. The shield can be used to protect technologists from sitting or lying "hot" patients. A fold down shelf is conveniently located to hold syringe shields, carriers and injectors, etc. Place shielding where and when it is needed. Wheel and swivel locks on each caster provide secure placement.

SPECIFICATIONS:

Dimensions: 24" w x 24.25" deep x 32.5" h (61 x 62 x 83 cm)
Shield: 22" w x 19" h (56 x 48 cm)
Folding Shelf: 24" w x 5" deep (61 x 12.7 cm)
Height Adjustable: 35" to 45" (89 to 114 cm)
Lead Shielding: 1" thick (2.5 cm)
Finish: Powder coat
Weight: 230 lb (104 kg)
Shipping Weight: 366 lb (166 kg)

042-522 Radiation Shield, Mobile,
Adjustable Height, Compact, 1" lead

MOBILE RADIATION SHIELD

Adjustable height shield



Rolling convenience with height adjustability

The Adjustable Height Mobile Radiation Shield puts shielding where it is needed. Roll into place, secure the wheels and adjust to the desired height. Panel can be height adjusted within a 4.5" range. Wheel and swivel locks on each caster provide secure placement.

SPECIFICATIONS:

Dimensions: 37.75" w x 29" depth x 40.5" h (95.9 x 73.7 x 102.9 cm)
Shield: 36" w x 19" h (91.4 x 48.3 cm)
Height Adjustable: 36" to 40.5" (91.4 to 103 cm)
Lead Shielding: 1" thick (2.5 cm)
Finish: Powder coat
Weight: 446 lb (202.3 kg)
Shipping Weight: 560 lb (254 kg)

042-519 Radiation Shield, Mobile,
Adjustable Height, 1" lead

CLEAR-LEAD™ MOBILE BARRIERS

Save Time – new sleek design is easy to clean!



042-580 Clear-Lead Mobile Barrier, Regular Window

- Meets all infection control criteria
- Sleek design is easy to clean
- Durable and shatter resistant protection
- Effortless maneuverability
- 0.5 mm lead equivalency window
- 0.8 mm lead opaque panel
- Superior stability

Introducing Clear-Lead™ – lead acrylic mobile x-ray barriers, featuring the first design change in thirty years. These barriers are for use around any imaging procedures using ionizing radiation, e.g., fluoroscopy, x-ray, CT, mammography.

The greatest advantage of these barriers over others on the market is the time they save behind the scenes. Their sleek design exposes no hardware, nooks or crannies while the opaque section of the shield is covered with plastic laminate – the result, a smooth, easy-to-clean finish.

Handle-grab edges and large casters make them effortless to maneuver. Clear-Lead™ Mobile Barriers offer large, durable and shatter resistant protection wherever it's needed.



042-582 Clear-Lead Mobile Barrier, Full Window

SPECIFICATIONS:

042-580 Clear-Lead Mobile Barrier, Regular Window

Dimensions: 75" h x 30" w (191 x 76 cm)
Window: 24" h x 30" w (61 x 76 cm)
Opaque Panel: 48" h x 30" w (122 x 76 cm)
Shielded Area: 72" h x 30" w (183 x 76 cm)

Shielding:

Window: 0.5 mm lead equivalency
Opaque Panel: 0.8 mm lead

Casters: Four hospital grade; two locking, two non-locking
Weight: 102 lb (46.5 kg)

042-582 Clear-Lead Mobile Barrier, Full Window

Dimensions: 75" h x 30" w (191 x 76 cm)
Window: 60" h x 30" w (152 x 76 cm)
Opaque Panel: 12" h x 30" w (31 x 76 cm)
Shielding Area: 72" h x 30" w (183 x 76 cm)

Shielding:

Window: 0.5 mm lead equivalency
Opaque Panel: 0.8 mm lead

Casters: Four hospital grade; two locking, two non-locking
Weight: 98 lb (44.5 kg)

042-580 Barrier, Clear-Lead™, Regular Window
Window size 24" h x 30" w (61 x 76 cm)

042-582 Barrier, Clear-Lead™, Full Window
Window size 60" h x 30" w (152 x 76 cm)

SHIELDED STORAGE CABINET

View contents without opening doors



The Lead-Lined Storage Cabinet is a must for any laboratory that handles radioactive materials. The lead shielding provides an extra measure of safety over the use of lead pigs alone. The .25" lead glass

sliding doors allow the viewing of materials without the extra exposure of having to open doors. Seeing before opening also decreases the time spent with the door open during material retrieval. Handsomely crafted with adjustable shelves, the white formica surface blends with any decor.

SPECIFICATIONS:

Dimensions: 15.5" w x 10" h x 7.5" depth (39 x 25 x 19 cm)
 I.D.: 14.125" w x 7.75" h x 5.75" depth (36 x 20 x 15 cm)
 Lead Shielding: .0625" thick (.15 cm)
 Lead Glass Window: .25" thick (.64 cm)
 Weight: 35 lb (16 kg)

154-090 Cabinet, Storage, .0625" lead

RECTANGULAR CONTAINER



The Rectangular lead and steel wall container is designed for storage and shipment of radioactive materials. The container has 1" lead shielding and offers adequate protection for most radioisotopes.


Hinged cover may be locked with padlock.

SPECIFICATIONS:

Dimensions:
 I.D.: 6" l x 6" w x 8" h (15 x 15 x 20 cm)
 O.D.: 8" l x 8" w x 11" h (20 x 20 x 28 cm)
 Lead Shielding: 1" thick (2.5 cm)
 Weight: 175 lb (80 kg)
 Shipping Weight: 180 lb (82 kg)

001-001 Container, Rectangular, 1" lead

NIPTONGS



These low-cost tongs are used to handle small radioactive or otherwise dangerous objects up to 1" diameter. The tongs have a 45° v cut groove on each jaw. The compression spring maintains a strong grip on the object until the tension is released by squeezing the finger bar. Niptongs are made of chrome-plated, high carbon steel with hardwood handles and are easily disassembled for decontamination and cleaning.

011-012 Niptongs, 12" (30 cm)

011-036 Niptongs, 36" (91 cm)

FORCEPS



Surgical forceps are excellent handling devices for small items in the hot lab. No need to physically touch vials, ampules, etc. Forceps are made of surgical grade steel.

066-533 Forceps, Curved, locking, 9.5" l (24.1 cm)

066-535 Forceps, Straight, locking, 9.5" l (24.1 cm)

066-536 Forceps, Curved, non-locking, 12.5" l (31.7 cm)

POWDER-FREE RADIATION ATTENUATING GLOVES

Reduced exposure, freedom of movement, natural latex free



Powder-Free Radiation Attenuating Gloves provide increased protection from direct x-ray beam and scattered radiation, and reduce harmful radiation exposure during any procedure requiring the use of fluoroscopy.

Just .007" thick, these gloves provide superior performance and dexterity over thicker, conventional radiation gloves. The reduced thickness allows for more flexibility and greater touch sensitivity while decreasing finger fatigue – all factors critical when working with delicate instruments. Applications include fluoroscopy, radioisotope handling, radiology and nuclear medicine.

Because these gloves are powder-free, lead-free and latex-free, they offer a reduced risk of natural rubber latex allergies and are an environmentally friendly alternative to leaded gloves.

Powder-Free Radiation Attenuating Gloves meet the following attenuation values:

58.7% attenuation at 60KVP HVL = 2.3 mm Al

49.9% attenuation at 80KVP HVL = 3.3 mm Al

44.6% attenuation at 100KVP HVL = 4.3 mm Al

40.6% attenuation at 120KVP HVL = 5.6 mm Al

All gloves are shipped sterile and ready to use in single pair packages.

208-065 Gloves, size 6-1/2, pr

208-070 Gloves, size 7, pr

208-075 Gloves, size 7-1/2, pr

208-080 Gloves, size 8, pr

208-085 Gloves, size 8-1/2, pr

208-090 Gloves, size 9, pr

LEAD GLASS GOGGLES



Comfortable protection

Designed to reduce radiation exposure to the eyes, these lead glass goggles are framed with soft, pliable vinyl and held on the head securely by an adjustable strap. Four direct ventilators prevent fogging. Eye shielding is provided by a 2" x 4.25" single sheet of fluoroscopic quality lead glass. The 4.2 density glass will effectively eliminate more than 95% of direct radiation produced by gamma rays. Glass provides 2.00 mm lead equivalency.

117-425 Goggles, Lead Glass

LEADED NEOPRENE GLOVES



These seamless, Lead-Loaded Neoprene Gloves are manufactured from an exclusive formulation. When handling materials emitting low-energy gammas they have proven highly effective and may reduce dosages up to 85%.

The high-density lead/neoprene shielding film provides the glove with good shielding properties while

still maintaining adequate flexibility & dexterity. The glove is easily cleaned with antiseptic, disinfectant and Radiacwash (featured separately) solutions. Ideal for I-125 and Tc-99m procedures. The gloves are provided in glove box size and can be cut or rolled down to the desired length.

SPECIFICATIONS:

Lead Equivalency: 0.10 mm
Nominal Thickness: 30 mils
Transmission: 50%
Half-Value Thickness: (0.32") 0.80 mm
Linear Absorption Coefficient: 0.87 mm

116-010 Gloves, Leaded Neoprene, pr

SOFT-SKIN COAT APRON



This comfortable, supple, protective apron is a pleasure to wear. The lead vinyl core is totally sealed in a tough, easy-to-clean, chemical, abrasion and aging resistant outer covering that is more pliable and lighter than standard lead rubber aprons.

The apron design differs from conventional aprons. There are no straps or buckles for support. Instead the apron is supported across the breadth of the shoulders, held snugly in any position the wearer assumes with the closures. The contour conforming principles assure the wearer of comfort and protection at all times.

SPECIFICATIONS:

Dimensions: 36" l x 24" w (91.4 x 61 cm)
Lead Equivalency: 0.5 mm
Color: Royal Blue (09)
Weight: 10 lb (4.6 kg)

103-701 Apron, Sof-Skin Coat

Lead-Lined Laboratory Furniture

Compact, attractive and offering benefits and options targeted at improving departmental efficiency and safety.



All Lead-Lined Cabinets:

- Are available with 0.25", 0.5" or 1" lead shielding
- Have key locking doors and drawers; card slots for labels
- Support up to 1500 lb
- Can be moved with a standard pallet jack
- Have built-in seismic anchor brackets
- Have lead shielding fully encased in steel - no exposed lead
- Have durable powder-coat finish and stainless steel tops
- Measure 36.5" high, 24" deep (except Phantom Cabinet)
- Can stand alone or in combination with other lead-lined or conventional cabinets (except Phantom Cabinet)

Biodex Lead-Lined Laboratory Furniture is compact, attractive and offers benefits and options targeted at improving department efficiency and safety. Choose the appropriate shielding thickness to satisfy your particular radiation safety needs. Mix and match cabinets to improve workflow and procedural efficiency. These 24" depth cabinets are sized to industry standards for neat configuration with existing furniture. Stainless tops with backsplash and lips help control spills and reduce cleanup and decontamination time. The durable powder-coat finish resists scratching, ensuring a clean, aesthetic appearance for years to come. Install any of these units in the smallest hot lab without special lifting or moving equipment - a simple pallet jack will do the job.

VERSATILE - All cabinets can be used individually, or in any combination to meet your storage, decay, and workbench requirements. Unified stainless steel countertops with 4" backsplash and 0.5" lip are available for multiple cabinet configurations. Tops can be ordered 6" deeper than cabinets to accommodate pipes and services.

SAFE - Completely shielded with lead on all six sides, most units are available with shielding thickness of 0.25", 0.5" and 1".

MODULAR - As needs change, so may department layouts. Biodex Lead-Lined Laboratory Furniture can stand alone or be combined with additional units at any time - without costly room alterations.

SECURE - All doors have key locks for security. All cabinets include a set of brackets for seismic anchoring.

BIODEX

www.biodex.com
1-800-224-6339


Int'l 631-924-9000

www.biodex.com/furniture

LEAD-LINED RADIOISOTOPE STORAGE CABINET



Designed for safely storing radioactive materials.



Cabinet drawing dimensions can be accessed at www.biodex.com

The Radioisotope Storage Cabinet is designed for safely storing radioactive materials. It features 12 key-locked drawers. Each drawer is easily removed for cleaning or decontamination. A card slot identifies contents.

SPECIFICATIONS:

Dimensions: 30.5" w x 24" depth x 36.5" h (77.5 x 61 x 92.7 cm)

Lead Shielding: .25" (.64 cm), .5" (1.3 cm) or 1" (2.5 cm) thick on all six sides

Drawer Dimensions: I.D.: 6" w x 11.88" depth x 4.5" h (15.2 x 30.2 x 11.4 cm)

Drawers: Key-locked

Countertop: Stainless steel with 4" (10.2 cm) backsplash and .5" (1.3 cm) spillproof lip

Finish: Powder coat

Weight:

244-110: 1069 lb (486 kg)

244-111: 1415 lb (642 kg)

244-112: 2558 lb (1160 kg)

244-110 Cabinet, Radioisotope Storage, .25" lead


244-111 Cabinet, Radioisotope Storage, .5" lead

244-112 Cabinet, Radioisotope Storage, 1" lead

LEAD-LINED DECAY CABINET



Can be used to decay sharps containers prior to disposal



Cabinet drawing dimensions can be accessed at www.biodex.com

The Decay Cabinet is designed for long and short term storage of decaying radioactive material. Two adjustable shelves support up to 100 lb each. The door is key-locked to prevent unauthorized access. The cabinet will accommodate sharps containers and other boxed waste prior to disposal. It can also be used to store flood sources.

SPECIFICATIONS:

Dimensions: 30.5" w x 24" depth x 36.5" h (77.5 x 61 x 92.7 cm)

Lead Shielding: .25" (.64 cm), .5" (1.3 cm) or 1" (2.5 cm) thick on all six sides

Shelf Dimensions: 26" w x 18" depth (66 x 45.7 cm), 100 lb (45.4 kg) capacity, adjustable height

Door: Key-locked

Countertop: Stainless steel with 4" (10.2 cm) backsplash and .5" (1.3 cm) spillproof lip

Finish: Powder coat

Weight:

244-140: 1010 lb (458 kg)

244-141: 1267 lb (575 kg)

244-142: 2125 lb (964 kg)

244-140 Cabinet, Decay, .25" lead

244-141 Cabinet, Decay, .5" lead

244-142 Cabinet, Decay, 1" lead

Note: Reverse door swing available.


Call for quote.

LEAD-LINED PREPARATION ENCLOSURE BASE CABINET



Built to support the Preparation Enclosure

This Cabinet is designed to support the Lead-Lined Preparation Enclosure. Full height, overlapping double doors with key locks open to an adjustable shelf with a 100 lb capacity. The cabinet may be used for decay and storage.



Cabinet drawing dimensions can be accessed at www.biodex.com

SPECIFICATIONS:

Dimensions: 36.5" w x 24" depth x 36.5" h (92.7 x 61 x 92.7 cm)

Lead Shielding: .25" (.64 cm), .5" (1.3 cm) or 1" (2.5 cm) thick on all six sides

Shelf Dimensions: 30" w x 18" depth (76.2 x 45.72 cm),

100 lb (45.4 kg) capacity, adjustable height

Doors: Key-locked

Countertop: Stainless steel with 4" (10.2 cm) backsplash and .5" (1.3 cm) spillproof lip

Finish: Powder coat

Weight:

244-190: 1063 lb (483.2 kg)

244-191: 1540 lb (699 kg)

244-192: 2433 lb (1104 kg)


244-190 Cabinet, Preparation Enclosure Base, .25" lead

244-191 Cabinet, Preparation Enclosure Base, .5" lead

244-192 Cabinet, Preparation Enclosure Base, 1" lead

LEAD-LINED PREPARATION ENCLOSURE


Connects to external ductwork



Adjustable 12" shield can fold forward to load large objects. The hood can also be completely closed and used for storage. Preparation Enclosure features built-in electrical outlets.

The Lead-Lined Preparation Enclosure is designed for applications that require handling gaseous radioactive materials. The interior provides ample floor space. A large lead glass window and halogen light allow safe and unobstructed viewing. The enclosure opening has an adjustable shield that creates access ports. A swing down shield covers the ports when not in use.

Gaseous materials are directed by a baffle (eliminating dead space) through a stainless steel chimney. The chimney is connected to external blowers (not supplied) and ductwork (not supplied). The blowers create a negative pressure preventing gas leaks.



Cabinet drawing dimensions can be accessed at www.biodex.com

SPECIFICATIONS:

Dimensions: 36" w x 24" depth x 30.5" h (91.4 x 61 x 77.5 cm)

Interior Floor Space: 31" w x 19" depth (78.7 x 48.3 cm)

Lead Shielding: .25" thick (.64 cm)

Adjustable Shield: 12" w x 10" h x .5" thick (30.5 x 25.4 x 1.3 cm)

Exhaust: 6" dia (15.2 cm) chimney, fixed upper and adjustable lower baffles. Blower and filter not included

Lighting: 110V Halogen lamp with two 25 watt bulbs, UL listed; 220V available upon request.

Lead Glass Window:

Dimensions: 34.5" w x 11.8" h x .75" thick (87.6 x 30 x 1.9 cm)

Density: 5.05 g/cm³

Finish: # 3 brushed, stainless steel

Weight: 733 lb (332 kg)


244-007 Preparation Enclosure, .25" lead

LEAD-LINED DECAY AND STORAGE CABINET



The Decay and Storage Cabinet performs two functions. Radioisotopes can be safely stored in drawers, while the cupboard section is for the storage of decaying material.

Drawers are easily removed for cleaning and decontamination. The decay section includes two heavy duty adjustable shelves.



Cabinet drawing dimensions can be accessed at www.biodex.com

SPECIFICATIONS:

Dimensions: 36.5" w x 24" depth x 36.5" h (92.7 x 61 x 92.7 cm)

Lead Shielding: .25" (.64 cm), .5" (1.3 cm) or 1" (2.5 cm) thick on all six sides

Drawer Dimensions: I.D.: 6" w x 11.88" depth x 4.5" h (15.2 x 30.2 x 11.4 cm)

Shelf Dimensions: 22" w x 18" depth (55.9 x 45.7 cm)

Door and Drawers: Key-locked

Countertop: Stainless steel with 4" (10.2 cm) backsplash and .5" (1.3 cm) spillproof lip

Finish: Powder coat

Weight:

244-160: 1103 lb (501.4 kg)

244-161: 1421 lb (645 kg)

244-162: 2500 lb (1134 kg)

244-160 Cabinet, Decay and Storage, .25" lead


244-161 Cabinet, Decay and Storage, .5" lead

244-162 Cabinet, Decay and Storage, 1" lead

LEAD-LINED WASTE CABINET



The Waste Cabinet is designed for storing "non-sharps" radioactive waste. A spacious chute with shielded cover allows waste to be dropped directly into a polyethylene container prior to decay and disposal.



Cabinet drawing dimensions can be accessed at www.biodex.com

SPECIFICATIONS:

Dimensions: 30.5" w x 24" depth x 36.5" h (77.5 x 61 x 92.7 cm)

Lead Shielding: .25" (.64 cm), .5" (1.3 cm) or 1" (2.5 cm) thick on all six sides

Chute: 6.5" dia (16.5 cm) with .5" thick (1.3 cm) lead shielded cover

Container: 22" h x 17" dia (56 x 43 cm) polyethylene, 12 gal capacity

Door: Key-locked

Countertop: Stainless steel with 4" (10.2 cm) backsplash and .5" (1.3 cm) spillproof lip

Finish: Powder coat

Weight:

244-150: 913 lb (415 kg)

244-151: 1282 lb (583 kg)

244-152: 2290 lb (1039 kg)

244-150 Cabinet, Waste, .25" lead

244-151 Cabinet, Waste, .5" lead

244-152 Cabinet, Waste, 1" lead

Note: Reverse door swing available.


Call for quote.

LEAD-LINED WASTE AND STORAGE CABINET



The Waste and Storage Cabinet performs two functions. Radioisotopes can safely be stored in drawers, while the cupboard section is for the storage of "non-sharps" radioactive waste.

Drawers are easily removed for decontamination. The cupboard section includes a shielded port and a 12 gallon polyethylene container.



Cabinet drawing dimensions can be accessed at www.biodex.com

SPECIFICATIONS:

Dimensions: 36.5" w x 24" depth x 36.5" h (92.7 x 61 x 92.7 cm)

Lead Shielding: .25" (.64 cm), .5" (1.3 cm) or 1" (2.5 cm) thick on all six sides

Drawer Dimensions: I.D.: 6" w x 11.88" depth x 4.5" h (15.2 x 30.2 x 11.4 cm)

Waste Section:

Chute: 6.5" dia (16.5 cm) with .5" thick (1.3 cm) lead shielded cover

Container: I.D.: 22" h x 17" dia (56 x 43 cm) polyethylene, 12 gal capacity

Door and Drawers: Key-locked

Countertop: Stainless steel with 4" (10.2 cm) backsplash and .5" (1.3 cm) spillproof lip

Finish: Powder coat

Weight:

244-170: 863 lb (391 kg)

244-171: 1439 lb (653 kg)

244-172: 2462 lb (1117 kg)

244-170 Cabinet, Waste and Storage, .25" lead

244-171 Cabinet, Waste and Storage, .5" lead

244-172 Cabinet, Waste and Storage, 1" lead

LEAD-LINED SINK AND WASTE CABINET



Two sections for hot and cold waste

The Sink and Waste Cabinet performs three functions. A stainless steel sink allows the convenience of running water in the hot lab. The space under the sink is used for cold storage. Separated from the sink section by a lead barrier, the waste section includes a shielded port that allows waste to be dropped into a polyethylene container for storage until decayed.



Cabinet drawing dimensions can be accessed at www.biodex.com

SPECIFICATIONS:

Dimensions: 30.5" w x 24" depth x 36.5" h (77.5 x 61 x 92.7 cm)

Lead Shielding: .25" (.64 cm), .5" (1.3 cm) or 1" (2.5 cm) thick on all six sides

Sink Section: 10" w x 14" depth x 10" h (25.4 x 35.6 x 25.4 cm)

Integral stainless steel sink with gooseneck faucet and wrist blades; suitable for cold storage

Waste Section:

Chute: 6.5" dia (16.5 cm) with .5" thick (1.3 cm) lead shielded cover

Container: I.D.: 20.5" h x 11.25" dia (52 x 28.6 cm), polyethylene, 7 gal capacity

Door: Key-locked

Countertop: Stainless steel with 4" (10.2 cm) backsplash and .5" (1.3 cm) spillproof lip

Finish: Powder coat

Weight:

244-130: 849 lb (386 kg)

244-131: 1253 lb (569.5 kg)

244-132: 2260 lb (1027.3 kg)

244-130 Cabinet, Sink and Waste, .25" lead

244-131 Cabinet, Sink and Waste, .5" lead

244-132 Cabinet, Sink and Waste, 1" lead

Notes: Reverse door swing available.

Foot Pedal operated sink is available.


Call for quote.

LEAD-LINED GENERATOR AND STORAGE CABINET



Larger generator drawer will accommodate the Mallinckrodt Generator with Elution Shield.

Trap doors on the top and sides of this drawer provide safe access to the generator. The cabinet has two shielded compartments on the right side for storing decaying generators prior to disposal.



Cabinet drawing dimensions can be accessed at www.biodex.com.

SPECIFICATIONS:

Dimensions: 36.5" w x 24" depth x 36.5" h (92.7 x 61 x 92.7 cm)

Lead Shielding:

Cabinet: .5" (1.3 cm) or 1" (2.5 cm) thick on all six sides

Drawer: .25" thick (.64 cm) on five sides

Dividers: .5" thick (1.3 cm)

Drawer Dimensions: I.D.: 15.5" w x 11.9" depth x 18.2" h (39.4 x 30.2 x 46.2 cm)

Doors and Drawer: Key-locked

Countertop: Stainless steel with 4" (10.2 cm) backsplash and .5" (1.3 cm) spillproof lip

Finish: Powder coat

Weight:

244-181: 1917 lb (869.5 kg)

244-182: 2815 lb (1277 kg)

244-181 Cabinet, Generator and Storage, .5" lead

244-182 Cabinet, Generator and Storage, 1" lead

LEAD-LINED UNIT DOSE CABINET




The Unit Dose Cabinet is designed for hot labs with limited space.

The upper left compartment holds two sharps containers to facilitate decay rotation. A small chute with removable shielded lid allows syringe disposal into the front container. A large rear port allows safe removal of the second decayed container. The upper right compartment is for bulk storage of unit dose ammo boxes. The shielded door provides access without additional exposure.

Side-by-side middle drawers can be used for storage of syringes, sources, pigs, radioisotopes and other small items requiring lead shielding.

Phantoms and flood sources can be stored in the bottom drawer which runs the full width of the cabinet.



Cabinet drawing dimensions can be accessed at www.biodex.com.

SPECIFICATIONS:

Dimensions: 30.5" w x 24" depth x 36.5" h (77.5 x 61 x 92.7 cm)

Lead Shielding: .25" (.64 cm), .5" (1.3 cm) or 1" (2.5 cm) thick on all six sides

Drawer Dimensions:

Bottom Drawer: I.D.: 23.5" w x 21.8" depth x 4.5" h (59.7 x 55.4 x 11.4 cm)

Left & Right Drawers: I.D.: 9.38" w x 17.88" depth x 4.5" h (23.8 x 45.4 x 11.4 cm)

Upper Doors: I.D.: 11.5" w x 17.62" depth x 11" h (29.2 x 44.7 x 27.9 cm)

Doors and Drawers: Key-locked

Countertop: Stainless steel with 4" (10.2 cm) backsplash and .5" (1.3 cm) spillproof lip

Finish: Powder coat

Weight:

244-120: 1156 lb (524.3 kg)

244-121: 1614 lb (732 kg)

244-122: 2328 lb (1056 kg)

244-120 Cabinet, Unit Dose, .25" lead

244-121 Cabinet, Unit Dose, .5" lead

244-122 Cabinet, Unit Dose, 1" lead

Related:

039-338 Monoject Container, 4 qt., 10/pkg

Fits 039-325, 039-326, 039-350 and 039-330


039-387 Monoject Container, 8 qt., 10/pkg

Fits 039-325, 039-326 and 039-330

LEAD-LINED PHANTOM CABINET



Shown with Decay Module



Cabinet drawing dimensions can be accessed at www.biodex.com

The Phantom Cabinet stores flood sources and phantoms* on their ends. The interior is sectioned by an adjustable divider. The Lead-Lined Phantom Cabinet cannot be ordered separately, it must be part of a multiple cabinet configuration.

**Note: Phantom Cabinet does not accommodate 24" circular phantoms.*

SPECIFICATIONS:

Dimensions: 11.6" w x 24" depth x 36" h (29.5 x 60.9 x 91.4 cm)
Lead Shielding: .25" thick (.64 cm)
Door: Key-locked
Finish: Powder coat
Weight: 319 lb (145 kg)

244-009 Cabinet, Phantom, .25" lead

Note: Must be ordered as part of a multiple cabinet configuration

LEAD-LINED REFRIGERATOR



Lead-Lined Refrigerator is an under-the-counter unit



Cabinet drawing dimensions can be accessed at www.biodex.com

The Lead-Lined Refrigerator has 4.4 cubic feet of storage space and is ideal for storage and refrigeration of radiopharmaceuticals, tagged biological and other radioactive materials.

SPECIFICATIONS:

Dimensions: 22.2" w x 20.8" depth x 34.2" h (56.4 x 52.9 x 86.9 cm)
Lead-Lined Refrigerator is an under-the-counter unit
Lead Shielding: .125" thick (.32 cm)
Capacity: 4.4 cu ft
Door: Key-locked
Power: 115 V, 60 Hz, 1.3 amp
Finish: Powder coat
Weight: 315 lb (143 kg)

244-003 Refrigerator, Lead-Lined


Related:

087-112 Converter, 230V

LEAD-LINED PET UNIT DOSE CABINET



Available with .25", .5" and 1" lead shielding.



Cabinet drawing dimensions can be accessed at www.biodex.com.

- Designed for PET hot labs with limited space
- Accommodates:
 - Compact L-Block Shield with Built-in Dose Calibrator Shield
 - PET Sharps Container Shield
 - Lead Brick Cave
 - Atomlab Dose Calibrators and most others
- Sliding shelves for:
 - PET shipping containers
 - Small items
- Lead shielded on all six sides
- Key-locked doors

Designed for PET hot labs with limited space, the PET Unit Dose Cabinet provides a space-efficient work area over a fully shielded storage cabinet.

The cabinet supports the 042-433 Biodex Compact L-Block with Built-in Dose Calibrator Shield, the 039-412 Sharps Container Shield, the 042-434 Lead Brick Cave, and accommodates all of our Atomlab Dose Calibrators and many others. The dose calibrator display unit mounts on a stand above the countertop to maximize work space.

The lower cabinet has key-locking doors, two sliding bottom shelves, and two sliding upper shelves. The bottom shelves will accommodate PET shipping containers. The top shelves conveniently store syringes, syringe shields, and other small items. This cabinet is completely shielded on all six sides with .25", .5" and 1" lead, and can stand alone or be grouped with other cabinets.

All cabinets in this product line are built to the industry standard height of 36.5". All units include a stainless steel countertop incorporating a 0.5" lip and 4" backsplash. When ordering multiple units for grouped configuration, a unified countertop may be ordered to provide a continuous work surface.

Upon request, Biodex will factory mount the L-Block and Sharps Shield in place eliminating the need for on-site lifting equipment. A simple pallet jack is all that is required to move the unit.

SPECIFICATIONS:

Dimensions: 36.5" w x 24" depth x 36.5" h (93 x 61 x 93 cm)

Lead Shielding: .25" thick (.64 cm)

Finish: Powder coat

Doors: Key-locked

Countertop: Stainless steel with 4" (10.2 cm) backplash and .5" (1.3 cm) spillproof lip

Weight Capacity: 1550 lb (703 kg)

Weight: 1240 lb (562 kg)

244-200 Cabinet, PET, Unit Dose, .25" lead

Does not accommodate Lead Brick Cave

244-205 Cabinet, PET, Unit Dose, .25" lead

Accommodates Lead Brick Cave 042-434

Note: The cabinet is also available with .5" and 1" lead shielding. Call for quote.

Related:

042-433 L-Block Shield, Compact, 1.5" lead

With built-in Dose Calibrator Shield

042-434 Lead Brick Cave, 3-wall, 2" lead

Fits 042-433 L-Block Shield

086-332 Dose Calibrator, Atomlab™ 500*Plus*, 100-240 VAC

Includes: Smart Display, ionization chamber, well counter, RS-232 port, vial/syringe dipper and well insert.

039-412 Sharps Container Shield, PET, 1" lead

Uses one 039-413 Sharps Container

039-413 Sharps Container, 3.2 qt., 30/pkg

Fits 039-412

PET UNIT DOSE TABLE

More Features, More Shielding



- *Designed for PET hot labs with limited space*
- *Lead shielding under calibration chamber*
- *Shipped Completely Assembled:*
 - *Compact L-Block Shield with Built-in Dose Calibrator Shield*
 - *PET Sharps Container Shield*
- *Shelf space for:*
 - *PET shipping containers*
 - *Small items*

Designed for PET hot labs with limited space, the PET Unit Dose Table includes added features to help maximize the work area and the protection. The table includes the 042-433 Biodex Compact L-Block with Built-in Dose Calibrator Shield and the 039-412 Sharps Container Shield.

The Compact L-Block with Dose Calibrator Shield features a large 8" x 8" x 4" lead glass window with adjustable window angle, 1.5" thickness lead shielding in front, and 1" thick lead in the base and in the chamber shield. It is designed to accommodate a chamber that is through mounted in the countertop. The chamber shield accommodates all Atomlab chambers and many others.

The Sharps Container Shield is constructed of steel and lined with 1" thick (2.5 cm) lead. It features a lockable sliding cover for container removal and a hinged top door for syringe disposal.

Shipped completely assembled, the Unit Dose Table eliminates the need for on-site lifting equipment. A simple pallet jack is all that is required to move the unit. Once the Dose Calibrator Chamber is placed in the shield, the chamber support shelf provides 1" thick lead shielding. The Dose Calibrator display unit mounts on a stand above the counter-top to further maximize work space.

A spacious bottom shelf can hold several PET shipping containers. The middle shelf provides convenient storage for syringe shields and other small items.

In applications where additional shielding is required, the optional 042-434 Lead Brick Cave is available.

SPECIFICATIONS:

042-448 PET Unit Dose Table:

Dimensions: 36.75" w x 24" depth x 36" h (93.5 x 61 x 91.5 cm)

Front legs incorporate adjustable levelers.

Lead Shielding: 1" thick (enclosed in chamber mounting shelf)

Finish: Powder coat

Weight Capacity: 1550 lb (703 kg)

Weight: 294 lb (133.3 kg)

042-433 Compact L-Block with Dose Calibrator Shield:

Dimensions: 18" w x 21.5" depth x 26" h (45.7 x 54.6 x 66 cm)

Lead Shielding:

Front: 1.5" thick (3.8 cm)

Base: 1" thick (2.5 cm)

Calibrator Shield: 1" thick (2.5 cm)

Calibrator Shield Inside Dimensions: 6.85" I. D. x 10.25" h (17.4 x 26 cm)

Lead Glass Window:

Dimensions: 8" w x 8" h x 4" thick (20.3 x 20.3 x 10.2 cm)

Density: 5.2 g/cm³

Finish: Powder coat

Weight: 590 lb (259 kg)

Detailed Specifications on page 20

039-412 High-Energy PET Sharps Container Shield:

Dimensions: 12" h x 8.75" dia (30.5 x 22.2 cm)

Lead Shielding:

Sides and Bottom: 1" thick (2.5 cm)

Rotating Cover: .875" thick (2.2 cm)

Hinged Door: .625" thick (1.5 cm)

Security: Key-locked

Weight: 175 lb (79.4 kg)

Shipping Weight: 222 lb (100.6 kg)

Detailed Specifications on page 18

042-448 Table, PET, Unit Dose

Includes: 042-433 Compact L-Block with Dose Calibrator Shield and 039-412 Sharps Container Shield

Related:

042-434 Lead Brick Cave, 3-wall, 2" lead

Fits 042-433 L-Block Shield

042-435 Lead Brick Cave Cover

Fits 042-434 Lead Brick Cave

086-332 Dose Calibrator, Atomlab™ 500Plus, 100-240 VAC

Includes: Smart Display, ionization chamber, well counter, RS-232 port, vial/syringe dipper and well insert.

039-413 Sharps Container, 3.2 qt., 30/pkg

Fits 039-412

DOSE DRAWING SYRINGE SHIELDS

For drawing a dose from a shielded vial



Dose Drawing Syringe Shield (glass non-replaceable)



Dose Drawing Syringe Shield with replaceable glass

Dose Drawing Syringe Shields reduce hand exposure when drawing doses from vial shields. The barrel of the syringe is constructed of optically clear high density (5.6) lead glass that offers complete 360° visibility and reduces radiation exposure from Tc-99m by more than 99%. At the end of the lead glass barrel is a .125" lead flange that provides additional hand shielding between a vial shield and the syringe shield.

The shield features quick and smooth syringe insertion with an O-ring seal and anti-roll cap. Dose Drawing Syringe Shields for safety syringes facilitate two thumbscrews to secure the syringe.

The shields are constructed with an elastomer material that acts as a shock absorber between the metal and glass, improving the resistance to breakage.

Dose Drawing Syringe Shields:

Accommodates luer and non-luer lock syringes

007-661 Syringe Shield, 3 cc

007-663 Syringe Shield, 5 cc and 6 cc

007-665 Syringe Shield, 10 cc and 12 cc

Dose Drawing Syringe Shields with Replaceable Glass:

Accommodates luer and non-luer lock syringes

007-691 Syringe Shield, 3 cc

007-693 Syringe Shield, 5 cc and 6 cc

007-695 Syringe Shield, 10 cc and 12 cc

Replacement Glass for Dose Drawing Syringe Shields:

Accommodates luer and non-luer lock syringes

127-691 Syringe Shield Replacement Glass, 3 cc

127-693 Syringe Shield Replacement Glass,
5 cc and 6 cc

127-695 Syringe Shield Replacement Glass,
10 cc and 12 cc

HIGH DENSITY LEAD GLASS VIAL SHIELD



- *Lead glass provides clear visibility*
- *Accommodates most vials*
- *Centering action holds vials securely*

The High Density (5.6) Lead Glass Vial Shield reduces hand exposure and offers complete 360° visibility. The lead glass vial shield is suitable for low-energy radioisotopes.

The shield has a removable cap that makes cleaning and needle insertion as simple as possible while maintaining a sleek attractive appearance.

Automatic centering action positions vials within the shield and holds them securely for extra safety and convenience.

SPECIFICATIONS:

Lead Equivalency: .12" (3 mm)

HVL for 99m-Tc: 10

Accommodates Vial Sizes: 5 thru 30 ml

Weight: 3 lb (1.4 kg)

001-075 Vial Shield, .12" lead equiv

ZEVALIN™ Y-90 REACTION VIAL SHIELD



Designed for the reaction vial supplied in the Y-90 Zevalin™ kit

- Lead and aluminum construction attenuates beta and bremsstrahlung radiation
- Inside height and diameter sized specifically for Y-90 reaction vial
- Pivoting aperture cover for quick and convenient access to vial septum
- Removable top and bottom to minimize exposure during vial transfer

Y-90 Zevalin mixing and dose drawing is performed using a 10 ml reaction vial supplied by IDEC in the Zevalin kit. This shield is designed to fit that specific vial, and to attenuate beta and bremsstrahlung radiation associated with yttrium-90. The shield's top and bottom are threaded and knurled to allow quick insertion and removal of vials with minimal hand exposure. The top cover incorporates a septum access aperture with a pivoting shield to further reduce exposure.

SPECIFICATIONS:

Dimensions: 2.2" dia x 3.2" h (5.6 x 8.1 cm)

Inside Dimensions: 1.06" dia x 2" h (4 x 5.1 cm)

Shielding:

Cylinder: 0.5" thick (1.3 cm) lead and 0.06" thick (1.5 mm) aluminum

Bottom End Cap: 0.5" thick (1.3 cm) lead and 0.06" thick (1.5 mm) aluminum

Top End Cap: 0.25" thick (0.64 cm) lead

Cover: 0.31" thick (0.79 cm) lead, 0.06" thick (1.5 mm) aluminum

Weight: 3.9 lb (1.8 kg)

053-607 Vial Shield, Y-90, Reaction, 10 ml

TUNGSTEN VIAL SHIELD



- Virtually unbreakable
- Safe handling of radioactive liquids

The Tungsten Vial Shield is designed to greatly reduce exposure to vials containing liquid radioisotopes. The shield is constructed of .19" (.48 cm) thick tungsten, equivalent to .38" (1 cm) lead at 150 keV.

Vials can be loaded from the top or bottom of the shield. A loss-proof slide injection port on the top allows easy access to the vial septum.

The tough tungsten will retain its shape under the roughest handling and is virtually unbreakable.

SPECIFICATIONS:

Dimensions: 1.44" dia x 2.8" h (3.7 x 7 cm)

I.D.: 1" x 2.25" (2.5 x 5.7 cm)

Lead Equivalency: .38" (1 cm)

Weight: 1.65 lb (0.75 kg)

053-805 Vial Shield, Tungsten

VIAL SHIELD



For heated preparation of Sulphur Colloid, MAG-3, Cardiolite, AcuTect and NeoTect

This Vial Shield is designed to aid in preparation of radiopharmaceuticals that require boiling.

Vents are located to minimize scatter leakage, boiling water can circulate freely around the vial, heating the solution rapidly and uniformly. The carrying handle makes it easy to lower and remove the vial from the boiling water bath.

The vial is constructed of lead .25" (.64 cm) thick. A 5.6 density flush mounted lead glass window provides protection and visibility. The radiation level for 25 mCi of Tc-99m is reduced to background.

SPECIFICATIONS:

Dimensions: 2" dia x 3.875" h (5 x 9.8 cm)

Lead Shielding: .25" thick (.64 cm)

Accommodates Vial Sizes: up to 1.5" dia x 3.125" h (3.8 x 7.9 cm)

Weight (including handle): 3 lb (1.4 kg)

001-236 Vial Shield, .25" lead

VIAL PIG



Compatible with the PET Shipping System

Vials containing PET or other high-energy radionuclides can be safely transported in this 10/30ml Vial Pig. The pig offers a minimum of 1" of lead shielding to accommodate the concentrated energy.

The Vial Pig can be used independently or with the Biodex PET Shipping System, which meets DOT II Type A packaging requirements.

SPECIFICATIONS:

Dimensions:

Exterior: 6.63" h x 4.15" dia (16.8 x 10.5 cm)

Interior: 2.76" h x 1.51" dia (7 x 4.4 cm)

Lead Shielding:

Sides and Bottom: 1" thick (2.5 cm)

Top: 1.75" (3.8 cm)

Weight: 21.3 lb (9.7 kg)

001-706 Pig, Vial, PET, 10/30 ml, 1" lead
For 30 ml vials

Includes: 001-707 Vial Pig Adapter, to accommodate 10 ml vials and three Absorbent Sheets

Related:

001-771 Sheets, Absorbent, 100/pkg

COLOR-CODED SYRINGE LABELS



Labels quickly identify the radiopharmaceutical


MDP 212-001 Purple	TI 212-009 Reflex Blue	CARDIOLITE 212-019 Pink
MAA 212-002 Light Orange	Unmarked 212-011 White	CHOLETEC 212-023 Orange
SC 212-003 Yellow	HIDA 212-012 Brown	Tc RBC 212-026 White
DTPA 212-004 Olive Green	HDP 212-013 Purple	MYOVUE 212-027 White
GH 212-005 Red	MIBI 212-014 White	Y-90 212-029 White
Tc04 212-006 Pantone Purple	DMSA 212-015 Green	In-111 212-030 White
PYP 212-007 Blue	MAG-3 212-016 Blue	Saline 212-031 White
Ga 212-008 Orange	I-131 212-017 White	FDG 212-032 White
Neurolite 212-028 White		

This color-coded system has been designed to quickly identify radiopharmaceuticals. The .5" (1.3 cm) diameter circles are colored, labeled with the radiopharmaceutical name, and easily affixed to the plunger top of syringes. The colors have been selected as a standard most representative of radiopharmaceutical manufacturers' packaging schemes.

Color-Coded Syringe Labels (500/roll):

212-001 Label, MDP
212-002 Label, MAA
212-003 Label, SC
212-004 Label, DTPA
212-005 Label, GH
212-006 Label, Tc04
212-007 Label, PYP
212-008 Label, Ga
212-009 Label, TI
212-011 Label, (blank)
212-012 Label, HIDA
212-013 Label, HDP
212-014 Label, MIBI
212-015 Label, DMSA
212-016 Label, Mag-3
212-017 Label, I-131
212-019 Label, Cardiolite
212-023 Label, Choletec
212-026 Label, Tc RBC
212-027 Label, Myoview
212-028 Label, Neurolite
212-029 Label, Y-90
212-030 Label, In-111
212-031 Label, Saline
212-032 Label, FDG

ZEVALIN™ In-111 UNIT DOSE PIG



The 001-789 Unit Dose Pig is designed to reduce exposure from gamma emitting radiopharmaceuticals, such as the In-111 component of a Zevalin™ treatment. The barrel of the pig is constructed of .5" thick (1.3 cm) lead encased in durable Lexan. The pig accommodates a 10 cc B-D syringe filled to capacity. The pig fits into Zevalin shipping container and most ammo cans.

SPECIFICATIONS:

Dimensions: 10.3" l x 2.3" dia. (26 x 5.8 cm)

Lead Shielding: .5" thick (1.3 cm)

Weight: 8.2 lb (3.7 kg)

001-789 Pig, Unit Dose, Zevalin In-111
Accommodates 10 cc syringe with or without needle, filled to 10 cc

Related:

001-284 Pig Rack, PET

001-771 Sheets, Absorbent, 100/pkg

U.S. Patent No. 6,822,253

U.S. Patent No. 6,963,073

U.S. Patent No. 7,019,317

ZEVALIN™ SHIPPING CONTAINER



The Zevalin™ Shipping Container is designed to transport up to four unit dose pigs of Zevalin In-111 and/or Y-90. The system consists of a DOT type shipping container and foam inserts.

SPECIFICATIONS:

Dimensions:

Container: 11.75" l x 11.75" w x 12.5" h

(29.8 x 29.8 x 31.8 cm)

Cubic Feet: ~1 cu ft (.03 cu meters)

Regulations:


- Meets DOT II Type A packaging requirements when shipping Zevalin™.

- Meets IATA Dangerous Goods Regulations, 53rd Edition Sections 5.0.4.3, 10.5, and 10.6.1 thru 10.6.3.5

Compliance reports for packaging are available upon request.



ZEVALIN™ Y-90 UNIT DOSE PIG



The 001-788 Unit Dose Pig is designed to reduce exposure from beta emitting radiopharmaceuticals, such as the Y-90 component of a Zevalin™ treatment. The barrel of the pig is constructed of .09" thick (2.3 mm) lead and .36" (9.1 mm) acrylic that attenuates beta emissions and errant bremsstrahlung. The lead and acrylic is encased in a durable Lexan. The pig accommodates a 10 cc B-D syringe filled to 9 cc. The pig fits into the Zevalin shipping container and most ammo cans.

SPECIFICATIONS:

Dimensions: 10.3" l x 2.3" dia. (26 x 5.8 cm)

Lead Shielding: .09" thick (2.3 mm)

Acrylic Shielding: .36" thick (9.1 mm)

Weight: 3 lb (1.4 kg)

001-788 Pig, Unit Dose, Zevalin Y-90
Accommodates 10 cc syringe with or without needle, filled to 9 cc

Related:

001-284 Pig Rack, PET

001-771 Sheets, Absorbent, 100/pkg

U.S. Patent No. 6,822,253

U.S. Patent No. 6,963,073

U.S. Patent No. 7,019,317

001-790 Shipping Container, Zevalin
Note: Unit Dose Pigs sold separately

Related:

001-788 Pig, Unit Dose, Zevalin Y-90
Accommodates 10 cc syringe with or without needle, filled to 9 cc

001-789 Pig, Unit Dose, Zevalin In-111
Accommodates 10 cc syringe with or without needle, filled to 10 cc

001-726 Tags, Wire Security, 25/pkg
Used to identify unauthorized access

001-721 Document Protector

Int'l 631-924-9000 • www.biodex.com

PRO-TEC® UNIT DOSE PIG

Safety and efficiency – from the pharmacy to the patient... and back!



- Attractive design
- Overlapped lead eliminates streaming
- Durable, high-impact Lexan encapsulates lead components
- Opens and closes with a single twist
- 0.25" lead shielding
- O-ring seal
- Fits into ammo cans (vertically and horizontally)
- Compatible with automatic washing equipment
- Accommodates the following syringes:
 - 3 cc BD
 - 3 cc BD Safety-Lok
 - 3 cc Monoject
 - 3 cc Monoject Safety
 - 5 cc BD
 - 6 cc Monoject
 - 10 cc BD

The Pro-Tec Unit Dose Pig accommodates the most commonly used conventional and safety-engineered syringes. The lead components are fully encapsulated in durable Lexan, making this unit rugged and easy to clean. A single twist to open or close reduces loading/unloading time. The overlap design eliminates streaming regardless of the dose's position inside the pig.

A replaceable O-ring protects against leakage. Complete encapsulation of lead components in high-impact Lexan protects the lead shielding from physical damage during handling. The durable Lexan will not be damaged by automatic washing systems. Smooth plastic surfaces make the application and removal of adhesive labels easy. This product is designed to comply with IATA and DOT II requirements when transported in an appropriate shipping container.



Designed to fit conventional syringes



Designed to fit safety syringes

Innovative and cost effective, the Pro-Tec® Unit Dose Pig will help improve the safety and efficiency of radiopharmaceutical handling procedures

Pro-Tec® shielding for safety syringes will help pharmacies and clinics adhere to ALARA principles and improve compliance with OSHA directives – without compromising efficiency.

SPECIFICATIONS:

Dimensions: 8.5" h x 1.94" dia. (22 x 4.9 cm)

Lead Shielding: 0.25" nominal thickness (0.64 cm)

Construction: Lead, fully encapsulated with polycarbonate on the outside and polypropylene on the inside

Color: Red, White, Blue

Weight: 3 lb (1.4 kg)

001-280 Unit Dose Pig, Pro-Tec®, Red

001-281 Unit Dose Pig, Pro-Tec®, White

001-282 Unit Dose Pig, Pro-Tec®, Blue

Note: Unit Dose Pigs are sold in lots of six

U.S. Patent No. 6,822,253

U.S. Patent No. 6,963,073

U.S. Patent No. 7,019,317

Related:

008-400 Wall Rack, Unit Dose Pig

001-283 Pig Rack, Pro-Tec®

001-754 Shipping Bag, Pro-Tec Pig

066-533 Forceps, Curved, Locking, 9.5" l (24.1 cm)

PRO-TEC® PIG SHIPPING BAG



Accommodates up to 11 Pro-Tec Pigs

The durable, nylon, waterproof shipping bag is designed to accommodate up to 11 Pro-Tec Pigs containing syringes. It is the connecting piece in the Biodex solution, ensuring safe, convenient handling of syringes from the pharmacy and back again. The system meets DOT 7A Type A packaging requirements.

SPECIFICATIONS:

Dimensions: 10.5" h x 5.75" w x 11.5 depth (26.7 x 14.6 x 29.2 cm)

Material:

Exterior: Nylon

Interior: Polyester

001-754 Shipping Bag, Pro-Tec® Pig

Related:

001-280 Unit Dose Pig, Pro-Tec®, Red

001-281 Unit Dose Pig, Pro-Tec®, White

001-282 Unit Dose Pig, Pro-Tec®, Blue

Replacement:

001-756 Placard, Type I, 100/pkg

001-757 Placard, Type II, 100/pkg

001-758 Covers, hook & loop fasteners, plastic, Sm, 90 x 130 mm, 50/pkg

001-759 Covers, hook & loop fasteners, plastic, Md, 130 x 40 mm, 50/pkg

001-779 Covers, hook & loop fasteners, plastic, Lg, 160 x 180 mm, 50/pkg

U.S. Patent No. 6,822,253

U.S. Patent No. 6,963,073

U.S. Patent No. 7,019,317



UNIT DOSE PIG WALL RACK

Improve lab safety, efficiency and organization with the new Unit Dose Pig Wall Rack..

See page 11 for details.

PET PIG RACK



Accommodates PET, Zevalin and BEXXAR pigs.

Designed for countertop use, the PET Pig Rack holds up to eight PET pigs safely and conveniently.

The bottom of the rack has hexagonal cutouts that match the bottom of a pig to ensure each pig is held securely. When inserted into the rack, the pig can be easily opened or closed with a single twist.

SPECIFICATIONS:

Dimensions: 12" w x 1.75" h x 6" d (30.5 x 4.4 x 15.2 cm)

Weight: 2.8 lb (1.3 kg)

001-284 Pig Rack, PET

Related:

001-785 Pig, Unit Dose, PET, 3/5 cc, .5" lead

001-788 Pig, Unit Dose, Zevalin Y-90

001-789 Pig, Unit Dose, Zevalin In-111

001-860 Pig, Unit Dose, BEXXAR

U.S. Patent No. 6,822,253

U.S. Patent No. 6,963,073

U.S. Patent No. 7,019,317

PRO-TEC® PIG RACK



Designed for countertop use, the Pro-Tec Pig rack holds up to ten Pro-Tec Pigs safely and conveniently. The bottom of the rack has hexagonal cutouts that match the bottom of a Pro-Tec Pig to ensure each pig is held securely. When inserted into the rack, the pig can be easily opened or closed with a single twist.

SPECIFICATIONS:

Dimensions: 12" w x 1.75" h x 6" d (30.5 x 4.4 x 15.2 cm)

Weight: 2.8 lb (1.3 kg)

001-283 Pig Rack, Pro-Tec®

Related:

001-280 Unit Dose Pig, Pro-Tec®, Red

001-281 Unit Dose Pig, Pro-Tec®, White

001-282 Unit Dose Pig, Pro-Tec®, Blue


U.S. Patent No. 6,822,253

U.S. Patent No. 6,963,073

U.S. Patent No. 7,019,317

BEXXAR™ SHIPPING SYSTEM

Transports Unit Dose Pigs for BEXXAR



- *Single twist thread to open and close*
- *No exposed lead, durable, high-impact Lexan shell with painted lead interior*

The BEXXAR Shipping System transports one Unit Dose Pig that houses a 60 cc syringe for the radiotherapy drug BEXXAR. The pig is designed to reduce exposure while shipping I-131. The body of the pig is constructed of lead shielding encased in a durable Lexan shell and will hold a 60 cc syringe filled to 30 cc. The lead thickness of the pig, combined with the lead in the shipping container has been calculated to allow less than 50 mR at the surface of an ammo can. The system holds 160 mCi and meets DOT II shipping requirements.

SPECIFICATIONS:

001-860 Unit Dose Pig

Dimensions: 10.30" l x 2.50" dia (26.2 x 6.35 cm)

Lead Shielding:

Top: .56" (1.4 cm)

Side: .26" (.66 cm)

Bottom: .50" (1.27 cm)

Weight: 6.5 lb (3 kg)

001-880 Shipping System, BEXXAR

Dimensions:

Container: Ammo Can exterior 12" l x 6" w x 13.5 h

(30.5 x 15.2 x 34.3 cm)

Weight (Combined): 31 lb (14.1 kg)



Regulations:

- Meets DOT II Type A packaging requirements when shipping up to 160 mCi of I-131 BEXXAR™
- Complies with USA DOT 7A Type A Radioactive Material Requirements
- Meets IATA Dangerous Goods Regulations, 53rd Edition Sections 5.0.4.3, 10.5, and 10.6.1 thru 10.6.3.5



001-880 Shipping System, Unit Dose, BEXXAR

Includes: One 001-860 Unit Dose Pig

Component:

001-860 Pig, Unit Dose, BEXXAR

Accommodates a 60 cc syringe

Related:

001-284 Pig Rack, PET

001-771 Sheets, Absorbent, 100/pkg

001-861 Discs, Charcoal Absorbent, 100/pkg

001-721 Document Protector

U.S. Patent No. 6,963,073



FOR THE RADIOPHARMACIST

BEXXAR Dose Preparation System consists of a specially designed Syringe Shield, Vial Shield and Stand. Together this system provides a safe and shielded environment while making it easy to prepare a BEXXAR dose.
See page 58



FOR THE TECHNOLOGIST

Syringe Dipper accommodates B-D and Monoject 60 cc syringes to assay the dose prior to administration. The dipper fits in most commercial dose calibrators with 2.5" x 10" well chamber interior.

See page 75



Graseby Injection Pump Shield provides .5"-thick lead shielding to attenuate I-131 during the slow infusion process. The pump shield accommodates Graseby units #3100 and #3400.

See page 13



I-131 SODIUM IODIDE SHIPPING SYSTEM

Transports Unit Dose Vial Pigs for I-131 Sodium Iodide



Transport up to a 30 ml vial of I-131 Sodium Iodide in the I-131 Shipping System. The complete system offers 1" lead shielding to accommodate up to 235 mCi of I-131 and meets DOT II Type A packaging requirements.

The I-131 Shipping System transports one Unit Dose Vial Pig that accommodates an 8 to 10 ml or 30 ml vial containing I-131 Sodium Iodide capsules. The vial pig is designed to reduce exposure while shipping diagnostic or therapy doses of I-131 in capsule form. The vial pig offers 1" of lead shielding to accommodate the concentrated energy.

The system meets DOT II Type A packaging requirements when shipping up to 235 mCi of I-131 capsules.

I-131 Shipping System for Vial Pigs consists of:

- Vial Pig
- Absorbent Sheets
- Shipping Container
- Vial Pig Adapter to accommodate 8 to 10 ml vials



NEW

TABLETOP SHIELD

For detailed product specifications, see page 22.

SPECIFICATIONS:

001-706 Vial Pig

Dimensions:

Exterior: 6.63" h x 4.15" dia (16.8 x 10.5 cm)

Interior: 2.76" h x 1.51" dia (7 x 3.8 cm)

Lead Shielding:

Sides and Bottom: 1" thick (2.5 cm)

Top: 1.75" (4.4 cm)

Weight: 21.3 lb (9.7 kg)

001-881 Shipping System, I-131

Dimensions:

Container: Ammo Can exterior 12" l x 6" w x 13.5 h

(30.5 x 15.2 x 34.3 cm)

Weight (Combined): 31 lb (14.1 kg)



Regulations:

- Meets DOT II Type A packaging requirements when shipping up to 235 mCi of I-131 Sodium Iodine
- Complies with USA DOT 7A Type A Radioactive Material Requirements
- Meets IATA Dangerous Goods Regulations, 53rd Edition Sections 5.0.4.3, 10.5, and 10.6.1 thru 10.6.3.5

Compliance reports for radioactive materials packaging are available by request or visit our website: www.biodex.com.

001-881 Shipping System, I-131, Vial

For 8 to 10 or 30 ml vials

Includes: One 001-706 Vial Pig and shipping container

Component:

001-706 Pig, Vial, PET, 10/30 ml, 1" lead

For 30 ml vials

*Includes: 001-707 Vial Pig Adapter, to accommodate 8 to 10 ml vials, and Three Absorbent Sheets**

Related:

001-726 Tags, Wire Security, 25/pkg

Used to identify unauthorized access

001-771 Sheets, Absorbent, 100/pkg

001-707 Adapter, Vial Pig, 10 ml

Allows 001-706 Vial Pig to accommodate 10 ml vials

001-721 Document Protector

COMPACT PET SHIPPING SYSTEM FOR ONE UNIT DOSE PIG




This Compact PET Shipping System transports one 3 cc or 5 cc syringe containing high-energy radionuclides such as FDG F-18. Syringes fit into the pig with or without an attached needle. Designed to conserve space and minimize weight, the fully loaded shipping container weighs less than 40 lb. An important feature is that the shipping container can be left at a convenient height while the pig (9 pounds) can be easily removed from the container. The pig is then placed behind an L-Block for dose loading and unloading. All Biodex PET L-Block Shields incorporate a hex-shaped plate that corresponds with the hex shape at the base of the pig. This design element allows one-handed loading and unloading of syringes.

The 001-785 Unit Dose Pig is encapsulated in high-impact Lexan and polypropylene, making the pig durable, easy to clean and compatible with automatic washing systems. A single twist opens or closes the pig, reducing handling time.

The system meets DOT II Type A packaging requirements when shipping up to 500 mCi (18.5 GBq) of FDG F-18.

Compact PET Shipping System for One Unit Dose Pig consists of:

- PET Unit Dose Pig
- Absorbent Sheets
- Shipping Container with lead shielding



001-785 PET Unit Dose Pig also sold separately.

The PET Unit Dose Pig is encapsulated in durable, high-impact Lexan and polypropylene:

- Single twist thread to open and close
- No exposed lead
- Compatible with automatic washing equipment

SPECIFICATIONS:

001-785 Unit Dose Pig

Dimensions: 10.2" h x 2.4" dia (26 x 6 cm)

Lead Shielding:

Body: .5" thick (1.3 cm)

Ends: 1.2" thick, bottom; 1.44" thick, top (3 cm and 3.6 cm)

Weight: 8.7 lb (4 kg)

001-786 PET Shipping System, Single

Dimensions:

Container: 11.75" l x 11.75" w x 12.5" h (29.8 x 29.8 x 31.8 cm)

Cubic Feet: ~1 cu ft (.03 cu meters)

Weight (Combined): 38.3 lb (17.5 kg)



Regulations:

- Meets DOT II Type A packaging requirements when shipping up to 500 mCi (18.5 GBq) of FDG F-18



- Meets IATA Dangerous Goods Regulations, 53rd Edition Sections 5.0.4.3, 10.5, and 10.6.1 thru 10.6.3.5

Compliance reports for radioactive materials packaging are available by request or visit our website: www.biodex.com.

001-786 Shipping System, PET, Single

For single Unit Dose Pig

Includes: 001-785 Unit Dose Pig, absorbent sheets and shipping container with lead shielding

Component:

001-785 Pig, Unit Dose, PET, 3/5 cc, .5" lead

Accommodates syringes with or without needle

Note: PET Unit Dose Pigs are sold in lots of three.

Related:

001-284 Pig Rack, PET

001-730 Cart, Transport, PET Shipping System

001-771 Sheets, Absorbent, 100/pkg

001-726 Tags, Wire Security, 25/pkg

Used to identify unauthorized access

001-721 Document Protector

U.S. Patent No. 6,586,758

U.S. Patent No. 6,822,253

U.S. Patent No. 6,963,073

U.S. Patent No. 7,019,317

COMPACT PET SHIPPING SYSTEM FOR TWO UNIT DOSE PIGS




This Compact PET Shipping System transports two 3 cc or 5 cc syringes containing high-energy radionuclides such as FDG F-18. Syringes fit into the pigs with or without an attached needle. Designed to conserve space and minimize weight, the fully loaded shipping container weighs only 55 lb. An important feature is that the shipping container can be left at a convenient height while the pig (9 pounds) can be easily removed from the container. The pig is then placed behind an L-Block for dose loading and unloading. All Biodex PET L-Block Shields incorporate a hex-shaped plate that corresponds with the hex shape at the base of the pig. This design element allows one-handed loading and unloading of syringes.

The 001-785 Unit Dose Pig is encapsulated in high-impact Lexan and polypropylene, making the pig durable, easy to clean and compatible with automatic washing systems. A single twist opens or closes the pig, reducing handling time.

The system meets DOT II Type A packaging requirements when shipping up to 160 mCi (5.92 GBq) in one pig and 235 mCi (8.70 GBq) in the other, totaling 395 mCi (14.62 GBq) of FDG F-18.

Compact PET Shipping System
for Two Unit Dose Pigs consists of:

- PET Unit Dose Pigs (two)
- Absorbent Sheets
- Shipping Container with lead shielding



001-785 PET
Unit Dose Pig also
sold separately.

The PET Unit Dose Pig is encapsulated in durable, high-impact Lexan and polypropylene:

- Single twist thread to open and close
- No exposed lead
- Compatible with automatic washing equipment

SPECIFICATIONS:

001-785 Unit Dose Pig

Dimensions: 10.2" h x 2.4" dia (26 x 6 cm)

Lead Shielding:

Body: .5" thick (1.3 cm)

Ends: 1.2" thick, bottom;

1.44" thick, top (3 cm and 3.6 cm)

Weight: 8.7 lb (4 kg)

001-787 PET Shipping System, Double

Dimensions:

Container: 11.75" l x 11.75" w

x 12.5" h (29.8 x 29.8 x 31.8 cm)

Cubic Feet: ~1 cu ft (.03 cu meters)

Weight (Combined): 55.2 lb (25.2 kg)



Regulations:

- Meets DOT II Type A packaging requirements when shipping up to 160 (5.92 GBq) and 235 mCi (8.70 GBq) of FDG F-18



- Meets IATA Dangerous Goods Regulations, 53rd Edition Sections 5.0.4.3, 10.5, and 10.6.1 thru 10.6.3.5

Compliance reports for radioactive materials packaging are available by request or visit our website: www.biodex.com.

001-787 Shipping System, PET, Double

For two single Unit Dose Pigs

Includes: Two 001-785 Unit Dose Pigs, absorbent sheets and shipping container with lead shielding

Component:

001-785 Pig, Unit Dose, PET, 3/5 cc, .5" lead

Accommodates syringes with or without needle

Note: PET Unit Dose Pigs are sold in lots of three.

Related:

001-284 Pig Rack, PET

001-730 Cart, Transport, PET Shipping System

001-771 Sheets, Absorbent, 100/pkg

001-726 Tags, Wire Security, 25/pkg

Used to identify unauthorized access

001-721 Document Protector

U.S. Patent No. 6,586,758

U.S. Patent No. 6,822,253

U.S. Patent No. 6,963,073

U.S. Patent No. 7,019,317

COMPACT PET SHIPPING SYSTEM FOR THREE UNIT DOSE PIGS




This Compact PET Shipping System transports three 3 cc or 5 cc syringes containing high-energy radionuclides such as FDG F-18. Syringes fit into the pigs with or without an attached needle. Designed to conserve space and minimize weight, the fully loaded shipping container weighs only 95 lb. An important feature is that the shipping container can be left at a convenient height while the pig (9 pounds) can be easily removed from the container. The pig is then placed behind an L-Block for dose loading and unloading. All Biodex PET L-Block Shields incorporate a hex-shaped plate that corresponds with the hex shape at the base of the pig. This design element allows one-handed loading and unloading of syringes.

The 001-785 Unit Dose Pig is encapsulated in high-impact Lexan and polypropylene, making the pig durable, easy to clean and compatible with automatic washing systems. A single twist opens or closes the pig, reducing handling time.

The system meets DOT II Type A packaging requirements when shipping up to 235 mCi (8.70 GBq) in one pig, and 160 mCi (5.92 GBq) in the second pig and 140 mCi (5.18 GBq) in the third pig, totaling 535 mCi (19.80 GBq) of FDG F-18.

Compact PET Shipping System for three Unit Dose Pigs consists of:

- PET Unit Dose Pigs (three)
- Absorbent Sheets
- Shipping Container with lead shielding



001-785 PET Unit Dose Pig also sold separately.

The PET Unit Dose Pig is encapsulated in durable, high-impact Lexan and polypropylene:

- Single twist thread to open and close
- No exposed lead
- Compatible with automatic washing equipment

SPECIFICATIONS:

001-785 Unit Dose Pig

Dimensions: 10.2" h x 2.4" dia (26 x 6 cm)

Lead Shielding:

Body: .5" thick (1.3 cm)

Ends: 1.2" thick, bottom; 1.44" thick, top (3 cm and 3.6 cm)

Weight: 8.7 lb (4 kg)

001-739 PET Shipping System, Triple

Dimensions:

Container: 11.75" l x 11.75" w x 12.5" h (29.8 x 29.8 x 31.8 cm)

Cubic Feet: ~1 cu ft (.03 cu meters)

Weight (Combined): 95 lb (43.1 kg)



Regulations:

- Meets DOT II Type A packaging requirements when shipping up to 235 (8.70 GBq), 160 (5.92 GBq) and 140 mCi (5.18 GBq) of FDG F-18
- Meets IATA Dangerous Goods Regulations, 53rd Edition Sections 5.0.4.3, 10.5, and 10.6.1 thru 10.6.3.5

Compliance reports for radioactive materials packaging are available by request or visit our website: www.biodex.com.

001-739 Shipping System, PET, Triple

For three single Unit Dose Pigs

Includes: Three 001-785 Unit Dose Pigs, absorbent sheets and shipping container with lead shielding

Component:

001-785 Pig, Unit Dose, PET, 3/5 cc, .5" lead

Accommodates syringes with or without needle

Note: PET Unit Dose Pigs are sold in lots of three.

Related:

001-284 Pig Rack, PET

001-730 Cart, Transport, PET Shipping System

001-771 Sheets, Absorbent, 100/pkg

001-726 Tags, Wire Security, 25/pkg

Used to identify unauthorized access

001-721 Document Protector

U.S. Patent No. 6,586,758

U.S. Patent No. 6,822,253

U.S. Patent No. 6,963,073

U.S. Patent No. 7,019,317

INTEGO™ VIAL SHIPPING SYSTEM


The PET shipping system for use with the Medrad® Intego PET infusion system.



Biodex has designed a NEW Vial Shipping System that employs a unique Tungsten Vial Shield for use with the Medrad® Intego™ PET Infusion System. Manufactured to Medrad's specifications, the Vial Shield transports a 30 ml Hospira vial.

To lift or lower the vial within the container or the Intego™ Infusion System, a detachable handle is provided. Accommodation is made for the handle to travel with the shipping container.

The system meets DOT II Type A packaging requirements when shipping up to 2.5 Ci (92.5 GBq) of FDG F-18.



SPECIFICATIONS:

001-708 Vial Shield with Lifting Handle
Weight: 15.2 lb (6.89 kg)

001-723 Intego Shipping Container

Dimensions:

Container: 11.75" l x 11.75" w x 12.5" h
(29.8 x 29.8 x 31.8 cm)

Cubic Feet: ~1 cu ft (.03 cu meters)

Weight: 48.8 lb (22.1 kg)

Weight (Combined): 64 lb (29 kg)



Regulations:

- Meets DOT II Type A packaging requirements when shipping up to 2.5 Ci (92.5 GBq) of FDG F-18



- Meets IATA Dangerous Goods Regulations, 53rd Edition Sections 5.0.4.3, 10.5, and 10.6.1 thru 10.6.3.5

Compliance reports for radioactive materials packaging are available by request or visit our website: www.biodex.com.

001-708 Shield, Vial, Tungsten, Intego™

For 30 ml Hospira vials.

Includes lifting handle.

001-723 Shipping Container, Intego™

Related:

001-730 Cart, Transport, PET Shipping System

001-726 Tags, Wire Security, 25/pkg

Used to identify unauthorized access

001-721 Document Protector

U.S. Patent No. 6,586,758

COMPACT PET SHIPPING SYSTEM FOR VIAL PIG




Transport up to a 30 ml vial of PET pharmaceutical in the PET Shipping System. The complete system offers 1.55" lead shielding to accommodate 2.5 curies of FDG F-18 and meets DOT II Type A packaging requirements.

The PET Shipping System for Vial Pigs is designed to transport a 10 ml or 30 ml vial containing high-energy radionuclides. Designed to conserve space and minimize weight, the entire system weighs only 50 lb. An important feature is that the shipping container can be placed at a convenient height while the pig is easily removed from the shipping case. For added safety and convenience, the vial pig can then be placed in the 042-466 Dose Drawing System for drawing doses from the vial.

The system meets DOT II Type A packaging requirements when shipping up to 2.5 Ci (92.5 GBq) of FDG F-18.

PET Shipping System for Vial Pigs consists of:

- PET Vial Pig
- Absorbent Sheets
- Shipping Container with lead shielding



SPECIFICATIONS:

001-706 Vial Pig

Dimensions:

Exterior: 6.63" h x 4.15" dia (16.8 x 10.5 cm)

Interior: 2.76" h x 1.51" dia (7 x 3.8 cm)

Lead Shielding:

Sides and Bottom: 1" thick (2.5 cm)

Top: 1.75" (4.4 cm)

Weight: 21.3 lb (9.7 kg)

001-724 PET Shipping System, Vial

Dimensions:

Container: 11.75" l x 11.75" w x 12.5" h
(29.8 x 29.8 x 31.8 cm)

Cubic Feet: ~1 cu ft (.03 cu meters)

Weight (Combined): 49.7 lb (22.5 kg)



Regulations:

- Meets DOT II Type A packaging requirements when shipping up to 2.5 Ci (92.5 GBq) of FDG F-18



- Meets IATA Dangerous Goods Regulations, 53rd Edition Sections 5.0.4.3, 10.5, and 10.6.1 thru 10.6.3.5

Compliance reports for radioactive materials packaging are available by request or visit our website: www.biotech.com.

001-724 Shipping System, PET, Vial*

For 10 or 30 ml vials

Includes: 001-706 Vial Pig and shipping container with lead shielding

Components:

001-706 Pig, Vial, PET, 10/30 ml, 1" lead

For 30 ml vials

Includes: 001-707 Vial Pig Adapter, to accommodate 10 ml vials, and Three Absorbent Sheets*

Related:

001-707 Adapter, Vial Pig, 10 ml

Allows 001-706 Vial Pig to accommodate 10 ml vials

001-730 Cart, Transport, PET Shipping System

001-771 Sheets, Absorbent, 100/pkg

001-726 Tags, Wire Security, 25/pkg


Used to identify unauthorized access

001-721 Document Protector

U.S. Patent No. 6,586,758

COMPACT PET SHIPPING SYSTEM FOR DOUBLE-ENDED PET PIG

One system for shipping and injecting FDG F-18



The Double-Ended PET Pig is an all-in-one solution for the safe transport and administration of 511 keV radionuclides, such as FDG F-18. The pig is constructed in three sections of .6" (1.5 cm) lead enclosed in durable Lexan. Depending on the stage of the injection process, only one section of the pig is removed. The injection itself is accomplished without removing the syringe from the pig. Simply open the administration port and push the tungsten plunger against the syringe plunger. During injection the pig is positioned on a stand. The Double-Ended Pig accommodates a 5 cc syringe and fits into its own compact shipping container. The system employs a unique design that decreases the weight and size of the container.

Here's how it works:

1. Remove the Double-Ended Pig from the compact shipping container. (See dia. A.)
2. Place behind an L-Block Shield.
3. Unscrew the top. (See dia. B.)
4. Remove the syringe and place into a dose calibrator. (See dia. C.)
5. Return the syringe to the pig and put the top back on. (See dia. D.)
6. Transport the pig to the patient injection area. (See dia. E.)
7. Place the pig into the Pig Cradle so the top section of the pig is over the long section of the cradle. (See dia. F.)
8. Remove the bottom of the pig and connect to your preferred injection device, butterfly, etc. (See dia. G.)
9. Open the plunger lock located on the top section by pulling the slide toward the edge of the pig. With a pen-like device, push the tungsten plunger to administer the dose. (See dia. G.)
10. Recap the syringe.
11. Put the bottom back on the pig.
12. Return the pig to the compact shipping container. (See dia. H.)

SPECIFICATIONS:

001-793 Double-Ended PET Pig

Dimensions: 9.5" h x 2.3" dia at maximum point (24 x 5.8 cm)

Shielding:

Sides: .6" thick (1.5 cm) lead

Bottom: 1.44" thick (3.7 cm) lead

Top: .875" thick (2.2 cm) tungsten

Weight: 9 lb (4 kg)

001-794 PET Shipping System, Double-Ended Pig

Dimensions:

Container: 11.75" l x 11.75" w x 12.5" h (29.8 x 29.8 x 31.8 cm)

Cubic Feet: ~1 cu ft (0.3 cu meters)

Weight: 36 lb (16.3 kg)



Regulations:

- Meets DOT II Type A packaging requirements when shipping up to 475 mCi (17.58 GBq) of FDG F-18.
- Meets IATA Dangerous Goods Regulations, 53rd Edition Sections 5.0.4.3, 10.5, and 10.6.1 thru 10.6.3.5.



Compliance reports for radioactive materials packaging are available by request or visit our website: www.biodex.com.

001-794 Shipping System, PET, Double-Ended Pig, Single
Includes one 001-793 Double-Ended PET Pig and Shipping Container with lead.

Component:

001-793 Pig, Double-Ended, PET, .6" lead
Accommodates 5 cc syringes with or without needle.

Related:

001-797 Pig Cradle
Accommodates 001-793.

001-771 Sheets, Absorbent, 100/pkg

001-721 Document Protector

U.S. Patent No. 6,586,758 B2


Other patents pending

PET DOSE DRAWING SYSTEM

For drawing FDG F-18 from a vial into a syringe



042-466 Dose Drawing System



042-467 Dose Drawing System with 2" L-Block

The Dose Drawing System is used to draw FDG F-18 doses from a vial. It consists of a specially designed Dose Drawing Syringe Shield, vial shield and stand.

When drawing FDG F-18, the vial shield and syringe shield remain on the stand. The stand allows the vial shield to rotate and has a fixed stop at the 45-degree downward angle. The vial shield is constructed of 1" lead encased in steel. The shield's top is threaded to allow quick insertion and removal of vials. The top cover incorporates a pivoting shield for septum access that allows the use of a vent needle, if desired, along with the needle from a 5 cc syringe.

The tungsten syringe shield mounts on a stainless steel rotating platform that slides the syringe into the vial and locks in position. The dose can then be drawn with forceps to minimize hand exposure. The syringe shield is constructed of .55" tungsten. A high-density lead glass window allows viewing up to the 5 cc mark on a 5 cc syringe.

The PET Dose Drawing System provides a safe and shielded environment while making it easy to draw a FDG F-18 dose.



If transported in a Biodex PET Vial Shipping System, the vial shield is placed in the stand.



Vial shield shipping cover is replaced with the drawing shield cover.



Syringe is placed in the Dose Drawing Syringe Shield. The vial shield rotates to 45-degree downward angle. Swing open the cover for vial septum access.



Syringe shield slides upward, allowing the syringe to puncture the vial septum. Dose is drawn.

SPECIFICATIONS:**007-984 Dose Drawing Syringe Shield**

Dimensions: 1.7" dia x 5.5" l (4.3 x 14 cm)

Shielding: .55" thick (1.3 cm) tungsten

Lead Glass Window: .25" thick (.64 cm), 5.6 density

Construction:

Exterior: Tungsten

End Cap Mechanism: Stainless steel

Weight: 4 lb (1.8 kg)

001-855 Dose Drawing Vial Shield

Dimensions:

Exterior: 6.1" h x 4" dia (15.5 x 10.2 cm)

Interior: 2.76" h x 1.51" dia (7 x 3.8 cm)

Lead Shielding:

Cylinder: 1" thick (2.54 cm)

Top End Cap: 1" thick (2.5 cm)

Cover: .75" thick (1.9 cm)

Construction:

Exterior: Steel

Interior: Molded lead

Top End Cap: Aluminum

Bottom End Cap: Plastic

Finish: Powder coat

Weight: 20 lb (9 kg)

042-455 Dose Preparation Stand

Dimensions: 12.6" w x 10.6" depth x 12.6" h (32 x 26.9 x 32 cm)

Construction: Stainless steel and steel

Finish: Stainless steel and powder coat

Weight: 21 lb (9.5 kg)

042-449 2" L-Block Shield

Dimensions: 14" w x 18" depth x 24.7" h (36 x 46 x 62 cm)

Lead Shielding:

Front: 2" thick (5.08 cm)

Base: 1" thick (2.5 cm)

Lead Glass Window:

Dimensions: 8" w x 8" h x 4" thick (20 x 20 x 10 cm)

Density: 5.2g/cm³

Finish: Powder Coat

Weight: 345 lb (156 kg)

Shipping Weight: 420 lb (191 kg)

042-426 Interlocking Lead Brick Cave

Dimensions:

I.D.: 14" w x 17.8" depth x 13.8" h (35.5 x 45.3 x 34.6 cm)

Lead Shielding: 2" thick (5 cm)

Finish: Paint

Weight: 532 lb (241 kg)

042-466 Dose Drawing System, PET

Includes: 007-984 Dose Drawing

Syringe Shield, 001-855 Dose Drawing

Vial Shield, 042-455 Dose Preparation Stand and

001-707 Vial adapter to accommodate 10 ml vials.

**042-467 Dose Drawing System, PET,
with L-Block Shield (042-449)****Components:**

007-984 Syringe Shield, Dose Drawing,
PET, 5 cc

001-855 Vial Shield, Dose Drawing,
PET, 30 ml

042-455 Stand, Dose Preparation

042-449 L-Block Shield, 2" lead
With 8" x 8" x 4" lead glass window

Related:

007-986 Syringe Shield, Dose Drawing,
PET, 3 cc

007-988 Syringe Shield, Dose Drawing,
PET, 2.5 cc

007-989 Syringe Shield, Dose Drawing,
PET, 5 cc

042-426 Lead Brick Cave, 3-wall, 2" lead
*Fits 042-428, 042-449 and 042-419 L-Block Shields
Accommodates 042-466 PET Dose Drawing System*

066-536 Forceps, Curved, non-locking,
12.5 l" (31.7 cm)

CYCLOTRON WORKBENCH



- Includes 2" lead storage safe
- Ideal for target decay
- Use in PET, nuclear medicine or radiation therapy departments

Sturdy as they come, this steel table can be used for just about any application requiring a strong, level platform. Ideal for holding heavy L-Block shields and caves, the surface is powder coated and the front legs feature adjustable levelers.

The middle shelf features a storage safe that is ideal for storing large quantities of high-energy radioisotopes and/or target decay. Shielded with a thickness of 2" of lead, the safe is encased in a powder-coated steel jacket and features an adjustable shelf. The lead-lined door is hung with heavy duty non-sagging hinges and is key-locked to prevent unauthorized access.

SPECIFICATIONS:**042-456 Cyclotron Workbench**

Dimensions: 36.75" w x 24" depth x 36" h (93.5 x 61 x 91.5 cm)

Front legs incorporate adjustable levelers

Finish: Powder coat

Weight Capacity: 1,600 lb (725.76 kg) includes 1,050 lb for safe

Weight: 1200 lb (544.32 kg)

Shipping Weight: 1250 lb (567 kg)

244-006 Lead Storage Safe, 2"

Dimensions: 17.4" w x 17" depth x 19" h (44.2 x 43.2 x 48.3 cm)

I.D.: 12" w x 12" depth x 12" h (30.5 x 30.5 x 30.5 cm)

Lead Shielding: 2" thick (5 cm)

Finish: Powder coat

Door: Key-locked

Weight: 1050 lb (476 kg)

042-456 Workbench, Cyclotron

Includes: Storage Safe, 2" lead (244-006)

BEXXAR™ DOSE PREPARATION SYSTEM



The Dose Preparation System is used to prepare the radioimmunotherapy drug BEXXAR. It consists of a specially designed syringe shield, vial shield and stand.

When preparing BEXXAR, the vial shield and syringe shield remain on the stand. The stand allows the

vial shield to rotate and has a fixed stop at the 45-degree downward angle. The vial shield is constructed of lead shielding encased in steel. The shield's top is threaded to allow quick insertion and removal of vials. The top cover incorporates a pivoting shield for septum access that allows the use of a vent needle along with the needle from the 60 cc syringe.

The syringe shield mounts on a rotating platform that slides the syringe into the vial and locks in position. The dose can then be drawn with forceps to minimize hand exposure. The syringe shield is constructed of .84" lead shielding encased in steel. A high-density lead glass window allows viewing up to the 50 cc mark on a 60 cc syringe. When the syringe is filled and assayed, it can be transferred to the BEXXAR Unit Dose Pig for safe transport.

The Dose Preparation System provides a safe and shielded environment while making it easy to prepare a BEXXAR dose.

SPECIFICATIONS:

007-978 Dose Preparation Syringe Shield

Dimensions: 3" dia x 5.6" l (7.6 x 14.2 cm)

Lead Shielding: .84" thick (2.1 cm)

Lead Glass Window:

Dimensions: .50" w x 4.25" l x .38" thick (1.27 x 10.80 x .97 cm)

Density: 5.6

Construction:

Exterior: Steel

End Cap Mechanism: Stainless steel

Finish: Powder coat

Weight: 13 lb (5.9 kg)

001-850 Dose Preparation Vial Shield

Dimensions:

Exterior: 6.35" h x 4" dia (16.1 x 10.2 cm)

Interior: 3.35" h x 1.17" dia (8.5 x 2.9 cm)

Lead Shielding:

Cylinder: 1.17" thick (3 cm)

Bottom: .81" thick (2.1 cm)

Top End Cap: 1" thick (2.5 cm)

Cover: .75" thick (1.9 cm)

Construction:

Exterior: Steel

Interior: Molded lead

Top End Cap: Aluminum

Bottom End Cap: Plastic

Finish: Powder coat

Weight: 21 lb (9.5 kg)

042-455 Dose Preparation Stand

Dimensions: 12.6" w x 10.6" d x 12.6" h (32 x 26.9 x 32 cm)

Construction: Stainless steel and steel

Finish: Stainless steel and power coat

Weight: 21 lb (9.5 kg)

042-460 Dose Preparation System, BEXXAR

Includes: 007-978 Dose Preparation

Syringe Shield, 001-850 Dose Preparation

Vial Shield and 042-455 Dose Preparation Stand

Component:

007-978 Syringe Shield, Dose Preparation, BEXXAR, 60 cc

001-850 Vial Shield, Dose Preparation, BEXXAR, 30 ml

042-455 Stand, Dose Preparation

Related:

066-536 Forceps, Curved, non-locking, 12.5 l" (31.7 cm)



Syringe shield mounts to the stand on a stainless steel rotating platform; syringe is placed in shield. Note vial shield in upright position.



The vial shield rotates to 45-degree downward angle. Swing open the cover for vial septum access.



Syringe shield slides upward allowing the syringe to puncture the vial septum. Dose is drawn.

LOW-ENERGY SHIELDED DECAY DRUM



Available with .125" or .25" shielding. Shown with optional Transport Dolly.

- *Completely shielded*
- *Convenient sliding access door*

The Shielded Decay Drum provides safe handling and storage of radioactive waste material.

A top-loading, sliding port provides fast, safe access. The drum includes three 30-gallon plastic liners.

With two Shielded Decay Drums on hand, a rotation system can be instituted for an even more practical storage and decay process.

The optional dolly simplifies transport of the decay drum.

SPECIFICATIONS:

Dimensions: 30" h x 19.5" dia (76 x 50 cm)

Lead Shielding: .125" (.3 cm) or .25" (.64 cm) lead

Top Sliding Port: 6.5" x 8.5" (16.5 x 22 cm)

Capacity: 30 gal

Weight:

039-284: 160 lb (72 kg) for .125"

039-244: 280 lb (128.8 kg) for .25"

Shipping Weight:

039-284: 238 lb (107 kg)

039-244: 353 lb (160 kg)

039-599 Transport Dolly

Weight: 14 lb (6.4 kg)

039-284 Decay Drum, .125" lead
Includes: Three 039-286 Plastic Liners

039-244 Decay Drum, .25" lead
Includes: Three 039-286 Plastic Liners

Related:

039-599 Dolly, Transport

Replacement:

039-286 Liner, Plastic, 30-gal capacity

HIGH-ENERGY SHIELDED DECAY DRUM



Shielded with .5" lead for high-energy isotopes.

- *Completely shielded*
- *Large drop port*

The Shielded Decay Drum provides the necessary shielding from high-energy radiopharmaceuticals, such as Ga-67, In-111 and I-131 products.

This heavy-duty drum has a handled lid with a 6.5" x 8" (16.5 x 20.3 cm) drop port. The drum includes three 30-gallon plastic liners.

With two Shielded Decay Drums on hand, a rotation system can be instituted for an even more practical storage and decay process.

The optional dolly simplifies transport of the decay drum.

SPECIFICATIONS:

Dimensions: 30" h x 19.5" dia (76 x 50 cm)

Lead Shielding: .5" thick (1.3 cm)

Drop Port: 6.5" x 8" (16.5 x 20.3 cm)

Weight: 480 lb (217.7 kg)

039-599 Transport Dolly

Weight: 14 lb (6.4 kg)

039-514 Decay Drum, 30" h, .5" lead
Includes: Three 039-286 Plastic Liners

Related:

039-599 Dolly, Transport

Replacement:

039-286 Liner, Plastic, 30-gal capacity

TEC-CONTROL™ CHROMATOGRAPHY SYSTEMS

For radiopharmaceutical quality control






Tec-Control™ Chromatography tests the radiochemical purity of specific Tc-99m-labeled radiopharmaceuticals. The accompanying chart shows which strips and solvents are required to perform each individual test. Some solvents must be purchased separately (see Sigma-Aldrich chart) due to hazardous material shipping restrictions.

Detailed instruction manuals are packaged with each strip container, although our Radiopharmaceutical QC Procedure Manual (151-000) explains paper chromatography in greater detail.

Radiopharmaceutical QC Manual

**MINIATURIZED
CHROMATOGRAPHY
PROCEDURES FOR
RADIOPHARMACEUTICALS:
2011 UPDATE
150-000**

BY
A. MICHAEL ZIMMER, PhD



150-000 Procedure Manual,
Radiopharmaceutical QC

This detailed manual explains Paper Chromatography, a QC method for evaluating the radiochemical purity of currently used Tc-99m-labeled radiopharmaceuticals. Procedures are quick and easy to use, a simple quality control solution for any nuclear medicine department.

TEC-CONTROL™ ALUMINUM BREAKTHRU KIT

Simple 3-step procedure



The Aluminum Breakthru Kit provides a rapid, easy and inexpensive way to test aqueous solutions, particularly pertechnetate generator eluate, for trace quantities of aluminum. Aluminum forms an intense red precipitate with the indicator paper, and the intensity of the color is directly proportional to the amount of aluminum in the solution. The USP allows a concentration of aluminum ion in an injection ≤ 10 micrograms per milliliter (10 $\mu\text{g}/\text{ml}$) in technetium 99m eluate prepared from Molybdenum 99 formed as a result of uranium fission.

Simple Procedure:

1. Place a drop of the eluate or solution to be tested on the indicator paper. The best procedure is to form a hanging drop using a 19-22G needle.
2. Place a drop of the standard aluminum solution on the indicator paper. Use the same size drop.
3. Compare the intensity of the red spot formed. If the eluate spot is less intense than the standard solution, the eluate contains less than 10 $\mu\text{g}/\text{ml}$ aluminum.

150-780 Chromatography Kit, Tec-Control
Aluminum Breakthru Kit includes:
*Aluminum standard, 5 ml, 10 $\mu\text{g}/\text{ml}$,
50 indicator strips and manual*

150-785 Chromatography Kit, Tec-Control
Aluminum Breakthru Kit includes:
*Aluminum standard, 5 ml, 5 $\mu\text{g}/\text{ml}$,
50 indicator strips and manual*

References:

Miniaturized Chromatographic Quality-Control Procedures for Tc-99m Radiopharmaceuticals; A. Michael Zimmer and Dan G. Pavel, Journal of Nuclear Medicine, Vol. 18/12, Dec. 1977, pg. 1230.
Technical Parameters Associated with Miniaturized Chromatography Systems; Raimund A. Taukulis, A. Michael Zimmer, Dan G. Pavel and Bhupendra A. Patel, University of Illinois Medical Center, Chicago, Illinois, Journal of Nuclear Medicine Technology, Vol. 7/1.

CHROMATOGRAPHY QC FOR THE FOLLOWING RADIOPHARMACEUTICALS:

RADIOPHARMACEUTICAL	SOLVENT(S) REQUIRED*	STRIPS REQUIRED
Aluminum Breakthru	150-781	150-782
Aluminum Breakthru	150-783	150-782
Bicisate (Neurolite™)	Ethyl Acetate 99.9%	150-130
Diphosphonate	Acetone (HPLC Grade), Distilled H ₂ O	150-001 & 150-005
Disofenin (Hepatolite™)	150-160, Distilled H ₂ O	150-122 & 150-125
DMSA	Acetone (HPLC Grade),	150-025
DTPA	Acetone (HPLC Grade), Distilled H ₂ O	150-001 & 150-005
Exametazine (Ceretek™)	Ethyl Acetate 99.9%	150-130
Glucoheptonate	Acetone (HPLC Grade), Distilled H ₂ O	150-001 & 150-005
HDP	Acetone (HPLC Grade), Distilled H ₂ O	150-001 & 150-005
MAA	Acetone (HPLC Grade),	150-001
MAG3 (Mertiatide™)	Acetone, (HPLC Grade), Chloroform 99.8% Tetrahydrofuran 99+%	150-951 & 150-952
MDP	Acetone (HPLC Grade), Distilled H ₂ O	150-001 & 150-005
Mebrofenin (Cholotec™)	150-160, Distilled H ₂ O	150-122 & 150-125
Pyrophosphate	Acetone (HPLC Grade), Distilled H ₂ O	150-001 & 150-005
Sestamibi (Cardiolite™ and Miraluma™)	Ethyl Acetate 99.9%	150-991
Sulphur Colloid	Acetone (HPLC Grade),	150-001
Tc-99m (reduced)	Acetone (HPLC Grade),	150-001
Tc-99m Monoclonal Antibodies	0.9% Saline	150-771
Tetrofosmin (Myoview™)	Ethyl Acetate 99.9%	150-971
In-111 Octreotide (Octreoscan™)	150-773 & 0.9% Saline	150-771 **
In-111 Monoclonal Antibodies	150-773 & 0.9% Saline	150-771 **
In-111/Y-90 Zevalin™	0.9% Saline	150-772
I-131 Monoclonal Antibodies And Neutrospec™	0.9% Saline	150-771

* See Sigma-Aldrich product matrix for solvents.

** Note: A Well Plate is suggested to perform QC on these radiopharmaceuticals.

Contact Fisher Scientific Corp. 800-766-7000 / www.fishersci.com / Item #14-245-71

Chromatography Strips

- 150-001** Strips, RED, 50/pkg
- 150-005** Strips, BLACK, 50/pkg
- 150-025** Strips, YELLOW, 50/pkg
- 150-122** Strips, ORANGE, 50/pkg
- 150-125** Strips, LIGHT BLUE, 50/pkg
- 150-130** Strips, GOLD, 50/pkg
- 150-771** Strips, DARK GREEN, 50/pkg
- 150-772** Strips, BLUE, 50/pkg
- 150-951** Strips, LIME, 50/pkg
- 150-952** Strips, PEACH, 50/pkg
- 150-971** Strips, TEAL, 50/pkg
- 150-991** Strips, PINK, 50/pkg
- 150-782** Strips, ALUMINUM, 50/pkg

Chromatography Solvents

- 150-160** Solvent, 20% Sodium Chloride, 30 ml
- 150-773** Solvent, DTPA, 5 ml
- 150-781** Solvent, Aluminum Standard, 5 ml, 10 µg/ml
- 150-783** Solvent, Aluminum Standard, 5 ml, 5 µg/ml

Related Items:

- 150-960** Developing Vials, 10 ml, 288/case
(used for all Tec-Control testing)
- 150-961** Developing Vials, 5 ml, 144/case
(used for Sestamibi & Tetrofosmin)
- 066-533** Forceps, Curved, locking, 9.5" l (24.1 cm)
- 066-535** Forceps, Straight, locking, 9.5" l (24.1 cm)
- 066-536** Forceps, Curved, non-locking, 12.5" l (31.7 cm)

Tec-Control Solvent Vendor:

Sigma-Aldrich Chemical Company
800-558-9160 / www.sigmaaldrich.com

Note: Customers outside the US should visit the Sigma-Aldrich web site to locate a regional office.

Solvent Description	Vendor Part #
Acetone HPLC Grade	27072-5
Ethyl Acetate 99.5% ACS Reagent	31990-2
Chloroform 99.8% ACS Reagent	31998-8
Tetrahydrofuran 99+% ACS Reagent	36058-9

PORTABLE SHIELDED ISOLATOR USP<797>

Certified Shielded Barrier Laminar Flow Isolator

- Small footprint
- Antechamber located under the Isolator
- Large, full-view shielded window
- Designed from the ground up for a Radiopharmacy
- Stainless Steel, inside and out
- Motorized height adjustable – 10 inch range
- Work sitting or standing with correct ergonomics
- USP<797> compliant overcomes expensive ventilation renovations
- Internal IV bar
- Sealed shielded dose calibrator chamber in the work surface
- Fits Atomlab 400 and 500 Dose Calibrators
- Use for blood work or other compounding or drawing requirements
- Mobile, on four heavy-duty locking casters
- Electric outlets in work area
- HEPA air filtration
- Lead shielded (.25" thick) for compounding sterile radiopharmaceuticals
- ISO Class 5 (Class 100) Isolator

The Cleanroom Solution

The Biodex Germfree Shielded Isolator provides the ultimate in both product and operator protection. The shielded isolator functions as a “glovebox” using negative pressure to meet NIOSH recommendations while complying with USP<797> regulations for use outside a cleanroom.

HEPA filtered unidirectional (laminar) air bathes the work area to protect the product from contamination and removes any particulate generated by sample manipulation. The operator is provided a complete shielded barrier from materials being manipulated in the work area. The Biodex Germfree Shielded Isolator meets or exceeds ISO 14644-1, ISO Class 5 (Class 100) air quality. Each shielded isolator undergoes rigorous physical testing to assure the unit meets performance requirements as validated. It is required that independent certification be performed before use.

What is a Barrier Laminar Flow Isolator?

In January 2004, USP<797> changed the way facilities that prepare compound sterile preparations (CSPs) approach their work. As the first official and enforceable requirement for CSPs, USP<797> deals with policies and practices for preparing CSPs. It applies to all facilities that prepare CSPs, including clinics, hospital care units, and main and satellite pharmacies. As facilities that prepare CSPs evaluate and change their procedures to meet USP<797> requirements, they are finding that barrier isolators provide an ideal alternative to more costly cleanrooms.

According to the Food and Drug Administration, a barrier isolator is “a decontaminated unit supplied with HEPA filtered air that provides uncompromised continuous isolation of its interior from the external environment, including surrounding cleanroom air and personnel.”

By installing a certified barrier isolator, one can obtain cleanroom conditions within a contained workspace. Barrier isolators provide



a ISO Class 5 (Class 100) environment for product preparation, with work occurring inside a closed, pressurized environment accessible only via sealed gloves that reach into the work area.

A Shielded Glovebox is not a certified barrier isolator. By choosing a certified barrier isolator, the requirement for a ISO Class 8 environment is eliminated. The Biodex Germfree Isolator is an excellent example of a well-designed shielded barrier Laminar Flow Isolator, meeting all standards and more.

Antechamber

- Sealed two-door air lock maintains complete environmental separation between work area and ambient conditions.
- HEPA filtered purge of trapped air lock air eliminates cross contamination between the work area and the room during both material ingress and egress.

Filtration


- Full framed, standard size mini-pleat HEPA's filter 100% of inlet and exhaust air from both the work area and the antechamber to provide a fully contained environment.
- HEPA filters are full coverage and front loading for easy replacement by a certifier.

Ergonomics

- Hydraulic assist height adjusting stand offers a full 10" range, allowing operators to sit or stand comfortably for extended time periods.
- Stainless steel sliding tray inside the antechamber air lock that pulls forward for loading and unloading items into the air lock, eliminating reaching strains. The sliding tray lifts up 7.5" when the interior antechamber door is opened.



Rapid Exchange Glove System



The GRx offers a fast, easy, efficient solution to the traditional glove changing methods. Gloves can now be changed in seconds, maintaining a clean environment.

- Two part sleeve/glove system allows the use of most types and sizes of commercially available gloves for better dexterity and tactility.
- Extra large, oval glove ports are placed with bottoms together to provide an anthropometrically correct configuration that accommodates a wide range of body types and increases range of movement.
- Glove ports have a 3" armrest to reduce operator fatigue.
- Large viewing panel is at an angle to reduce glare and operator strain.

Sharps Container

- One sharps container in a locking shielded holder, located under the isolator in a .25" lead shield.
- Sharps tube is a straight 2.5" diameter tube to facilitate quick drop of the largest syringes with an easily removable shielded seal/stopper to maintain work area containment.

Safety

- Lock-out handle requires key for access to work area and sharps container.
- Digital pressure readout with low pressure alarm for work area.
- Inward face velocity is 95-100 LFM at glove port opening, to protect operator/product during massive breach of containment.
- Glove changes can be made without breaking containment.
- Locking casters are standard, seismic anchors are available.

Cleaning

- Front viewing panel is top hinged and self supporting for easy access to the work area for cleaning and equipment ingress and egress.
- All filter diffuser / guard panels are removable for easy cleaning.
- Straight sides and back maximize work area to accommodate the many types and shapes of equipment and dose calibrators.
- All corners in work area, antechamber and work surface are easily reached and cleaned.

Configuration Options

- Overall work area pressure is negative with the option to externally exhaust.
- Optional Thimble Exhaust Transition Kit is available for connection to external exhaust systems
- Fits Atomlab and other dose calibrators

SPECIFICATIONS:

Dimensions:

Overall Dimensions: 36" w x 32.5" depth x 79" h (91.4 x 82.5 x 200.7 cm)

Designed to fit through standard door openings and elevators.

Height Adjustment: 79" – 89" (200.7 x 226 cm)

Work Area Dimensions: 35" w x 25" depth x 28" h (88.9 x 63.5 x 71 cm)

Viewing window: 30" w x 23" h (76.2 x 58.4 cm), 1.8"

(46 mm) thick leaded acrylic is equal to 2 mm lead equivalency

Construction: Stainless steel with a pharmaceutical grade finish

Lead Shielding: .25" thick (.64 cm) on back, sides, bottom and front

Fluorescent Lights: High efficiency; externally mounted to minimize heat build up. Separate lighted power ON/OFF switch.

Motor/Blower: High capacity with speed control to extend HEPA filters life. Separate lighted power ON/OFF switch.

Power: 115 VAC, 60 Hz (Optional 220 VAC, 50-60 Hz is available), 15 AMP line, running amperage = 6 AMP. Includes ten foot hospital grade power cord with molded grounded plug. Sealed outlet in work area.

Weight: 1600 lb (725.7 kg)

Warranty: Two years parts; one year labor

190-215 Shielded Isolator

Related:

214-129 Chair, Ergonomic

190-217 GRX – Rapid Exchange Glove System

Includes one pair of isolator sleeves, two outer rings, six inner rings and all necessary O-rings

190-219 Total Customer Care Plan

Includes installation and training (see web for details)

190-220 Thimble Exhaust Transition Kit

For connection to external exhaust systems

Replacement:

039-413 Sharps Container, 3.2 qt., 30/pkg

190-218 Hypalon Sleeves (pr)


For details on the Total Customer Care Plan, visit us on the web at www.biodec.com/isolator.

ATOMLAB™ 400 DOSE CALIBRATOR

Designed for facilities receiving unit doses including PET and BETA



The Atomlab™ 400 provides fast, accurate radionuclide activity measurements with performance that easily complies with the most stringent regulatory requirements.



- *Pre-programmed for 88 most commonly used radionuclides*
- *Large, easy-to-read backlit LCD*
- *Small footprint economizes workspace*
- *Ultra-fast response*
- *Automatic range selection; ranges up to 40 Curies of Tc-99m or 10 Curies of F-18*
- *Displays in Curies or Becquerels*
- *Remote Ionization Chamber*
- *Self-Diagnostic Software*
- *Desktop or wall mount display*
- *Two-year warranty*
- *RS-232 bi-directional serial communications port*

The Atomlab™ 400 provides fast, accurate radionuclide activity measurements with performance that easily complies with the most stringent regulatory requirements.

The unit is simple to operate. There is a routine list of ten pre-programmed isotopes plus another seven that are user selected from the library. The library contains 88 isotopes listed alphabetically, including Y-90 and Sr-89. Four isotopes are displayed at a time.

Activity is displayed on a LCD panel in either Curie or Becquerel units. Background correction is performed at the touch of a button. Range selection is automatic.

Activity measurements are performed by a microprocessor-controlled electrometer located within the detector assembly of the ionization chamber. The chamber is shielded with .25" (6.3 mm) lead. It can be located up to eight feet away from the display unit. Chamber bias is generated within the display unit by an electronic high voltage supply, eliminating the need for expensive battery changes.

The RS-232 port enables the Atomlab 400 Dose Calibrator to communicate with most commercially available nuclear medicine management systems.

SPECIFICATIONS:

Isotope Selection Keys: Ten pre-programmed – Tc-99m, Tl-201, Co-57, Cs-137, I-131, In-111, Ga-67, Xe-133, I-123, and Mo-99; seven additional keys for user-set isotopes; two new isotope keys and a full alphabetical list of 88 isotopes.

Activity Range: 0.01 μ Ci to 40 Ci (.0004 MBq to 1500 GBq) of Tc-99m or 10 Ci of F-18.

Energy Range: 25 keV to 3 MeV photons

Response Time: One to two seconds for doses greater than 200 uCi; three seconds for doses greater than 20 uCi; 50-100 seconds below 20 uCi of Tc-99m with default threshold, threshold adjustable to reduce counting time

Detector Linearity: $\pm 1\%$ or 0.2 μ Ci, whichever is greater

Electrometer Linearity: $\pm 1\%$ or 0.2 μ Ci, whichever is greater

Electrometer Accuracy: $\pm 1\%$ or 0.2 μ Ci, whichever is greater

Overall Accuracy: $\pm 3\%$ or 0.3 μ Ci, whichever is greater; overall accuracy is affected by such factors as the accuracy of the specific source calibration, geometric variations due to sample volume or configuration, detector linearity, electrometer accuracy and readout accuracy

Repeatability: $\pm 0.3\%$ above 1 mCi short term (24 hr); 1% long term (one yr); exclusive of background

Digital Calibration Dial: Four-digit LED dial display with increment/decrement keys to change the value; range is from 0.0 to 999.9

Detector: Well-type pressurized ionization chamber, with Argon fill gas; well opening 2.75" (7 cm), well depth 10.25" (26 cm)

Chamber Gas Pressure: 149KPa gauge (21.6 psig) at 20 degrees C or 250KPa absolute (36.3 psia) at 20 degrees C. IATA regulation 3.2.2.4 Exempts Gases of Division 2.2 from Dangerous Goods Regulations when transported at pressure less than 200KPa gauge (29 psig) at 20 degrees C. Device is shipped standard goods.

Detector Shielding: .25" (6.3 mm) lead on all sides except top well opening; supplementary shielding available

Chamber Bias: 355 \pm 5 volts

Environmental Operating Conditions:

Temperature: 0-40° C

Humidity: 0-90% rh, non-condensing

Power Requirements: 100 to 240 VAC, 0.6 – 0.3 amps, auto switching; APS Power Supply (APS22ES-150160), for medical use.

Line Frequency: 50/60 Hz; detachable line cord; built-in EMI filter and transient suppression

Detector and Interface Cables: 8' (243 cm) long, six conductor cables (two carry power, two for chassis ground, two carry serial data for digital I/O)

Display Unit:

Dimensions: 6.75" w x 6" depth x 5" h (17.1 x 15.3 x 12.7 cm)

Weight: 3.6 lb (1.64 kg); desktop or wall mountable

Detector Unit:

Dimensions: 6" dia x 15.5" h (15.24 x 39.37 cm)

Well I.D.: 2.75" dia x 10.5" h (7 x 26.7 cm)

Well I.D. with Liner: 2.5" dia x 10.25" h (6.35 x 26 cm)

Lead Shielding: .25" thick (6.3 mm)

Weight: 35 lb (16 kg)

Approvals: ETL to UL 60601-1 and cETL to CAN/CSA C22.2 No. 601-1-M90

IEC 60601-1, IEC 60601-1-4 and IEC 60601-1-2 and CE marked

Warranty: Two-year



For physics tests, sample reports, testimonials and warranty information, visit us on the web
www.biodex.com/dosecalibrators

New low pressure chamber. Atomlab Dose Calibrators are shipped air or ground as standard goods.

Which dose calibrator is right for you? Visit www.biodex.com/dosecalibrators for comparison chart.



The Atomlab™ 400 display can be mounted on a wall or placed on a desktop.

086-335 Dose Calibrator, Atomlab™ 400,
100-240 VAC
*Includes: Smart Display, ionization chamber,
RS-232 port, vial/syringe dipper and well insert.*

Related:

086-338 Dose Calibrator Shielding Rings,
Interlocking, 2.25" lead
*For additional protection from
high energy activity*

086-423 Moly Shield, Vial, .3" lead

086-435 Moly Shield, Syringe, .3" lead

086-509 Lineator

086-334 Cable, European to Wall Outlet

Replacement:

086-242 Vial/Syringe Dipper

086-241 Well Insert

An industry exclusive two-year warranty is standard.

Contact Biodex to schedule a one-on-one online demonstration
 Call 1-800-224-6339, or Email: aalessi@biodex.com

Int'l 631-924-9000 • www.biodex.com

ATOMLAB™ 500 DOSE CALIBRATOR

One dose calibrator that can be used for a wide variety of nuclear medicine, PET and radioimmunotherapy applications.



Proven performance for fast, accurate measurements. One dose calibrator for all your requirements.



- Easy to use, large color touch screen display with intuitive menus
- Automatic range selection; ranges up to 100 Curies of Tc-99m or 25 Curies of F-18
- Pre-programmed for 88 most commonly used radionuclides; any 12 can be conveniently touch selected
- Displays in Curies or Becquerels
- Small footprint economizes workspace
- Ultra-fast response
- Robust software and extensive functionality
- Report and label printers available
- The USB ports also allow software upgrades via USB memory devices.
- Desktop or wall mount display
- Communicates with most commercially available NM management systems via bi-directional serial communications port
- USB ports offer the ability to accommodate a USB mouse and printing devices.
- Two-year warranty



Standard Apps for Atomlab 500 include:

- Automated Quality Assurance Apps
 - Constancy and Expanded Constancy
 - Linearity and Auto Linearity
 - Accuracy
 - Geometry
- Nuclear Pharmacy Apps
 - Future dose computation
 - Volume determination
 - Inventory control of 25 samples, correcting volume, activity and moly concentration
- Multiple Detector Apps
 - Manages multiple ionization chambers
- Wipe Test Counter Apps
 - Upgradeable at any time to include a wipe test counter

The Atomlab 500 provides fast, accurate radionuclide activity measurements with performance that easily complies with the most stringent regulatory requirements. The system consists of a new low pressure ionization chamber with redesigned seal, electrometer with extraordinary linearity and an auto-ranging touch screen color display. Now one dose calibrator can be used for a wide variety of nuclear medicine, PET and radioimmunotherapy applications. Additionally, there are advanced, but easy-to-use programs for nuclear pharmacy, radiochemistry and radiochromatography.

Activity measurements are performed by the microprocessor controlled electrometer located within the chamber assembly. The chamber is shielded with .25" (6.3 mm) lead. It can be located up to eight feet away from the display unit. Chamber bias is generated by an electronic high voltage supply, eliminating the need for expensive battery changes.

Every element of the design and technical development will increase dose accuracy, department productivity and regulation compliance. The attractive and intuitive human interface guides the user through each operation. Software can easily be updated via the Biodex website or by using a convenient memory card. The touch-screen display can rest on a bench or mount on the wall of a hot lab, hot cell or laminar flow hood.


In addition to powerful self diagnostics, the Atomlab 500 includes an exclusive chamber monitoring technology to assure longer life and accuracy. Integrated pressure and temperature sensors feedback data so that the influence of gas pressure change will not effect an accurate reading.

Operation

The system is easy to use. There are 12 isotope selection touch keys pre-programmed for the most commonly used radionuclides. Any of those keys can be reprogrammed by the user for a desired isotope. There are 88 isotope-specific dial values listed in the library. Dial values can easily be changed if required.

Activity is displayed on the touch screen color display in either Curie or Becquerel units. Background correction is performed at the touch of a button. Range selection is automatic, from .01 microcurie to 100 Curies of Tc-99m or 25 Curies of F-18.


Accuracy



Atomlab Dose Calibrators have consistently proven to be highly accurate. Biodex and chamber manufacturer Sun Nuclear Corporation have participated in the isotope program sponsored by National Institute of Standards and

Technology. Each month a certified isotope is received from National Institute of Standards and Technology and is measured in the Atomlab™ Dose Calibrator, producing direct traceability to National Institute of Standards and Technology.


Quality Assurance




The Atomlab 500 has been designed to make life easier. The extensive selection of quality assurance applications streamlines and simplifies hot lab administration requirements. The system stores and decay corrects multiple reference sources and compares the measured activity to the calculated activity for the daily constancy test.

Linearity tests can be performed in the traditional manual method or by a fully automated program that allows for readings from a source to be taken, and automatically recorded at specified intervals. The system will graph the results.


The attenuation tube test for linearity can be performed using software that will guide the user through the procedure, store the values and make all calculations.



Constancy Activity




Linearity Decay



Geometry Test

Data Download



The Atomlab™ Data Manager is available as an option. The Windows™ based utility allows wipe test and dose calibrator QA results to be downloaded using a USB/serial converter. The results can be viewed and printed from the data manager software as required. In addition, information stored in the data manager can be exported into Microsoft® Excel or to department management systems.

ATOMLAB™ 500 DOSE CALIBRATOR

An extensive selection of quality assurance applications streamlines and simplifies hot lab administration requirements



The Atomlab™ 500 display can be mounted on a wall or placed on a desktop.

Commercial Nuclear Pharmacy

The Atomlab 500 Dose Calibrator features a nuclear pharmacy “Dose Calculation Screen” to meet the needs of a commercial nuclear pharmacy. All the information needed to draw doses efficiently is readily available on one screen. The feature easily performs pre- and post-decay calculations, volume calculations for specific times and isotope changes all with minimal screen touches. No calculators needed! Atomlab Dose Calibrators can read up to 100 Curies of Tc-99m, thereby eliminating the need for an aliquot preparation, a significant time saver every time a generator is milked.

Industrial Nuclear Pharmacy

The Atomlab 500 provides inventory control for 25 samples, storing and correcting the volume, activity, and moly concentration. The system will perform both volume and future dose calculations. In addition to inventory management, the Atomlab 500 provides quality assurance and record keeping functions. The inkjet printer allows hard copy records to be produced for all functions. The label printer allows the user to print labels for the syringe or vial.

Radiochromatography

The radiopharmaceutical quality control program is exceptional. The Atomlab 500 performs all counting and calculations for paper chromatography tests, computing the percentages of free pertechnetate, hydrolyzed reduced Tc-99m and labeled radiopharmaceuticals.

Radiochemistry

Up to seven ionization chambers, or six ionization chambers and a well counter can be connected via daisy chain to a single display. The activity in each detector can be selected and viewed from the single display.

Communications

RS-232 port and two USB ports to communicate in real time with the most commercially available nuclear medicine management systems, connect to external monitor or upload software upgrades.

Sample Labels

Future Dose

Future Dose Record	
Patient: _____	Procedure: _____
Id #: _____	Kit: _____
Prepared By: _____	Lot #: _____
Isotope: Tc-99m	
Dose Ordered: 18.00 mCi at 10:45 EST 06/01/2010	
Calc Current Activity: 18.67 mCi at 10:28 06/01/2010	
Administered By: _____	

CAUTION RADIOACTIVE MATERIAL	Radiopharmaceutical: Tc-99m
	Dose Ordered: 18.00 mCi at 10:45 EST 06/01/2010
	Calc Current Activity: 18.67 mCi at 10:28 06/01/2010

CAUTION RADIOACTIVE MATERIAL	Radiopharmaceutical: Tc-99m
	Dose Ordered: 18.00 mCi at 10:45 EST 06/01/2010
	Calc Current Activity: 18.67 mCi at 10:28 06/01/2010

Future Dose with Dose in Chamber

Future Dose Record	
Patient: _____	Procedure: _____
Id #: _____	Kit: _____
Prepared By: _____	Lot #: _____
Isotope: Tc-99m	
Calc Admin Dose: 106.3 mCi at 11:30 EST 06/01/2010	
Measured Activity: 119. µCi at 10:33 06/01/2010	
Administered By: _____	

CAUTION RADIOACTIVE MATERIAL	Radiopharmaceutical: Tc-99m
	Calc Admin Dose: 106.3 mCi at 11:30 EST 06/01/2010
	Measured Activity: 119. µCi at 10:33 06/01/2010

CAUTION RADIOACTIVE MATERIAL	Radiopharmaceutical: Tc-99m
	Calc Admin Dose: 106.3 mCi at 11:30 EST 06/01/2010
	Measured Activity: 119. µCi at 10:33 06/01/2010

SPECIFICATIONS:

Display: LCD Touch Panel 6.5" x 5", function keys are displayed for the operation being performed

Connectors: RJ-12 for well cable

Power: This system uses XP Power Supply for Medical Use, Model #PDM60US15

Line Voltage: 100 to 240 VAC, auto selectable by the power supply, 1.5-0.75 amps

Line Frequency: 50/60 Hz, detachable line cord, built-in EMI filter and transient suppression

Auxiliary Port: Two USB ports, one RS-232 port

Memory: Stores: Inventory and QA tests

Isotope Selection Keys: Twelve pre-programmed – Tc-99m, Tl-201, I-123, I-131, Cs-137, Co-57, Xe-133, Ga-67, In-111, F-18, Y-90, Mo-99; 25 user-defined isotopes and a full alphabetical list of 88 isotopes.

Activity Range: : 0.01 uCi to 100 Ci (.0004 MBq to 3700 GBq) of Tc-99m or 25 Ci of F-18

Energy Range: 25 keV to 3 MeV photons

Response Time: One to two seconds for doses greater than 200 uCi; three seconds for doses greater than 20 uCi; 50-100 seconds below 20 uCi of Tc-99m with default threshold; threshold adjustable to reduce counting time

Detector Linearity: $\pm 1\%$ or 0.2 μCi , whichever is greater

Electrometer Linearity: $\pm 1\%$ or 0.2 μCi , whichever is greater, up to 40 curies of Tc-99m, $\pm 1.5\%$ up to 100 curies of Tc-99m

Electrometer Accuracy: $\pm 1\%$ or 0.2 μCi , whichever is greater

Overall Accuracy: $\pm 3\%$ or 0.3 μCi , whichever is greater; overall accuracy is affected by such factors as the accuracy of the specific source calibration, geometric variations due to sample volume or configuration, detector linearity, electrometer accuracy and readout accuracy

Repeatability: $\pm 0.3\%$ above 1 mCi short term (24 hr); 1% long term (one yr)

Digital Calibration Dial: Four-digit dial with increment/decrement keys to change the value; range is from 0.0 to 999.9

Detector: Well-type pressurized ionization chamber, with Argon fill gas; well opening 2.75" (7 cm), well depth 10.25" (26 cm). Up to seven chambers can be serially connected to one display.

Chamber Gas Pressure: 149KPa gauge (21.6 psig) at 20 degrees C or 250KPa absolute (36.3 psia) at 20 degrees C. IATA regulation 3.2.2.4 Exempts Gases of Division 2.2 from Dangerous Goods Regulations when transported at pressure less than 200KPa gauge (29 psig) at 20 degrees C. Device is shipped standard goods.

Detector Shielding: .25" (6.3 mm) lead on all sides except top well opening; supplementary shielding available

Chamber Bias: 355 ± 5 volts

Environmental Operating Conditions: Temperature: 0-40° C; Humidity: 0-90% rH, non-condensing

Power Requirements: 100 to 240 VAC, 0.38 – 0.15 amps, auto switching; XP Power Supply (PDM60US15), for medical use.

Line Frequency: 50/60 Hz; detachable line cord; built-in EMI filter and transient suppression

Detector and Interface Cables: 8' (243 cm) long, six conductor cables (two carry power; two chassis ground; two carry serial data for digital I/O)

Display Unit:

Dimensions: 9.5" w x 12" depth x 12" h (24.1 x 30.5 x 30.5 cm)

Weight: 6.3 lb (2.9 kg); desktop or wall mountable

Detector Unit:

Dimensions: 6" dia x 15.5" h (15.24 x 39.37 cm)

Well I.D.: 2.75" dia x 10.5" h (7 x 26.7 cm)

Well I.D. with Liner: 2.5" dia x 10.25" h (6.35 x 26 cm)

Lead Shielding: .25" lead (6.3 mm)

Weight: 35 lb (16 kg)

Approvals: ETL to UL 60601-1 and cETL to CAN/CSA C22.2 No. 601-1-M90 IEC 60601-1, IEC 60601-1-4 and IEC 60601-1-2 and CE marked

Warranty: Two-year



For physics tests, sample reports, testimonials and warranty information, visit us on the web
www.biodex.com/dosecalibrators

New low pressure chamber. Atomlab Dose Calibrators are shipped air or ground as standard goods.

Which dose calibrator is right for you? Visit www.biodex.com/dosecalibrators for comparison chart.

Contact Biodex to schedule a one-on-one online demonstration
 Call 1-800-224-6339, or Email: aalessi@biodex.com



086-341, Optional label printer for Atomlab 500 and Atomlab 500 Plus

086-330 Dose Calibrator, Atomlab™ 500, 100-240 VAC
Includes: Smart Display, ionization chamber, RS-232 port, vial/syringe dipper and well insert.

Related:

086-333 Software, Atomlab™ 500 Data Manager

086-336 Chamber, Dose Calibrator

075-594 Chamber, Wipe Test

086-338 Shielding Rings, Interlocking, 2.25" lead
For additional protection from high energy activity

086-509 Lineator

086-435 Moly Shield, Syringe, .3" lead

086-243 Copper Dipper, Vial/Syringe

086-423 Moly Shield, Vial, .3" lead

086-341 Printer, Label, Dymo
Includes one roll 086-343 labels; two rolls 086-344 label

086-339 Printer, Ink Jet (report)

086-334 Cable, European to Wall Outlet

Replacement:

086-242 Vial/Syringe Dipper

086-241 Well Insert

086-343 Label, Blank, Lg, 300/roll
(For Dymo Printer, 086-341)

086-344 Label, "Radioactive", Sm, 260/roll, 2/pkg
(For Dymo Printer, 086-341)

An industry exclusive two-year warranty is standard.

EASY UPDATE

Your Atomlab 500 Dose Calibrator is upgradable. You can easily install software updates via the Biodex website or by using a convenient USB memory device.

ATOMLAB™ WIPE TEST COUNTER

Eliminate the tedium of wipe testing with the Atomlab™ Wipe Test Counter

Finally, a wipe test counter that's easy to use, easy to understand and dependable.

UPGRADE-ABILITY

The Atomlab™ Wipe Test Counter can be upgraded to include a fully functioning Atomlab™ Dose Calibrator simply by adding a new chamber. The "Smart Display" detects which chambers are connected and instantly reconfigures the screen to the appropriate icons.



- Easy-to-use, large touch screen display with intuitive menus
- 64 Channel MCA
- Adjustable wide window and individual isotopes
- 2x2 NaI drilled-well detector
- Remote shielded well
- Energy spectrums with individual ROI
- Ability to help identify isotopes causing contamination
- User-specific wipe locations and trigger levels
- Wipes that exceed trigger levels are immediately recognized: flashes on screen and prints in red
- Detailed wipe reports including cpm and dpm
- Wipe testing results stored
- Upgradable at any time to a dose calibrator by adding an ionization chamber
- Optional printer for hard-copy archives
- Two-year warranty
- Communicates with most commercially available NM management systems via bi-directional serial communications port



Standard Apps for Atomlab Wipe Test Counter include:

- Automated Quality Assurance Apps
 - Full Width at Half Max (FWHM)
 - Chi Square
 - Minimum Detectable Activity (MDA)

Finally, a wipe test counter that's easy to use, easy to understand, fast and dependable. A color touch-screen display utilizing intuitive software eliminates the tedium of wipe testing. Simply perform a daily calibration and background count, then count the wipe for each predetermined location. Trigger levels can be set for any isotope at any location including 200 dpm for iodine. In seconds the system will determine if the location is above or below the user defined trigger level.

When performing a wipe test, the full spectrum is displayed. A wide window that includes the isotope energies expected in a particular department is set by the user. The efficiencies of the isotopes selected for the window can be either factory defaults or user determined using an integrated detector efficiency program. Individual isotope ROIs along with the wide window can be set. This feature helps identify the isotope(s) causing contamination.

Up to 50 wipe locations can be entered as a restricted area, unrestricted area, sealed source or package. The results are displayed in dpm, cpm, μCi or kBq.

The system consists of a lead shielded 2" x 2" sodium iodide (NaI) well detector and a 64 channel multi-channel analyzer. The displayed energy range (spectrum) is 0-800 KeV, which is typically found in nuclear medicine departments.

The wipe counter is designed to meet or exceed all NRC (10 CFR 35.70, 10 CFR 20.1906 and 10 CFR 35.2067) and state wipe test requirements. There are automated programs for the quality assurance functions: calibration, FWHM, chi-square and minimum detectable activity (MDA).

Wipe test results and QA test data can be stored in memory and printed at any time.

SPECIFICATIONS:**DISPLAY**

Display: LCD Touch Panel 6.5" x 5", function keys are displayed for the operation being performed

Connectors: RJ-12 for well cable

Power: This system uses XP Power Supply for medical use, Model #PDM60US15

Line Voltage: 100 to 240 VAC, auto selectable by the power supply, 1.5 – 0.75 amps

Line Frequency: 50/60 Hz, detachable line cord, built-in EMI filter and transient suppression

Auxiliary Port: Two USB ports, one RS-232 port

Memory: Stores wipe, calibration, background, high voltage, isotope specification, isotope efficiency, Chi-Square testing results, technologist list, wipe locations list and latest MDA calculation. Results can be displayed and printed.

Preset Radionuclides: 27 including Tc-99m, Co-57, Cs-137, Ga-67, Tl-201, I-123, I-125, I-131, In-111, F-18

DETECTOR CABLE

Length: Standard 8 ft. (243 cm); custom cable lengths available up to 20 meters

Conductors: Six total – 2 for power, two for chassis ground, two for serial data

Connectors: RJ-12

WELL COUNTER

Detector: 2" x 2" NaI (Tl) integral line scintillation detector with a 0.75" dia x 1.44" depth well (1.9 x 3.7 cm)

Style: Remote Detector

Channels: 64

MCA: Integral to Well Counter

Spectral Resolution: FWHM 10%

Count Rate: (Maximum) 30,000 cps

Connectors: RJ-12 for power and signal

Lead Shielding: 0.5" (1.2 cm) integral lead shield

Optional Calibration Source: 0.1 µCi Cs-137 Rod Source

Optional Lead Shield: 0.5" (1.2 cm) lead shield, slides over integral shield for total of 1.0" of shielding – includes cover

Cover: Fits optional additional lead shield, has 0.25" (.6 cm) lead shielding

PHYSICAL DATA

Display Unit:

Size: 9.5" w x 12" h x 12" d (24.1 x 30.5 x 30.5 cm)

Weight: 4.2 lb (1.9 kg)

Well Counter:

Size: 6" d x 11" h (15.24 x 27.9 cm)

Weight: 29 lb (3.2 kg)

Approvals: ETL listed to UL 60601-1 and CAN/CSA C22.2 No. 601-1M90, IEC 60601-1, IEC 60601-1-4 and IEC 60601-1-2 and CE marked

Warranty: Two-year



The Atomlab™ Wipe Test Counter display can be mounted on a wall or placed on a desktop. Chamber shown with optional Lead Shield.

086-331 Atomlab™ Wipe Test Counter

Includes: Smart Display, well counter and RS-232 port.

Related:

086-333 Software, Atomlab 500 Data Manager

086-336 Chamber, Dose Calibrator

075-596 Lead Shield, Wipe Test Chamber

086-339 Printer, Ink Jet (report)

086-334 Cable, European to Wall Outlet

063-139 Rod Source, Cs-137, Calibrated, 0.1 µCi

006-350 Wipe Test Kits, 500/pk

086-342 Well Liners, Disposable, 100/pk


An industry exclusive two-year warranty is standard.

EASY UPDATE

Your Atomlab Wipe Test Counter is upgradable.

You can easily install software updates via the Biodex website or by using a convenient USB memory device.

Sample reports available at
www.biodex.com/wipetestcounter



Printed reports

Contact Biodex to schedule a one-on-one online demonstration
Call 1-800-224-6339, or Email: aalessi@biodex.com

Int'l 631-924-9000 • www.biodex.com

ATOMLAB™ 500Plus DOSE CALIBRATOR

Dose Calibration and Wipe Testing combined...it's a complete Hot Lab Management System.



*Atomlab™ 500Plus
Dose Calibrator and
Wipe Test Counter
Finally...one solution
for all your molecular
imaging needs.*



The Atomlab 500Plus Dose Calibrator and Wipe Test Counter have brought it all together - science, technology and application.

- Easy-to-use, large, color touch-screen display with intuitive menus
- Automatic range selection; ranges up to 100 Curies of Tc-99m or 25 Curies of F-18
- Pre-programmed for 88 most commonly used radionuclides; any 12 can be conveniently touch selected
- Displays in Curies or Becquerels
- Small footprint economizes workspace
- Ultra-fast response
- Robust software and extensive functionality
- Report and label printers available
- USB ports offer the ability to accommodate a USB mouse and printing devices.
- The USB ports also allow software upgrades via USB memory devices.
- Communicates with most commercially available NM management systems via bi-directional serial communications port
- 64 Channel MCA
- Adjustable wide window and individual isotopes

- 2x2 NaI drilled-well detector
- Remote shielded well
- Energy spectrums with individual ROI
- Ability to help identify isotopes causing contamination
- User-specific wipe locations and trigger levels
- Wipes that exceed trigger levels are immediately recognized: flashes on screen and prints in red
- Detailed wipe reports including cpm and dpm
- Wipe testing results stored
- Desktop or wall mount display
- Two-year warranty



Standard Apps for Atomlab 500Plus include:

- Unit offers all the Apps of both the Atomlab 500 and Wipe Test Counter

Combine the Atomlab 500 Dose Calibrator with the Atomlab Wipe Test Counter and create a complete, efficient and cost effective radioactivity measurement system... the Atomlab 500Plus.

Consider the software that is supplied with the Atomlab 500Plus. It's comprehensive, easy-to-use and feature rich. The software mirrors the way you think and work. It guides when necessary, but does not burden the advanced user. The touch-screen and easy-to-follow on-screen prompts mean you won't need "cheat sheets."

When required, you can be confident that every operation is captured and documented. That documentation makes compliance a breeze.

Your Atomlab 500 Dose Calibrator is upgradable. You can easily install software updates via the Biodex website or by using a convenient memory card.

Finally, one solution for all your molecular imaging needs.


SPECIFICATIONS:

Atomlab™ 500 Dose Calibrator


See page 69 for complete listing.

Atomlab™ Wipe Test Counter


See page 71 for complete listing.



Wipe Test Passed



Wipe Test Failed



QA Display
(Dose Calibrator)



The Atomlab™ 500Plus display can be mounted on a wall or placed on a desktop.

086-332 Dose Calibrator, Atomlab™ 500Plus,
100-240 VAC

Includes: Smart Display, ionization chamber,
well counter, RS-232 port, vial/syringe dipper
and well insert.

Related: Dose Calibrator

086-333 Software, Atomlab 500 Data Manager

086-336 Chamber, Dose Calibrator

075-594 Chamber, Wipe Test

086-338 Shielding Rings, Interlocking, 2.25" lead
For additional protection from
high energy activity

086-509 Lineator

086-435 Moly Shield, Syringe, .3" lead

086-243 Copper Dipper, Vial/Syringe

086-423 Moly Shield, Vial, .3" lead

086-341 Printer, Label, Dymo

Includes one roll 086-343 labels; two rolls
086-344 label

086-339 Printer, Ink Jet (report)

086-334 Cable, European to Wall Outlet

Replacement:

086-242 Vial/Syringe Dipper

086-241 Well Insert

086-343 Label, Blank, Lg, 300/roll
(For Dymo Printer, 086-341)

086-344 Label, "Radioactive", Sm, 260/roll, 2/pkg
(For Dymo Printer, 086-341)

Related: Wipe Test Counter

063-139 Rod Source, Cs-137, Calibrated, 0.1 µCi

075-596 Lead Shield, Wipe Test Chamber

006-350 Wipe Test Kits, 500/pkg

086-342 Well Liners, Disposable, 100/pk

An industry exclusive two-year warranty is standard.

New low pressure chamber. Atomlab Dose Calibrators are shipped air or ground as standard goods.


Which dose calibrator is right for you? Visit www.biodex.com/dosecalibrators for comparison chart.

Contact Biodex to schedule a one-on-one online demonstration
Call 1-800-224-6339, or Email: aalessi@biodex.com

Int'l 631-924-9000 • www.biodex.com

LINEATOR

Simulates eight different source strengths — using only one source



- *Simplifies compliance with NRC and state requirements*
- *Checks linearity without sample decay or fractioning*
- *Covers diagnostic or therapeutic quantities in a single pass*

The Lineator is a simple device to accurately and reliably verify the linearity of your dose calibrator. Test results are available in minutes, without waiting days for decay, making it feasible to perform a linearity test more often. Early identification can prevent problems before they occur.

The Lineator is a set of five tubes. Four tubes are used to perform the test and an interchangeable tube is used depending on the dynamic range needed to cover.

To perform a linearity test, insert a source of Tc-99m of the maximum activity to be measured into the central tube, then place the tube into the chamber of your dose calibrator and count. The remaining lead-lined tubes are placed, one at a time, concentrically over the central tube and counted individually or in combination. The readings are then normalized with predetermined factors, and the degree of linearity can be seen virtually at a glance.


Using only one source, the Lineator can simulate up to eight different source strengths. Each outer tube absorbs a portion of the source radiation and reduces the effective source activity seen by the dose calibrator.

SPECIFICATIONS:

Weight: 6 lb (3 kg)

086-509 Lineator

MOLY ASSAY SHIELD



Model 086-435 shown

SPECIFICATIONS:

086-435 Moly Assay Shield, Syringe

Dimensions: 9.125" h x 2.125" dia (23.2 x 5.4 cm)

Lead Shielding: .3" thick (7.6 mm)

086-423 Moly Assay Shield, Vial

Dimensions: 3.5" h x 2" dia (8.9 x 5.1 cm)

Lead Shielding: .3" thick (7.6 mm)

086-435 Moly Assay Shield, Syringe, .3" lead

086-423 Moly Assay Shield, Vial, .3" lead

VIAL/SYRINGE DIPPER




This rugged, Vial/Syringe Dipper has a comfortable handle and it will hold 1 cc to 10 cc syringes or up to a 30 ml vial. Included with every Atomlab Dose Calibrator, the Vial/Syringe Dipper will also fit any well chamber with 2.5" x 10" interior dimension. The material used will not cause attenuation problems and is very resistant to breaking.

086-242 Dipper, Vial/Syringe

SYRINGE DIPPER

Accommodates B-D and Monoject 60 cc syringes



The BEXXAR Syringe Dipper will hold a 60 cc syringe in the correct position in most commercial dose calibrators with 2.5" x 10" well chamber interior. The material used will not cause attenuation problems and is very resistant to breaking.

008-015 Dipper, Syringe, 60 cc

COPPER SYRINGE DIPPER



- *For use with I-123 and In-111*
- *Eliminate variations*

Designed for use with I-123 and In-111, the Copper Dipper removes variation in readings caused by attenuation differences from different materials and thicknesses used in syringes and vials.

SPECIFICATIONS:

Dimensions: 10.25" h x 1.63" dia
(26 x 4.1 cm)
Weight: .75 lb (.34 kg)

086-243 Copper Dipper, Vial/Syringe

DOSE CALIBRATOR SHIELDING RINGS

Especially suited for PET applications



086-338 Dose Calibrator Shielding Rings are 2.25" lead.

The Dose Calibrator Shielding Rings offer an additional 2.25" of lead shielding around the remote chamber for working with 511 keV radionuclides, such as FDG F-18.

SPECIFICATIONS:

086-338 Dose Calibrator Shielding Rings
Dimensions: 10.6" dia x 14.9" h
(27 x 37.9 cm)

I.D.: 6.1" dia (15.5 cm)
Lead Shielding: 2.25" thick (5.7 cm)
Weight: 359 lb (162.5 kg)
Shipping Weight: 388 lb (175.7 kg)

086-301 Dose Calibrator Shielding Rings
Dimensions: 11.625" l x 11.625" w x 16.4" h (29.5 x 29.5 x 41.6 cm)

I.D.: 11.25" dia x 16" h (28.6 x 40.6 cm)
Lead Shielding: 2.25" thick (5.7 cm)
Weight: 455 lb (206 kg)
Shipping Weight: 484 lb (219.5 kg)

086-338 Shielding Rings,
Interlocking, 2.25" lead
(Fits Atomlab Dose Calibrators 400,
500 and 500*Plus*.)
For additional protection from high energy activity

086-301 Dose Calibrator Shielding Rings,
Interlocking, 2.25" lead
(Fits Atomlab Dose Calibrators 100,
100*Plus*, 200 and 300.)
For additional protection from high energy activity

WELL INSERT



The Well Insert is included with the purchase of any Atomlab Dose Calibrator. The durable, clear Plexiglas insert is designed to protect the chamber from contamination and can be easily removed for cleaning. The insert will fit any well chamber with 2.5" x 10" interior dimension.

Keep a spare on hand for use while the other insert is being decontaminated.

086-241 Well Insert

ATOMLAB™ 950 THYROID UPTAKE SYSTEM

Accuracy, versatility and mobility

Complete, mobile, self-contained system:

- PC with high capacity hard drive, monitor, keyboard, trackball mouse and printer
- Windows® XP platform
- 1024-channel software-controlled multi-channel analyzer
- Software programs for thyroid uptake, wipe test, Schilling test, bioassay, hematology, manual MCA
- 2" x 2" NaI(Tl) detector with collimated shield (meets IAEA specifications)
- Distance measurement rod (detector-to-patient)
- Separate connections for probe and optional well counter for fast, convenient changeover
- Mobile stand with variable-height counterbalanced arm (for seated or supine patients)
- Industry exclusive two-year warranty

Comprehensive Testing and Utility Software:

Thyroid Uptake Program

- Supports multiple time-stamped uptake measurements
- Auto decay correction (or recount of standard)
- User-defined uptake protocols
- On-screen spectrum acquisition and analysis
- Reports include normal ranges, notes/comments, facility/physician/technologist

Wipe Test Program

- Complies with current regulatory standards
- Automatic subtraction of background activity
- Stores (and reports) individual wipe data on multiple wipe locations in any defined "area" - can be customized for package/shipment wipes

Schilling Test Program for use with a variety of commercial kits

- Program follows manufacturer's instructions for procedures and calculations

Bioassay Program for individual patients or employee groups

- Report choices include single bioassay report, summary of multiple bioassays on same patient, and summary report on multiple employees (gives most current assay result for each individual, for all isotopes of interest)

Hematology Program for red blood cell survival, blood/plasma volume, ERPF and GFR

Manual MCA Program for performing sample counts and analysis beyond the scope of the standardized programs

- Expandable library of 23 commonly used isotopes – pre-programmed with LLD, ULD, regulatory alarm levels, gain
- High-resolution graphics
- Keyboard entry of alphanumeric data, trackball mouse for selecting screens, menus, and options

System Administration / Quality Assurance Program

- Automatic daily calibration using Cs-137 source, with report
- Automatic high voltage adjustment, with verification
- Automatic Chi-Square test, with report
- Automatic communications test
- Administration summary QA report


All programs provide clear, concise reports for referring physicians, insurance providers, patient records and a database for physician and technologist identification.

Atomlab 950 software offers exceptional programming for Thyroid Uptakes, Wipe Testing, Bioassay, Schilling Tests, RBC Survival, Blood Volume, and Manual Count procedures.



The Atomlab 950 Thyroid Uptake System is an advanced multi-purpose spectrum analysis instrument designed for diverse nuclear medicine applications. Uptake studies, bioassays, Schilling tests, hematology tests, wipe tests, and other user-defined tasks are accomplished with speed and precision using this fully integrated computer-controlled instrument and its comprehensive selection of application software. The system's multi-channel analyzer has 1024 channels, with separate input connections for the probe and optional well counter. Engineered for mobility, durability, and operational efficiency, the Atomlab 950 handles clinical tests, safety compliance tasks, and system administration procedures quickly and accurately. Task-specific software provides step-by-step guidance throughout all procedures, automatically performs calculations, stores patient information and test results, and outputs clear, concise reports. The "Procedure Definition" screen is used to set up a facility's standardized test protocols by selecting the most appropriate program, editing program parameters to suit specific facility requirements, or by creating an entirely new procedure. When program design or modification is completed, the new program is stored for future use. A toolbar "help" button provides detailed system assistance.

The self-contained Atomlab 950 system is configured on a mobile platform with locking casters and a base that measures only 30" by 31.5" (76 x 80 cm). An upper shelf supports the computer, monitor, keyboard, and trackball mouse. A lower shelf supports the color inkjet printer. The base of the stand incorporates a shelf to hold the optional well counter. A 2" x 2" NaI(Tl) detector with collimator articulates on a counterbalanced arm. Twenty-two inches (56 cm) of vertical travel allows the probe to be positioned up to 55" (140 cm) from the floor to accommodate seated or supine patients. The probe swings more than 180° on the horizontal plane, and extends outward up to approximately 34" (86 cm) from the support column. This design makes positioning for uptake studies simple and comfortable for both patient and technologist.



Atomlab 950 Uptake Report displays all data sets with decay correct results and time intervals.

SPECIFICATIONS:

- Industry exclusive two-year warranty

MEDICAL SPECTROMETER HARDWARE

(call or visit the web at www.biodex.com for current specifications for computer)

Multi-Channel Analyzer

Channels: 1024

Inputs: Probe and well

Spectral Resolution: FWHM 10%

Count Rate: (Maximum) 100,000 cps

Count Rate Stability: 99%

Gross Count Rate Linearity: Within 5% up to 100,000 cps

Pulse Height Linearity: Within 2% (independent of detector)

Connectors: Signal (BNC); high voltage (MHV)

Power Supply: Regulated from 775-1225 VDC at 2 mAmps

Detector High Voltage Adjustment: Automatic H.V. adjustment for both probe & well; uses 10 μ Ci Cs-137 as the calibration source

MEDICAL SPECTROMETER SOFTWARE:

Programs: Thyroid Uptake, Wipe Test, Bioassay, Schilling Test, Hematology, Administration/ QA, Manual MCA

Radionuclides:

Factory Programmed: Au-198, Ba-133, Co-57, Co-58, Co-60, Cr-51, Cs-137, Fe-59, Ga-67, Hg-197, I-123, I-125, I-131, In-111, Ir-192, K-42, Na-24, Pd-103, Se-75, Sr-85, Tc-99m, Tl-201, Yb-169.

User Set: Unlimited user defined isotopes, setting ROI, half life, name, efficiency and gain

OTHER HARDWARE:

Probe: 2" x 2" NaI (TI) integral line scintillation detector with tube base

Uptake Stand:

Dimensions: 42" l x 31" w x 62" h (106.7 x 78.7 x 157.5 cm)

Collimated Shield: Flat field collimator meeting IAEA specifications

Arm: Counterbalanced, two section arm, moves 22.5" vertically and extends 29" horizontally from stand's vertical column

Casters: 3" maxi-lok

Weight: 296 lb (134.3 kg)

Shipping Weight: 430 lb (195 kg)

Optional:

187-246 Well Counter:

Detector: 2" x 2" NaI (TI) integral line scintillation detector with a .75" diameter x 1.44" deep well (1.9 x 3.7 cm)

Lead Shielding: 1" thick (2.5 cm)

Cover: .125" thick (.32 cm)

Connectors: Signal (BNC); High Voltage (MHV)

Weight: 54 lb (24.5 kg)

187-256 Well Counter:

Detector: 2" x 2" NaI (TI) integral line scintillation detector with a .75" diameter x 1.44" deep well (1.9 x 3.7 cm)

Lead Shielding: 2" thick (5 cm)


Cover: .125" thick (.32 cm)

Connectors: Signal (BNC); High Voltage (MHV)

Shipping Weight: 166 lb (75.2 kg)

Certification: ETL Listed to UL 60601-1 and CAN/CSA C22.2

No. 601.1-M90, EN 60601-1 and CE marked



Atomlab 950 simplifies daily calibration



187-140 Thyroid Uptake System,
Atomlab™ 950, Mobile, PC, 115 VAC
Mobile System includes:

- Monitor, standard keyboard, trackball, Windows® and Atomlab 950PC software, Printer and 1024 channel multi-channel analyzer
- 2" x 2" tube assembly and base
- Mobile support stand with collimator

187-130 Thyroid Uptake System,
Atomlab™ 950, Tabletop, PC, 115 VAC
Tabletop System includes:

- Monitor, standard keyboard, trackball, Windows® and Atomlab 950PC Software, Printer and 1024 channel multi-channel analyzer
- 2" x 2" tube assembly and base
- Tabletop stand with collimator

187-145 Thyroid Uptake System,
Atomlab 950, Mobile, PC, 230 VAC

187-135 Thyroid Uptake System,
Atomlab 950, Tabletop, PC, 230 VAC

Optional:

187-246 Well Counter, 1" lead (2.5 cm)

187-256 Well Counter, 2" lead (5 cm)

Related:

063-139 Rod Source, Cs-137, Calibrated, 0.1 μ Ci

063-100 Rod Source Set

Includes: Ba-133, Cs-137, Co-57, Na-22, Mn-54, Co-60 and Cd-109

101-103 Check Source, Cs-137, 10 μ Ci*

Uncalibrated, 1" dia x .25" thick (2.5 x .64 cm)

**Recommended Check Source for calibration of probe and well*

043-365 Thyroid Uptake Neck Phantom
Includes bottle carrier, capsule holder and 12 polyethylene bottles

086-342 Well Liners, Disposable, 100/pk



*For sample reports, testimonials
and warranty information,
visit us on the web www.biodex.com*

Int'l 631-924-9000 • www.biodex.com

ATOMLAB™ 950 THYROID UPTAKE SYSTEM COMPREHENSIVE TESTING AND UTILITY SOFTWARE

All programs provide clear, concise reports for referring physicians, insurance providers, patient records and a database for physician and technologist identification.

THYROID UPTAKE PROGRAM

- Supports multiple time-stamped uptake measurements
- Auto decay correction (or recount of standard)
- User-defined uptake protocols
- On-screen spectrum acquisition and analysis
- Reports include normal ranges, notes/comments, facility/physician/technologist

The Atomlab 950 provides a choice of four pre-programmed uptake study protocols, with provision to modify any procedure to suit specific requirements. Atomlab software guides the user through each step of any defined procedure to count the standard, lab background, patient thyroid, and patient background, and then automatically computes uptake percentage. The system allows multiple uptakes at varied intervals for each patient. Automatic decay correction calculates elapsed time between dose count and actual patient thyroid count. It is therefore not necessary to perform thyroid counts in exact hourly increments from times of dose counts. Results are calculated either from decay-corrected dose counts or by recounting a standard. Reports show user-defined normal ranges for comparison, and have provision for entering individualized notes and comments.

Preset Thyroid Uptake protocols:


1. **Automatic Count Time with Automatic Isotope Decay** – Atomlab Software calculates count time to satisfy a level of accuracy established by the technologist. High count rate offers the advantage of short count time when counting the capsule. Because the system calculates decayed activity of the dose, recounting a standard is unnecessary.
2. **Automatic Count Time with Recounting a Standard** – The system calculates the appropriate counting time based on the count rates and desired accuracy. A standard is maintained to recount for the decayed activity.
3. **Manual Count Time with Automatic Isotope Decay** – The technologist sets the count time for four separate counts used to calculate uptake – generally 60 seconds each. The system monitors elapsed time and calculates decayed value for administered dose.
4. **Manual Count Time with Recounting a Standard** – The technologist sets count time for four separate counts, then recounts standard and lab background, and counts the thyroid, and patient background – generally 60 seconds each.

WIPE TEST PROGRAM

- Complies with current regulatory standards
- Automatic subtraction of background activity
- Stores (and reports) individual wipe data on multiple wipe locations in any defined “area” - Can be customized for package/shipment wipes

Designed for use with an optional well counter, this program was designed to comply with current government regulations for counting and reporting results of wipes taken in contamination surveys. The program automatically calculates net contamination levels after subtracting measured background activity, and reports results in user-specified units of dpm, μCi , or kBq. Lower limit of detection (LLD) is automatically determined during the background count. Clicking on the “LLD Estimate” button gives the technologist a quick estimate of LLD to determine whether counting times are long enough. Specific wipe test areas may be named, added to the system database, and edited in much the same way as patient names and IDs. They may also be designated as “restricted”, “unrestricted”, or “sealed source”. The user enters the isotopes of interest for any specific survey area, and the system counts and reports results accordingly. A detailed spectrum analysis may be

performed (on demand) for any specific wipe. Each designated “area” may include up to ten separate wipe locations (or items) for testing and documenting in that area’s report (i.e. Area Name - “Camera Room”; Location - “counter top”, “table”, “cabinet door”, etc.). A “Wipe Preferences” screen allows the user to specify detector, the detector’s geometric efficiency, default count times, and the activity units to be used in reports.



Wipe test program screen.

MEASURED in (dpm)		
Isotope Name	Background	LLD
Ba-133	104.591	116.413
Co-57	23.862	47.046
Co-137	511.204	568.594
Wide Window	336.822	132.063

Lower limit of detection display.

BIOASSAY PROGRAM

For individual patients or employee groups


- Report choices include single bioassay report, summary of multiple bioassays on same patient, and summary report on multiple employees (gives most current assay result for each individual, for all isotopes of interest).

This program is used to assay for any isotope for which counting efficiency has been determined. Typical applications are monitoring activity in the thyroid of staff members who are exposed to I-123, I-125, or I-131, and for long-term monitoring and reporting of patients treated therapeutically with I-131. Three report formats are provided: The Individual Bioassay Report documents a single bioassay. An Individual Bioassay History Summary lists multiple assay results on the same individual. The Employee Summary Report lists the most current count for every individual staff member. The Bioassay Program allows the contamination trigger level to be set for the facility. The system determines and displays the lower limit of detection (LLD) for each isotope of interest after counting background, allowing the technologist to change parameters to adjust the LLD. Counts per minute are converted to dpm, μCi , or kBq using stored isotope efficiency values.

MANUAL MCA PROGRAM

For performing sample counts and analysis beyond the scope of the standardized programs.

- Expandable library of 23 commonly used isotopes – pre-programmed with LLD, ULD, regulatory alarm levels, gain
- High-resolution graphics
- Keyboard entry of alphanumeric data, trackball mouse for selecting screens, menus, and options



Manual program count screen.

In addition to providing “standardized” programs for routine test and survey procedures, the 950 has a “manual” program to accommodate other counting and spectrum analysis tasks. The user can select from three counting methods: Preset Time, Preset ROI Counts, and Continuous Counting (counts until STOP is pressed up to a maximum of 4095 seconds). Repetitive counts can be used to build a multiple count and/or multiple isotope reports. After any count, the Spectral Analysis option can be used to produce a Spectrum Analysis Report. The monitor shows when a test is in progress, and graphically displays the spectrum, time, counts, and cpm.

THE ADMINISTRATION/QUALITY ASSURANCE PROGRAM

- Automatic daily calibration with Cs-137 source, with report
- Automatic high voltage adjustment, with verification
- Automatic Chi-Square test, with report
- Automatic communications test
- Administration summary QA report
- Automatic high voltage adjustment
- Isotope efficiency calculations
- Multiple reports to document system performance
- Complete spectral analysis
- Isotope editing and entering



Chi-Square report.



Administration report.

There are a multitude of Administration functions, each designed to minimize the time normally spent maintaining accurate use of the system. There is an automatic set of the high voltage for each detector, and also an automatic calibration program using a Cs-137 source. The system will advise the user when a daily calibration has not been performed. The Chi-Square test program automatically counts the source 10 times and performs the Chi-Square calculations. The results are displayed and can be printed in a Chi-Square report.

The isotope editing program allows the user to enter custom isotopes or variations of the preset isotopes, varying the windows, gains, and half-lives, plus the user can set the counting efficiency for each isotope. There is a built-in calculator for calculating geometric efficiency. There is also a function that calculates the current activity of a calibrated source.

Site information is entered in the administration program containing the facility name, address, and a list of the technologists and physicians in the department. Reports can be printed showing the administration program's parameters.

VENTI-SCAN™ IV

Radioaerosol Delivery System




The Venti-Scan™ IV provides greater protection and uses a special disposable kit for greater trapping efficiency and unrestricted breathing.

- **System design reduces set-up time:**
 - kit automatically locks into position
 - one step oxygen connection to dedicated external port
 - precise injection port alignment
 - push button disposal of used kit
- Fully enclosed lead shielding from top to bottom
- Uses HEPA filter for increased trapping efficiency and resistance-free breathing
- Trapping efficiency greater than 99.9%
- Small baffle design ensures homogeneous distribution
- Mean particle size = 0.50 microns
- Lightweight and portable, weighs only 8 lb (3.7 kg)

Sharp, Clear Images - Easy Breathing

The Venti-Scan™ IV Radioaerosol Delivery System features a small baffle within the nebulizer to produce an optimal particle size, resulting in a sharp image, quickly. In addition, the kit includes a pleated hydrophilic HEPA filter which traps moisture. This makes it ideal for radioaerosol studies by impeding the exhaled radioaerosol particles from passing through. Tried and tested, the pleated contour increases surface area to decrease breathing resistance, making it virtually resistance-free with exceptional trapping efficiency.



Venti-Scan IV Radioaerosol Delivery System provides homogeneous distribution

Note the film study shown above. Images taken from a Venti-Scan IV study using our small baffled nebulizer kit are crisp and clear.

Designed For Efficiency

The Venti-Scan IV is designed to make performing a study more convenient for the technologist while providing superior images. When the kit is inserted into the Venti-Scan IV canister it automatically locks securely into position, assuring all port alignments. Oxygen connection is a simple attachment to a dedicated external port. The injection site on the kit is precisely angled to align with the canister port. They are positioned perfectly for a bull's eye every time. The new system also offers a quick, safe disposal method. Unplug the oxygen hose, invert the canister over a shielded waste container and push the release button to free the contaminated kit. This minimizes handling and exposure.



Venti-Scan IV shown mounted on IV Pole

Full Technologist Protection

The Venti-Scan IV shield is an enclosure providing lead-shielded protection from top to bottom. The Venti-Scan IV Disposable Kit includes everything needed for a single study including a comfortable, natural contour mouthpiece, HEPA filter, nose clip and disposal bag. The system uses clean-bore straight path tubing (superior to corrugated) to ensure that particles cannot get trapped in any internal ridges that typically cause clumping. The top of the canister has a shielded sliding port to accommodate the Venti-Pak Accessory Kit for ventilator-assisted patients. An IV pole mount is included with the shield for convenient positioning and administration.


Patients and technologists have always been comfortable and confident with the Venti-Scan. And, when the patient is comfortable, the procedure goes smoothly, without interruption. The end result is a superior study.




1 Insert Venti-Scan™ IV Kit in canister



2 Inject Tc-99m DTPA



3 Connect to O₂ supply; position patient comfortably for resistance-free breathing



177-075 Venti-Pak for ventilator assisted patients.

Call for a FREE** 30-Day Trial Evaluation.

When you purchase your first 75 kits - keep the canister... compliments of Biodex. (Value \$495.00)

Budgeted Savings

A recent cost comparison of departments averaging 20 studies per month resulted in an annual savings between \$2,000 and \$5,000. Venti-Scan disposables are already the most inexpensive kits available and can be even more time-effective by establishing a blanket order.

ONE PHONE CALL...ONE P.O.#...ONE YEAR...

Custom tailored to your department's volume, a blanket order assures that the disposables are there at the same time each month - when you need them.

Venti-Scan™ IV Radioaerosol Administration System:

177-090 Venti-Scan™ IV

Includes: Shielded canister with IV pole mount

Replacement:

177-091 Radioaerosol Kit, Disposable, for Venti-Scan™ IV ®*

Includes: 12" (30.5 cm) tubing, small particle delivery system with mouthpiece, HEPA filter, nose clip and disposal bag

177-092 Radioaerosol Kit, Disposable, for Venti-Scan™ IV ®*

Includes: 24" (61 cm) tubing, small particle delivery system with mouthpiece, HEPA filter, nose clip and disposal bag

Related:

177-075 Venti-Pak (adapter kit for ventilator assisted patients) ®*

CALL FOR FREE SAMPLE**

PEER PERSPECTIVE

"I love the Biodex Venti-Scan IV Radioaerosol System. It is very safe to use for both the technologist and the patient. I especially like that the system is fully enclosed by lead shielding. It is user friendly and works well with our patients. I highly recommend it."

- Saswati Roy, Chief Technologist, St. Francis Hospital, Milwaukee, WI

"We have been using the Biodex Venti-Scan IV Radioaerosol System in our facility since January 2005. We find the system is very compact, fits securely to an ordinary I.V. pole and is easy to load and dispose. The small size makes for easy decay and storage. Best of all, our patients do very well with this unit and the scans come out very nice. Thanks, Biodex, for a great product that makes routine patient exams much better!"

- Barbara J. Newman RT(R) (M) CNMT, Baylor Regional Medical Center, Plano, TX

"I've used the Venti-Scan IV Radioaerosol Delivery System for the last four years and I've been very pleased with the results. When used with the I.V. pole mount, this system makes performing lung scans quite easy."

The Venti-Scan IV delivers excellent aerosol images of the lungs. It produces images with outstanding count rates, which tends to make the aerosol image go very fast. This is good for us as technologists - and good for patients as well because it means they spend less time breathing and less time on the camera. The images for most patients have very negligible hot spots and that helps with physician interpretations. When imaging is complete, the kits are easy to dispose. There is no struggle to remove them from the shield, and being small, they take up little space in decay storage.

I have used several other lung ventilation kits over the years but the Venti-Scan IV seems to provide the best image quality. All in all, I highly recommend the Biodex Venti-Scan IV Radioaerosol System."

- Mark R. Blake RT(N), CNMT, UHHS Richmond Heights Hosp., Richmond Heights, OH

For physics tests us on the web at www.biodex.com/lungventilation

® DENOTES SINGLE USE ONLY

* NOTE: Sold in cases of 25.

** Offer applies to the United States only.

Int'l 631-924-9000 • www.biodex.com

AEROTECH™ I RADIOAEROSOL SYSTEM

The original Cadema Radioaerosol System



Shown with kit# 177-124

- *Enhanced Imaging*
- *Multiple views from one diagnostic study*
- *Completely portable and respirator compatible*
- *Completely Shielded*


Aerosol inhalation for radioimaging studies is an important tool for use in the localization and diagnosis of lung disease. AeroTech I puts information in the hands of the user to assist in diagnosis of a variety of lung diseases.

Enhanced Imaging

The AeroTech I delivers an even distribution of radioaerosol throughout the lungs, reaching the respiratory bronchioles and alveolar ducts and sacs. Excellent aerosol deposition allows the acquisition of multiple views from each diagnostic study. The end result is high quality imaging for confidence in patient diagnosis.

Respirator Compatible

For patients on a respirator, AeroTech I is easily adapted to fit on line using a single accessory tube.



177-325 Venti-Pak for ventilator assisted patients.

© DENOTES SINGLE USE ONLY

* NOTE: Sold in cases of 25.



Designed for Patient Comfort

Aerosol inhalation studies allow the patient to breathe normally throughout the procedure, minimizing the likelihood of patient noncompliance.

Superior Safety

Shield construction, disposable components and rapid aerosol delivery time combine to minimize radiation exposure to the patient and technologist. AeroTech I helps meet ALARA radiation protection objectives.

Flexibility

AeroTech I is available with a choice of designs. The 177-324 Delivery System generates a smaller MMAD particle than the standard 177-124 Delivery System. Use the 177-324 model for procedures where delivery time is not significantly restricted.

177-124 AeroTech™ I Disposable Radioaerosol Kit®*

LATEX FREE
Includes 24" (61cm) Tubing, Standard Nebulizer Delivery System with Mouthpiece, Bacteria Filter, Nose Clip and Disposal Bag

177-324 AeroTech™ I Disposable Radioaerosol Kit, with small particle delivery system®*

LATEX FREE
Includes 24" (61 cm) Tubing, Small Particle Nebulizer Delivery System with Mouthpiece, Bacteria Filter, Nose Clip and Disposal Bag

177-095 AeroTech™ I Shield ©

Related:

177-325 Venti-Pak, 5/pkg (adapter kit for ventilator assisted patients) ©

Welcome to the Biodex Cleanroom



Sometimes, it's not just what goes into a product that makes it great. It's also what stays out. That's why the Biodex Cleanroom, dedicated to assembly of our lung ventilation kits, is maintained at the highest standards.

The Biodex Cleanroom uses a system where conditioned and filtered air is forced downward by ceiling diffusers and exhausted at floor level. This allows achievement of three to five pounds of pressure, limiting the number of airborne particulate on products being assembled. Atmospheric conditions including temperature, humidity and airborne particles within a specific area, are monitored and charted daily to ensure no variation. Classified as an ISO 07 - Class 10,000 Cleanroom, as defined by Federal Standard 209E, this is an exceptionally clean environment.

"We take a lot of pride in keeping our Cleanroom operating at peak efficiency," explains the Consumer Assembly Supervisor at Biodex. "Every employee that enters this room adheres to strict guidelines and is required to wear clean room garments to protect the product and minimize the particulate in the environment. They are trained in appropriate operating procedures from putting on a hair net, frock, gloves and face mask, to proper assembly of each individual product. Nothing is left to chance."

As for having to "suit up" when entering and leaving the cleanroom on each shift, as well as being in an area where interaction with other Biodex personnel is somewhat limited, the staff really doesn't mind.

**It's a very clean job, but
somebody has to do it.**

BIODEX

www.biodex.com

1-800-224-6339

Int'l 631-924-9000

PULMONEX II® XENON SYSTEM

Double the breathing capacity of other systems

- Resistance-free breathing with 25-liter capacity
- Complete .125" lead shielding in upper cabinet for patient and operator safety
- Easy access to replace Drierite and Soda-Lime cartridges
- Effortless mobility for easy patient positioning
- Stylish steel cabinet
- Convenient built-in stainless steel tray for holding disposables, xenon gun, syringe, etc.

The Pulmonex II® Xenon System is the best choice for the performance of all regional ventilation studies. It's safe, simple to operate and affordable.

Resistance-Free Breathing

The injected bolus of xenon will reach the patient exactly when desired. Oxygen may be added to the system any time during the study with the press of a button. An in-line cartridge containing Soda-Lime absorbs CO₂, preventing acidosis. Large breathing passages, two 10-liter breathing bags (air-in and air-out) and motor-assisted airflow combine to provide resistance-free breathing.

Simple to Operate

All three steps of a Pulmonex II study (start up, equilibrium imaging and washout) are controlled by a single valve handle on the front panel. The valve directs the motor-assisted flow of gases throughout the system. A manually adjusted 15-minute timer initiates all functions, then automatically shuts down the system to complete the study after patient and system washout. With controls conveniently located on the front panel, the user can operate the system and observe the patient and gamma camera from one position.

Panel controls are clearly marked for each mode of the procedure with large viewing windows to make it easy to monitor the patient's breathing. A trim, clean design, large handles and total mobility permit easy positioning of the system for studies in both seated and supine positions.



Easy access to Soda-Lime and Drierite cartridges, supplied with the system. Xenon Convenience Kits, including pre-filled Soda-Lime and Drierite cartridges for single use only, are also available



Pulmonex II is a completely closed system that provides maximum, reliable test results with minimal effort

Designed for Safety

Internal systems of the Pulmonex II are shielded for patient and operator safety. The system features two built-in gas traps that operate with a blower fan. Exhaled xenon is pulled through activated charcoal housed within two .125" lead shielded "U"-shaped traps. The double traps extend the life of your charcoal and provide a lengthy migration path for xenon effluent, allowing greater decay and absorption before exhaustion. A cartridge containing Drierite serves as a moisture absorber for air passing into the trap. The charcoal trap can then more effectively remove xenon effluent after each study. Airflow regulation of the trap pump assures complete patient and system washout. Averaging 30-50 studies per month, the charcoal trap will last approximately one year; charcoal traps are easily replaced.

A disposable bacteriostatic filter, used in conjunction with a disposable mask or mouthpiece, prevents system contamination.

SPECIFICATIONS:

Dimensions: 48.5" h x 22" depth x 20.5" w (123.2 x 55.9 x 52.1 cm)
 Motor: UL approved, 12 volt DC
 Electrical Requirements: 115 vac, 1 amp or 230 VAC, 0.5 amp, 50/60 Hz
 Shipping Weight: 375 lb (172.5 kg)
 Certifications: ETL and cETL Listed
 Warranty: One year



*For warranty information,
 visit us on the web www.biodex.com*

132-503 Xenon System, Pulmonex II®,
 Double-Trap, 115 VAC

Related:

136-755 Xenon Trap Monitor
130-900 Pulmonex Kit, Free Breathing Hose

Replacement:

132-319 Charcoal Trap*

For any Pulmonex II

132-555 Cartridge, Refillable,
 Soda-Lime & Drierite

**50% discount on charcoal trap with blanket order commitment
 Call for details.*

** Offer applies to the United States only.*

NEW**AUTOMATIC XENON DISPENSER**

*Assures quick, accurate and easy delivery of
 xenon directly to the patient*



The Automatic Xenon Dispenser attaches to the front of the Pulmonex System without tools. Simply load the xenon vial into the supplied plunger assembly, and at the precise moment you want the delivery of xenon, simply press the dispense button. Xenon is automatically delivered to the patient.

150-315 Dispenser, Xenon, Automatic
Includes 48" Flexible Tubing with Luer-Lock Adapter

NEW**XENON TRAP MONITOR**

Mounts directly to the Pulmonex II Xenon System for continuous monitoring

The Xenon Trap Monitor continuously monitors trap effluent during a xenon study, meeting compliance requirements in most states. An exhaust hose connects from the Pulmonex II exhaust port to the Xenon Trap Monitor's intake port, making it the perfect complement to the Pulmonex II® Xenon System to ensure a safe environment.

The Xenon Trap Monitor is simple to operate; all controls are located on the front panel. Counting results are displayed with both audible and visual signals to indicate when the xenon trap exhaust port exceeds the threshold of 99 pCi/mL. The unit can also check background levels and perform a self-test for proper operation using a check source.

The Xenon Trap Monitor includes all the hardware necessary to mount the unit directly to the Pulmonex II Xenon System.

SPECIFICATIONS:

Dimensions: 6" w x 7" l x 4" h (15.2 x 17.8 x 10.2 cm)
 Lead Shielding: .5" thick (1.3 cm)
 Input: 22 mm Hose Adapter
 Detector: Halogen Quenched GM Tube, 1.1 cm dia, 2 mg/cm² mica window
 Voltage: 500 volts regulated
 Buttons: On/Off, Next/Mute
 Display: 4 Digit LED, 4 RGB Function LEDs (Self Test, Background, Check Source or Count)
 Readings: pCi/ml or counts
 Speaker: Internal, beeps in alarm mode
 Background Count Time: 1 minute
 Power: 18 volts, UL approved, external; 115 VAC power adapter
 Weight: 5.25 lb (11.55 kg)
 Certification: ETL and cETL Listed
 Warranty: One year



136-755 Monitor, Xenon Trap
Mounts to Pulmonex II Xenon Systems

Related:

101-103 Check Source, Cs-137, 10 µCi
Uncalibrated, 1" dia x .25" thick (2.5 x 64 cm)

XENON CONVENIENCE KITS™*

Face Mask or Mouthpiece Kits - your choice, your convenience.

- *Air-Cushioned Face Mask Kits with clear ultra-flex expandable tubing from 6" to 24", or*
- *Mouthpiece Kits with Direct Dose Administration adapter for leak-proof Xenon delivery*

Biodex Xenon Convenience Kits™ bring together all the components needed to complete a single xenon study. Xenon can be administered via direct dose administration, that conveniently luer locks without the use of a needle; or by injection port, which requires a syringe needle. Easy to use, simple to order and disposable, Biodex Xenon Convenience Kits™ are more than a bargain, they are a sensible and time-saving

investment. In fact, our complete kits can prolong the life of your xenon charcoal traps by ensuring that the pre-filled Drierite cartridges are fresh for every scan while the pre-filled Soda-Lime (CO₂ absorber) cartridges eliminate possible breakdown of granules that can lead to a clogged system and hinder patient air flow.



XENON CONVENIENCE KITS™

AIR-CUSHIONED™ Face Mask:

- 132-680** Convenience Kit, Face Mask ®
Includes: Bacteria filter
- 132-681** Convenience Kit, Face Mask ®
Includes: Bacteria filter and ultra-flex tubing (shown)
- 132-781** Convenience Kit, Face Mask ®
Includes: Bacteria filter, ultra-flex tubing, Drierite & Soda-Lime cartridges



XENON CONVENIENCE KITS™

AIR-CUSHIONED™ Face Mask w/Injection Port:

- 132-690** Convenience Kit, Face Mask with Injection Port ®
Includes: Bacteria filter (shown)
- 132-691** Convenience Kit, Face Mask with Injection Port ®
Includes: Bacteria filter and ultra-flex tubing
- 132-784** Convenience Kit, Face Mask with Injection Port ®
Includes: Bacteria filter, ultra-flex tubing, Drierite & Soda-Lime cartridges

NEW XENON CONVENIENCE KIT™

AIR-CUSHIONED™ Face Mask w/Luer-Lock Injection Port

- 132-684** Convenience Kit, Face Mask with Luer-Lock Injection Port ®
Includes: Bacteria filter, ultra-flex tubing, elbow and 90° luer-lock injection port

CALL FOR FREE SAMPLE**

® DENOTES SINGLE USE ONLY

* NOTE: Sold in cases of 25.

** Offer applies to the United States only.

XENON CONVENIENCE KITS™*


Direct Dose Administration for leak-proof xenon delivery



XENON CONVENIENCE KITS™

AIR-CUSHIONED™ Face Mask
w/Direct Dose Administration:

- 132-699** Convenience Kit, Face Mask with Direct Dose®
Includes: Bacteria filter
- 132-692** Convenience Kit, Face Mask with Direct Dose®
Includes: Bacteria filter and ultra-flex tubing
- 132-793** Convenience Kit, Face Mask with Direct Dose®
Includes: Bacteria filter, ultra-flex tubing, Drierite & Soda-Lime cartridges (shown)



XENON CONVENIENCE KITS™

with Mouthpiece and Direct Dose Administration Adapter:

- 132-770** Convenience Kit, Mouthpiece®
Includes: Bacteria filter and nose clip
- 132-774** Convenience Kit, Mouthpiece and Administration Adapter®
Includes: Bacteria filter, nose clip and ultra-flex tubing (shown)
- 132-771** Convenience Kit, Mouthpiece and Administration Adapter®
Includes: Bacteria filter, nose clip, ultra-flex tubing, Drierite & Soda-Lime cartridges

XENON CONVENIENCE KIT™

AIR-CUSHIONED™ Face Mask and Direct Dose Administration Adapter:

- 134-772** Convenience Kit, Face Mask and Administration Adapter®
Includes: Bacteria filter and ultra-flex tubing

CALL FOR FREE SAMPLE**


® DENOTES SINGLE USE ONLY

* NOTE: Sold in cases of 25.

** Offer applies to the United States only.

Int'l 631-924-9000 • www.biodex.com

AIR-CUSHIONED FACE MASK*



The air-cushioned face mask has a pre-filled air-cushion (medium inflation pressure) that molds to the contour of the patient's face providing a tight, leak-proof and comfortable seal. Cushion pressure is easily adjusted by inserting a standard syringe without needle in the two-way valve located on the mask bottom. The flexibility of cushion pressure allows optimum surface contact for every patient. Transparency of the entire mask allows continuous visual identification of patient's vital signs.

AIR-CUSHIONED™ Face Mask
without Injection Port:

132-685 Face Mask without Injection Port, Adult ☉



AIR-CUSHIONED™ Face Mask
with Injection Port:

132-695 Face Mask with Injection Port, Adult ☉

AIR-CUSHIONED FACE MASK*

with Direct Dose Administration




In addition to the patient comfort provided by the soft, pliable air cushion of the disposable face mask, the direct dose feature reduces patient anxiety because no needle is used. The syringe luer locks tightly to the direct dose tube, away from the patient's face. The one-way valve prevents the gas from reversing through the mask luer connector, allowing the technologist to disconnect immediately after injecting.

AIR-CUSHIONED™ Face Mask
with Direct Dose Administration:

132-698 Face Mask with Direct Dose
Administration, Adult ☉

FACE MASK HARNESS



This traditional Face Mask Harness is made of soft rubber and can be adjusted for any head size. The square rubber base is comfortably positioned on the back of the patient's head and the harness tails are brought around to the front. The small holes fit snugly on the hook ring and hold the mask firmly in place.

139-677 Face Mask Harness, Adult

☉ DENOTES SINGLE USE ONLY

* NOTE: Sold in cases of 25.

TRU-FIT DISPOSABLE MOUTHPIECE



Contoured to the natural shape of the mouth, this mouthpiece can be held gently, but securely — eliminating the need to clamp down or the possibility of “popping” out of the patient’s mouth and interrupting the study. The soft plastic material is transparent so that any obstruction can be immediately detected by the technologist. Pleasantly scented with vanilla, this mouthpiece makes the procedure a little more comfortable for the patient.

130-551 Mouthpiece, Disposable,
Vanilla Scented ®


DISPOSABLE NOSE CLIP



Economical, this nose clamp is all plastic for single patient use. Used with any mouthpiece.

130-100 Nose Clip, Disposable, 100/pkg ®

DISPOSABLE BACTERIA FILTER*



Bacteria filters are used to reduce the possibility of cross contamination. The single-use filter is placed in line between a delivery tube and disposable mouthpiece or face mask. Electrostatically charged filter media is 99.9% effective in bacteria/virus retention while maintaining low breathing resistance. Two filters can be piggy-backed together for potential high-risk studies.

132-750 Bacteria Filter ®

BLANKET ORDERS

*You may qualify to receive a 50% discount** on a replacement Charcoal Trap for your Pulmonex just by placing a Blanket Order.*

- There is no binding contract. You can cancel at any time without penalty or change the quantities throughout the life of the blanket.
- You set the quantities and delivery intervals — based on your department needs.
- A single purchase order establishes a tailored blanket for up to one year.
- A Blanket Order exempts you from price increases on blanketed items.
- Disposables are there when you need them. No more rush orders, expensive overnight shipping charges or equipment downtime.
- Maximize your storage space — by using ours.

***Offer applies to United States only.*


® DENOTES SINGLE USE ONLY

* NOTE: Sold in cases of 25.

** Offer applies to the United States only.

Int'l 631-924-9000 • www.biodex.com

CORRUGATED TUBING



Corrugated Tubing is scored and capped at 6" intervals for easy cutting and firm attachment to fittings.

139-680 Tubing, Corrugated, 100 ft/roll ☉

TUBING SPLICE



130-639 Tubing Splice, 22 mm Male, 10/pkg ☉

Tubing Splice is used to connect sections of 22 mm corrugated tubing to individual mask and filter.

"Y" CONNECTOR



185-302 "Y" Connector, luer plug ☉

139-102 "Y" Connector, solid plug ☉

☉ DENOTES SINGLE USE ONLY

DRIERITE

Moisture absorber - new blue and white blend



139-101 Drierite, 1 lb

139-104 Drierite, 4.5 lb

132-772 Drierite, Disposable Cartridge ☉

Drierite serves as a moisture trap for the air going into the charcoal trap of the Pulmonex Xenon System. Drierite is blue and white when dry and turns pink when it has absorbed maximum moisture. Available in two convenient sizes, Drierite is packaged in an air-tight container with a wide-access mouth.

SODA-LIME

CO₂ absorber



130-019 Soda-Lime, 1 lb

130-020 Soda-Lime, 4.5 lb

132-773 Soda-Lime, Disposable Cartridge ☉

Soda-Lime, a highly efficient CO₂ absorber, prevents patients from re-breathing carbon dioxide and subsequent acidosis. Available in two convenient sizes, Soda-Lime is packaged in an air-tight container with a wide-access mouth.

FREE-BREATHING PULMONEX® HOSE KIT

When unrestricted breathing is critical



The Pulmonex Xenon System has negligible internal breathing resistance, yet excess breathing resistance can be due to improper hose apparatus to the patient, particularly if long, small diameter corrugated hoses are used.

This inexpensive, high quality, “free breathing” hose kit includes all required components and adapters to replace disposables on older Pulmonex styles (#130-500).

The kit contains two 130-901 “Clean-Bore” hoses (36" l x 1 1/8" dia), two 130-904 hose clamps, and a 130-902 non-rebreathing anesthesia valve (inlet and outlet 7/8" O.D.), which permits the use of standard disposable bacteria filters and masks while providing unrestricted airway.

130-900 Pulmonex Kit, Free-Breathing Hose

Replacement:

130-901 Hose, “Clean-Bore”, 36" l

130-902 Valve, Anesthesia

Hans Rudolph, non-rebreathing

130-903 Adapter, Hose

.875" x 1.125" dia (2.2 x 2.9 cm)

130-904 Clamp, Hose

XENOVAP™ ADMINISTRATION SYSTEM



For ventilator-assisted patients

The XenoVAP™ Administration System is an inexpensive, disposable device used to administer Xenon-133 and to collect the expired gas on ventilator-assisted patients. The system accommodates any standard tracheotomy tube (15 mm I.D.).

The ventilator-assisted patient is removed from their regular device and placed on this system for short term ventilator assistance during the procedure. The patient's breathing is first stabilized and then Xe-133 is injected into the system. The VAP is squeezed and held until sufficient counts are collected by the gamma camera system or until the patient must breathe. The pump is squeezed and released, forcing the patient to breathe through the system to collect an “equilibrium phase” image. When sufficient data has been collected, the valve is opened and the patient breathes room air and Xe-133 is eliminated into a collection bag. The gamma camera system should immediately begin collecting data for the “washout phase” image.

When the study is complete, the entire rebreathing system is disposed of in accordance with NRC regulations.


System contents include:

- Ventilator-Assist Pump (VAP) with 2037 cc volume capacity
- Endotracheal Tube Connection (15 mm I.D.) accommodates standard tracheotomy tube
- Injection Site for administration of Xe-133
- Oxygen Inlet to collection bag
- CO₂ Absorber (pre-filled with Soda-Lime and dust filter)
- 35-liter Collection Bag

060-139 Administration System, XenoVAP™ ®

XENON-133 REBREATHING SYSTEM

Disposable pre-packaged complete kits



Shown with **AIR-CUSHIONED™** Face Mask with Direct Dose Administration



System 060-133 shown with mouthpiece

SECTION 8


The Xenon-133 Rebreathing System provides the clinician with a simple, safe and inexpensive method of administering Xenon-133 to perform perfusion steady state and washout studies. A 35 liter bag provides ample volume for patient maintenance and collection of expired Xenon-133.

This disposable system is constructed of inert, non-permeable plastic, thereby precluding xenon absorption or transfer through the device's components. The system is designed to accept any xenon cylinder or syringe/gun administration system.

Choose the mouthpiece or the face mask with direct dose; each comes as a complete disposable system including a pre-filled Soda-Lime absorber cartridge to prevent messy measuring and to save time. Whichever system you prefer, the xenon gas is administered directly to patient for single breath static imaging. Thereafter, the patient rebreathes to equilibrium.

Perfusion studies are easily accomplished by placing the patient on the system and injecting the dissolved xenon/saline intravenously while the patient holds their breath.

Washout is accomplished by opening one valve, thereby admitting external air to the intake side. Upon completion of patient washout, the system is sealed by means of an appropriate clamp, such as a hemostat or kelley clamp. The user removes the system to a storage area for decay or to a hood for venting. The unit is then repackaged in the resealable bag provided and disposed of according to NRC regulations.



060-133 Xenon-133 Rebreathing System, Mouthpiece ②
Kit includes: Mouthpiece (130-551)

060-137 **AIR-CUSHIONED™**
Xenon-133 Rebreathing System, Adult ②
Kit includes: Air-Cushioned™ Face Mask with Direct Dose Administration (132-698)

② DENOTES SINGLE USE ONLY

ROD SOURCES



SPECIFICATIONS:

Dimensions: 2.96" l x 0.47" dia (76 x 11.9 mm)
Nominal Total Activity: 0.1 μ Ci

063-138 Rod Source, Ba-133

063-139 Rod Source, Cs-137

063-137 Rod Source, Co-57


063-140 Rod Source, Ge-68 simulated F-18

063-100 Rod Source Set

*Set includes: Ba-133, Cs-137, Co-57,
Na-22, Mn-54, Co-60 and Cd-109*

Not available for export.

SPOT MARKER SOURCE*



A Spot Marker Source is easily taped to the patient for purposes of orientation while performing a camera study.

SPECIFICATIONS:

Co-57 Spot Marker Sources
Dimensions: 1" dia x .25" thick
(25 x 6 mm) clear lucite disk
Spot: .125" dia (3 mm)
Activity:
099-289: 50 μ Ci
099-291: 100 μ Ci

099-289 Marker Source, Spot, Co-57, 50 μ Ci

099-291 Marker Source, Spot, Co-57, 100 μ Ci

**A photocopy of your NRC or Agreement State License must accompany orders for radioactive sources and must clearly indicate your authority to possess the source being ordered.*

Not available for export.

FILLABLE POINT MARKER SOURCE



Save time for camera-to-patient alignment

SPECIFICATIONS:

Dimensions: 1" dia x .5" thick
Volume: 2 cc

043-274 Marker Source, Fillable Point, 4/set

Marker Sources are easily filled with the same radionuclide to be used in imaging procedures. Small clear plastic receptacle, 1" diameter x .5" thick, has a centered channel to contain channel to contain 0.2 cc of the nuclide. A nylon screw-plug tightens against an O-ring completing a tight, safe seal. Marker is re-usable or nuclide can be decayed or removed.

CHECK SOURCE



Instrument functionality is easily assessed with a Cs-137 Check Source. The activity is 10 μ Ci. No license is necessary.

101-103 Check Source, Cs-137, 10 μ Ci
(uncalibrated), 1" dia x .25" thick (2.5 x .64 cm)
Not available for export.

LIVER MARKER/RULER



The Liver Marker/Ruler is designed to mark, outline, and measure the liver while performing routine studies. When placed over the area of interest, a


measurement is taken by counting the space between holes. The marker is made of vinyl-coated leaded rubber for flexibility and easy cleaning. Each unit is 5.5" long and 2" wide with holes 1 cm apart. Two units are supplied as a kit.

SPECIFICATIONS:

Dimensions: 5.5" l x 2" w (14 x 5 cm)
Holes: 9 mm dia
Lead Equivalency: 1 mm

123-500 Liver Marker/Ruler

PEN POINT MARKER SOURCES*



The Pen Point Marker Source contains Co-57 in a ceramic matrix at the end of a 9.26" (23.5 cm) anodized aluminum rod. The pen-shaped rod screws into a brass cap which shields the active point. The Pen Point Marker Source is used in tracing the outlines of anatomical features on a patient. The trace appears almost instantly on the camera display. Contained activity is supplied as a nominal value $\pm 15\%$.

SPECIFICATIONS:

Dimensions: 9.26" l x .25" dia (235 x 6.35 mm)
Nuclide: Co-57
Source: 2 mm bead
Nominal Activity:
063-700 100 μCi
063-701 200 μCi
Source Cover: Brass cap

063-700 Marker Source, Pen Point, 100 μCi

063-701 Marker Source, Pen Point, 200 μCi

*A photocopy of your NRC or Agreement State License must accompany orders for radioactive sources and must clearly indicate your authority to possess the source being ordered.

Not available for export.

RADIOACTIVE RULERS*



These slender plastic rulers contain alternate 1 cm sections of Co-57 (10 μCi) and cold material. They are placed on the area to be measured, and the camera base line and window are adjusted for Co-57. To imprint the film, ten counts are sufficient. The ruler's image appears on the scintiphoto, allowing definition of anatomical location and/or organ size by counting the number of visible "hot" and "cold" 1 cm areas.

SPECIFICATIONS:

124-230 Rigid Radioactive Ruler

Nuclide: Co-57

Nominal Activity: 160 μCi

"Cold" Areas: Blank

Length: 17 cm (effective length 15 cm)

Shipping Weight: 1 lb (0.5 kg)

124-231 Flexible Radioactive Ruler

Nuclide: Co-57

Nominal Activity: 460 μCi

"Cold" Areas: Lead

"Hot" Areas: 23 @ 20 μCi each

Length: 47 cm (effective length 45 cm)

Shipping Weight: 3 lb (1.5 kg)

Radioactive Rulers:

124-230 Rigid Ruler, 160 μCi

124-231 Flexible Ruler, 460 μCi

*A photocopy of your NRC or Agreement State License must accompany orders for radioactive sources and must clearly indicate your authority to possess the source being ordered.

Not available for export.

RECTANGULAR COBALT-57 FLOOD SOURCES*



All Flood Sources are shipped with the patented radshield.

Cobalt-57 flood sources are used to test the response uniformity of gamma cameras. This is a critical quality control procedure to ensure that the camera response is uniform over the total head area(s). Cobalt-57 is uniformly dispersed in an epoxy matrix fully sealed in a high integrity ABS encapsulation. Radionuclide purity is greater than 99.9% (combined Co-56/Co-58 is less than 0.08% at source reference date). All flood sources are supplied with a uniformity test statement, wipe test certificate, handling and storage information and a custom decay calendar. All flood sources are shipped with the patented RadShield. RadShield provides comparable shielding to traditional hard-cases but at half of the weight.

SPECIFICATIONS:

Cobalt-57 Flood Sources are manufactured and 100% tested to have a coefficient variation of less than $\pm 1.0\%$ and integral non-uniformity of less than $\pm 2.5\%$. SPECT imaging procedures require less than $\pm 1\%$ coefficient of variation.

Rectangular Cobalt-57 Flood Sources

Dimensions: 25.4" l x 17.9" w x .7" thick (65 x 45 x 1.8 cm)

Active Dimensions: 23.9" l x 16.5" w x .22" thick (61 x 42 x .56 cm)

Cardiac Rectangular Flood Source

Dimensions: 17.1" l x 10.9" w x .7" thick (43.4 x 27.7 x 1.8 cm)

Active Dimensions: 15.5" l x 9.25" w (39.4 x 23.5 cm)

Rectangular Flood Sources, Cobalt-57:

Includes: Shielded storage case

043-830 Flood Source, Rectangular, 5 mCi

043-840 Flood Source, Rectangular, 10 mCi

043-855 Flood Source, Rectangular, 15 mCi

043-845 Flood Source, Rectangular, 20 mCi

043-860 Flood Source, Cardiac, Rectangular, 10 mCi

Not available for export.

FREE DISPOSAL WHEN YOU PURCHASE A NEW FLOOD SOURCE FROM BIODEX. (Applies to United States Only)

*A photocopy of your NRC or Agreement State License must accompany orders for radioactive sources and must clearly indicate your authority to possess the source being ordered.

DOSE CALIBRATOR VIAL REFERENCE SOURCES*



The daily calibration of dose calibrators is recommended to ensure accurate and reproducible instrument response. Calibration is easily achieved and maintained by the use of long-lived reference sources.

These sources are solid cast epoxy, 20 ml active volume in the 27 ml Vial E. They are calibrated

with $\pm 5\%$ accuracy at the 99% confidence level, NIST traceable. Dose Calibrator Reference Sources are registered with the U.S. Food and Drug Administration Center for Devices and Radiological Health and the U.S. Nuclear Regulatory Commission.

Each source includes a certificate of calibration, a leak test certificate, and a radiation safety and handling sheet. The source is packaged in an individual lead shield that is color coded and vinyl covered to eliminate exposure to the lead. Dose Calibrator Reference Sources are available individually or as an economical set.

Vial Dose Calibrator Reference Sources:

All sources calibrated to $\pm 5\%$

063-562 Source, Ba-133, 250 μCi

063-350 Source, Co-60, 50 μCi

101-356 Source, Cs-137, 200 μCi

063-261 Source, Co-57 simulated Tc-99m, 5 mCi

063-720 Source, Co-57 simulated Tc-99m, 10 mCi

063-586 Reference Source Set

Set includes: Co-57, 5 mCi;

Cs-137, 200 μCi ; Ba-133, 250 μCi

Not available for export.

Reference:

1. Guide for Preparation of Application for Medical Programs, U.S. Nuclear Regulatory Commission, Regulatory Guide 10.8 and American National Standard, Calibration and Usage of "Dose Calibrator" Ionization Chambers for the Assay of Radionuclides, ANSI N-42. 13-1986.

DOSE CALIBRATOR SYRINGE REFERENCE SOURCES*



The daily calibration of your dose calibrator is recommended to ensure accurate and reproducible instrument response. Calibration, using long-lived standards, should be performed in a manner that most closely represents how you use your dose calibrator. The Dose Calibrator Syringe Reference Source was designed for imaging facilities that obtain their radiopharmaceuticals in unit dose syringes.

The sources are solid cast epoxy, 3 ml active volume in a 5 cc "mock" syringe. They are calibrated within $\pm 5\%$ accuracy at 99% confidence level, NIST traceable.

Each source includes a certificate of calibration, leak test certificate, and radiation safety and handling sheet. The source is packaged in an individual lead shield that is color coded to the source. Dose Calibrator Syringe Reference Sources are available individually or as an economical set.

Syringe Dose Calibrator Reference Sources:

All sources calibrated to $\pm 5\%$

063-361 Source, Ba-133, 250 μCi

063-360 Source, Cs-137, 200 μCi

063-362 Source, Co-57, 5 mCi

063-364 Source, Ge-68, Simulated F-18, 500 μCi

063-363 Syringe Source Set

Set includes: Ba-133, 250 μCi ;

Co-57, 5 mCi; Cs-137, 200 μCi

Not available for export.

FREE DISPOSAL WHEN YOU PURCHASE A NEW REFERENCE SOURCE FROM BIODEx. (Applies to United States Only)

*A photocopy of your NRC or Agreement State License must accompany orders for radioactive sources and must clearly indicate your authority to possess the source being ordered.


JASZCZAK SPECT PHANTOM

ECT Phantom for PET and SPECT



Deluxe SPECT Phantom is shown

Inserts for PET and Spect Phantoms



Cold Rods

Cold Spheres

The Jaszczak SPECT Phantom provides consistent performance information for any SPECT or PET system. Multiple performance characteristics of camera-based SPECT systems are evaluated from a single scan of the phantom.

On-axis and off-axis transverse line spread function may be easily measured without removing the cover plate. Measurements of full-width-half (or tenth) maximum can be readily determined, either in air or in water.

The Phantom is used for:

- System performance valuation of: Collimator, Artifacts, Calibration, and Reconstruction Parameters
- Acceptance testing
- Routine quality, assurance and control
- Evaluation of center-of-rotation error
- Evaluation of non-uniformity artifact
- Evaluation of changes of radius-of-rotation on spatial resolution
- Evaluation of reconstruction filters on spatial resolution
- Evaluation of attenuation and scatter compensation
- Single slice volume sensitivity
- Total system volume sensitivity
- Lesion detectability

SPECT Phantoms are available in two models. The Deluxe Phantom is used for high resolution cameras. The Standard Phantom is used for lower resolution cameras.

SPECIFICATIONS:

Cylinder Interior Dimensions: 8.5" dia x 7.32" h (21.6 x 18.6 cm)

Cylinder Wall Thickness: 0.125" (3.2 mm)

Volume: 6.9 L

Volume With Inserts: 6.1 L

Cold Rod Insert Height: 3.46" h (8.8 cm)

Height of Spheres From Base Plate: 5" h (12.7 cm)

043-750 SPECT Phantom, Deluxe

Cold Rod Dimensions: 4.8 mm, 6.4 mm, 7.9 mm, 9.5 mm,

11.1 mm, 12.7 mm

Solid Sphere Diameters: 9.5mm, 12.7 mm, 15.9 mm, 19.1 mm, 25.4 mm, 31.8 mm

Shipping Weight: 15 lb (6.9 kg)

043-762 SPECT Phantom, Standard

Cold Rod Dimensions: 6.4 mm, 7.9 mm, 9.5 mm, 11.1 mm,

12.7 mm, 16.0 mm

Solid Sphere Diameters: 12.7 mm, 15.9 mm, 19.1 mm,

25.4 mm, 31.8 mm, 38 mm

043-750 Phantom, SPECT, Deluxe

043-762 Phantom, SPECT, Standard

Related:

043-763 Phantom Insert, Hollow Spheres

043-730 Phantom Insert, Triple Line

043-777 Phantom Insert, Cardiac

Phantom inserts are featured on page 107.

HOLLOW SPHERE INSERTS



- *Designed for use in all circular and elliptical SPECT cylinders*
- *Simulates hot and cold spherical "lesions"*
- *Quantitative evaluation of spatial resolution/object size, attenuation and scatter effects*
- *Evaluation of quantitative ECT reconstruction methods*

SPECIFICATIONS:

Set: Six hollow spheres (each individually removable and fillable)

Diameter: I.D.: 9.9 mm, 12.4 mm, 15.6 mm, 19.7 mm, 24.8 mm, and 31.2 mm

Volume of Spheres: 0.5 ml, 1.0 ml, 2.0 ml, 4.0 ml, 8.0 ml, and 16.0 ml

043-763 Phantom Insert, Hollow Spheres

Related:

043-765 Phantom, SPECT, Flangeless

043-750 Phantom, SPECT, Deluxe

043-762 Phantom, SPECT, Standard

043-740 Phantom, Lung-Spine


JASZCZAK FLANGELESS DELUXE SPECT PHANTOM

ECT phantom for SPECT, meets ACR requirements



Flangeless SPECT Phantom

Inserts for PET and Spect Phantoms



Cold Rods

Cold Spheres

The flangeless SPECT Phantom and the flangeless PET Phantom meet the requirements set by ACR. The appropriate phantoms provide consistent performance information for any

SPECT or PET system. The SPECT Phantom is the same as the Deluxe SPECT Phantom except that the flanges have been removed. The small SPECT Phantom is an ACR recommended phantom for small field of view dedicated cardiac SPECT systems. Multiple performance characteristics of camera-based SPECT systems are evaluated from a single scan of the phantom.

On-axis and off-axis transverse line spread function may be easily measured without removing the cover plate. Measurements of full-width-half (or tenth) maximum can be readily determined, either in air or in water.

The Flangeless SPECT Phantom is used for high resolution cameras, complying to the ACR requirements.

SPECIFICATIONS:

SPECT and PET Phantoms:

Cylinder Interior Dimensions: 8" dia x 7.32" h (20.4 x 18.6 cm)

Volume: 6.4 L

Cold Rod Insert Height: 3.46" h (8.8 cm)


Cold Rod Diameters: 4.8, 6.4, 7.9, 9.5, 11.1 and 12.7 mm

Height of Spheres From Base Plate: 5" h (12.7 cm)


Solid Sphere Diameters: 9.5, 12.7, 15.9, 19.1, 25.4 and 31.8 mm

ESSER FLANGELESS DELUXE PET PHANTOM

ECT phantom for PET, meets ACR requirement



Flangeless PET Phantom™



Flangeless PET Phantom Lid™

Image of PET Phantom Lid

PET Phantom Lid Containers

SPECT Phantom, Small

Cylinder Interior Dimensions: 5.5" dia x 5.9" h (14 x 15 cm)

Cold Rod Insert Height: 1.57" h (4 cm)

Cold Rod Diameters: 4.8, 6.4, 7.9, 9.5, 11.1, and 12.7 mm

Height of Spheres From Base Plate: 3.1" h (7.85 cm)

Solid Sphere Diameters: 6.4, 9.5, 12.7, 15.9, 19.1 and 25.4 mm

Flangeless Esser PET Phantom Lid™

Refillable thin-walled cylinders: 8,12,16,25 (x3)mm

Solid cylinder (Teflon®): 25mm

Cylinder height: 1.5 in

Lid Closure: Bayonet-Style with Lock Screw

043-765 Phantom, SPECT, Flangeless

043-759 Phantom, SPECT, Flangeless, Small

043-772 Phantom, PET, Flangeless

Related:


043-763 Phantom, Hollow Spheres Set (6)

043-730 Phantom Insert, Triple Line

043-777 Phantom Insert, Cardiac

Phantom inserts are featured on page 99.

SPECT ANTHROPOMORPHIC TORSO PHANTOM



Frontal view – shown with optional Cardiac Insert

- Evaluation of cardiac and lung ECT data acquisition and reconstruction methods
- Evaluation of non-uniform attenuation and scatter compensation methods

The Anthropomorphic Torso Phantom is used for the evaluation of non-uniform attenuation and scatter compensation methods. The phantom consists of a large, body-shaped cylinder with lung, liver and spine inserts. The phantom simulates the anatomical structures of radioactivity distributions for the upper torso of average to large male/female patients. Lung inserts can be filled with Styrofoam® beads and water to simulate lung tissue density.

When used with the optional Cardiac Insert, cardiac ECT data acquisition and reconstruction methods may also be evaluated.

SPECIFICATIONS:

Dimensions:

O.D.: 10.25" anterior-posterior x 15" lateral (26 x 38 cm)

I.D.: 9.5" anterior-posterior x 14.2" lateral (24 x 36 cm)

Wall Thickness: 0.37" (9.5 mm)

Volumes:

Left Lung (w/o Styrofoam® beads): ~ 0.9 L

Right Lung (w/o Styrofoam® beads): ~ 1.1 L

Left Lung (w/ Styrofoam® beads): ~ 0.36 L

Right Lung (w/ Styrofoam® beads): ~ 0.44 L

Liver: ~ 1.2 L

Background: ~ 10.3 L

Cylinder w/ Lung-Spine Insert: ~ 7.4 L

043-795 Phantom, Anthropomorphic Torso

Related:

043-777 Cardiac Insert

Phantom inserts are featured on page 99.



Top view – shown with optional Cardiac Insert

SPECT LUNG-SPINE PHANTOM



Lung-Spine Phantom shown with Cardiac Insert in place and Elliptical Cylinder

- Evaluation of cardiac and lung ECT data acquisition and reconstruction methods
- Evaluation of non-uniform attenuation and scatter compensation methods

The Lung-Spine Phantom consists of two chambers that are shaped to simulate the lungs. The chamber can be filled with Styrofoam® beads and water that mimics the lung tissue. When filled with Styrofoam® beads and a radioactive solution, the lung chambers simulate lung tissue with density of ~ 0.3 gm/cm³ and with any desirable radioactivity concentration. The Lung-Spine Phantom can be used with the optional Cardiac Insert (as shown) to realistically simulate the attenuation coefficients of any radioactivity uptake in various tissue in the human upper torso. Cardiac Insert is available separately.

SPECIFICATIONS:

Inside Diameter Elliptical Shape:

Dia Along Major Axis: 12.2" (30.5 cm)

Dia Along Minor Axis: 8.7" (22.1 cm)

Inside Height: 7.3" (18.6 cm)

Volume:

Empty Cylinder: ~ 9.4 L

Left Lung (w/o Styrofoam® beads): ~ 0.9 L

Right Lung (w/o Styrofoam® beads): ~ 1.1 L

Left Lung (w/ Styrofoam® beads): ~ 0.36 L

Right Lung (w/ Styrofoam® beads): ~ 0.44 L

Shipping Weight: 11 lb (4.9 kg)

043-740 Phantom, Lung-Spine

Related:


043-763 Phantom Insert, Hollow Spheres

043-730 Triple Line Insert

043-777 Cardiac Insert

Phantom inserts are featured on page 99.

TRIPLE LINE INSERT



- *Center-of-rotation error evaluation*
- *Evaluation of changes of radius-of-rotation on spatial resolution*
- *Spatial resolution measurement in air and in water, if mounted in cylinder*
- *Quantitative evaluation of reconstruction filters and scatter compensation method*

The Triple Line Insert is used to produce three 1 mm diameter parallel lines of tracer material spaced 7.5 cm apart. The locations of the fillable tubes are based on the recommendations in the NEMA Standards Publication for Performance Measurements of Scintillation Cameras, 1986.

Radioactive tracer liquid can be inserted into the line sources through surgical grade, stainless steel valves located at the ends of each line tube.

The cylinder can be filled with water to simulate the surrounding attenuating medium.

Quantitative measurements of on-axis and off-axis reconstructed line source resolutions can be performed in air by placing the triple line insert directly on the scanning bed.

The triple line insert provides accurate, reproducible images to quantitatively evaluate the effects of errors in center-of-rotation and radius-of-rotation on scanners. Using the insert, the influence of the type of reconstruction filter on SPECT spatial resolution measurements can be evaluated.

SPECIFICATIONS:


Useful Height of Line Sources: 2.76" (7 cm)
 Diameter of Insert: 7.3" (18.6 cm)
 Diameter of Line Sources: ~1 mm
 Spacing of Line Sources: 2.95" (7.5 cm)
 Shipping Weight: 3 lb (1.3 kg)

043-730 Phantom Insert, Triple Line

Related:

043-765 Phantom, SPECT, Flangeless
043-772 Phantom, PET, SPECT, ACR
043-750 Phantom, SPECT, Deluxe
043-762 Phantom, SPECT, Standard
043-740 Phantom, Lung-Spine

CARDIAC INSERT



- *Evaluation of cardiac ECT data*
- *Evaluation of attenuation and scatter*
- *Simulates normal and abnormal myocardial uptake*
- *Solid inserts simulate transmural and non-transmural cold abnormalities*
- *Fillable inserts simulate transmural and non-transmural cold or hot abnormalities*

This insert provides a multi-function simulation of the left ventricle, and can be used to evaluate SPECT imaging of cold defects within the "myocardium." Two solid acrylic sectors (45 and 60 degrees) are supplied with the insert, each one cm thick and two cm long. These non-filling defects may be placed at various positions within the "ventricle wall", either anteriorly or posteriorly. The long axis of the "ventricle" is adjustable from 30 to 60 degrees from the long axis of the cylinder. Four fillable defects are also included.

SPECIFICATIONS:

"Ventricle" Overall Dimensions: 3.7" l x 2.4" dia (9.3 x 6.1 cm)
 "Ventricle" Volume: ~ 60 ml
 "Myocardium" Thickness: 0.4" (1.0 cm)
 "Myocardium" Volume: ~ 110 ml
 Solid Defect Set (three pieces):
 1. 60° x 2 cm (h) x 10 mm (thick)
 2. 45° x 1.53 cm (h) x 10 mm (thick)
 3. 60° x 2 cm, with 5 mm wall thickness (non-transmural defect)
 Fillable Defect Set (four pieces):
 1. 180° x 2 cm (h) x 10 mm (thick) / Vol ~ 13 ml
 2. 90° x 2 cm (h) x 10 mm (thick) / Vol ~ 5.4 ml
 3. 45° x 2 cm (h) x 10 mm (thick) / Vol ~ 3.8 ml
 4. 45° x 2 cm (h), with 5 mm thick chamber / Vol ~ 1.4 ml


043-777 Phantom Insert, Cardiac

Includes: Defect Set

Related:

043-765 Phantom, SPECT, Flangeless
043-772 Phantom, PET, SPECT, ACR
043-750 Phantom, SPECT, Deluxe
043-762 Phantom, SPECT, Standard
043-740 Phantom, Lung-Spine
043-795 Phantom, Anthropomorphic Torso

HOFFMAN 3-D BRAIN PHANTOM™



Precise simulation of activity distribution for SPECT, PET and MRI

The Hoffman 3-D Brain Phantom provides the anatomically accurate three dimensional simulation of the radioisotope distribution found in the normal brain. The Phantom allows quantitative and qualitative study of the three dimensional effects of scatter attenuation as they would appear in Iodine-123-IMP or Iodine-123-HIPDM imaging with single photon emission computer tomography or fluorine-FDG-F18 imaging with positron emission computed tomography. The phantom simulates the 4:1 uptake ratio in the gray and white matter, normal in these studies. Ventricles that are normally void of radioactivity are present.

The phantom is comprised of sturdy plastic and a single fillable chamber that eliminates the necessity of preparing different concentrations of radioisotope. Nineteen independent plates stack neatly within the cylindrical phantom for easy disassembly and assembly. The user can easily add his own custom defects to simulate clinical abnormalities.

The Phantom can be filled with the appropriate radioactive material or contrast material for SPECT, PET or MRI applications.

Each of 19 inserts is made up of five thinner slices. Two slices 0.03" thick interspersed in 0.6" thick slices to create a composite slice.

SPECIFICATIONS:

Cylinder Dimensions: I.D.: 6.9" h x 8.2" dia (17.5 x 20.8 cm)

Fillable Volume: ~ 1.2 L

Shipping Weight: 23 lb (10.4 kg)

043-790 Phantom, Hoffman 3-D Brain

NEMA PET SCATTER PHANTOM™



- Complies with NEMA 2007 Standard
- Acceptance testing
- Determine the imaging systems relative sensitivity to scatter radiation
- Measure the effects of dead-time and the effects of random events generated at different levels of activity of the line source

SPECIFICATIONS:

Cylinder Outside Dimension: 20.3 cm dia x 70 cm long

Cylinder Hole Size: 6.4 mm

Cylinder Hole Offset: 4.5 cm

Line Source Dimensions: 5 mm O.D. x 80 cm long

Line Source Inside Diameter: 3.2 mm

Shipping Weight: 52 lb (23.5 kg)

043-768 Phantom, PET Scatter, NEMA 2007

Performance Measurements of Positron Emission Tomographs, NEMA Standards Publication No. NU2, National Electrical Manufacturers Association (NEMA), Washington, DC - 2007

NEMA PET SENSITIVITY PHANTOM™



- Complies with NEMA 2007 Standard
- Ideal for PET camera sensitivity

SPECIFICATIONS:

Five internally stacked concentric aluminum tubes – all 700 mm in length

1st Tube Inside Diameter: 3.9 mm Outside Diameter: 6.4 mm

2nd Tube Inside Diameter: 7.0 mm Outside Diameter: 9.5 mm

3rd Tube Inside Diameter: 10.2 mm Outside Diameter: 12.7 mm

4th Tube Inside Diameter: 13.4 mm Outside Diameter: 15.9 mm

5th Tube Inside Diameter: 16.6 mm Outside Diameter: 19.1 mm

6th Innermost Tube (a fillable polyethylene tube)

Inside Diameter: 1 mm Outside Diameter: 3 mm

Shipping Weight: 3 lb (1.3 kg)

043-769 Phantom, PET Sensitivity, NEMA 2007

Performance Measurements of Positron Emission Tomographs, NEMA Standards Publications No. NU2, National Electrical Manufacturers Association (NEMA), Washington, DC, 2007

PET-CT PHANTOM™



PET-CT Phantom™

The PET-CT Phantom™ includes internal structures (three rods and six spheres) which, when imaged with both modalities, can demonstrate how accurately the two image sets are aligned.

In addition, a single sample of radioactive water is attenuated by water, bone and CT contrast material (as well as air only) to determine how accurately the CT-based PET attenuation correction works.

The Phantom is used for:

- Acceptance testing of PET/CT and SPECT/CT systems
- Routine quality evaluation of PET/CT and SPECT/CT systems
- Evaluation of new image fusion software
- Evaluation of new attenuation correction algorithms
- Aluminum tubes are for registration
- The outer 2" OD micro cylinder is for comparing attenuation region to non attenuation region
- The 6" ring is for contrast solution
- Research

SPECIFICATIONS:

Main Cylinder:

- Interior length of phantom: 180 mm
- Fillable spheres (5) inner diameter: 10 mm, 13 mm, 17 mm, 22 mm, and 28 mm.

Distance from sphere plane to inside wall: 70 mm

Volume of empty cylinder: 9.7 liters

Main Cylindrical insert dimension:

- Outside diameter: 51 mm
- Length: 180 mm

Top Cylinder:

- Cylinder outside diameter: ~5.1 cm
- Cylinder inside diameter: ~4. cm
- Cylinder inside height: ~8.2 cm
- Cylinder outside height: ~12.0 cm
- Volume of empty cylinder: 408 cm³

Three Aluminum Tubes:

- One 5 inch long: ~1.7 cc
- Two 7 inch (ea): ~ 2.5 cc

Stepped Bone Ring:

- Pre-filled with liquid bone composition, not to be opened
- The volumes for the bone ring are:
- Outer volume: 15.6 cu inch: ~256 cc
- Inner volume: 6.7 cu inch: ~110 cc

043-771 Phantom, PET-CT

NEMA 2007/IEC 2008 PET PHANTOM



- *Complies with NEMA 2007 Standard*
- *Simulation of whole-body imaging using PET and camera-based coincidence imaging techniques*
- *Evaluation of reconstructed image quality in whole-body PET and camera-based coincidence imaging*
- *Determination of the coincidence count rate characteristics in brain and cardiac imaging*
- *Evaluation of the relationship between true coincidence count rate and radioactivity*
- *Determination of the address errors caused by address pile up*
- *Evaluation of the count loss correction scheme*

SPECIFICATIONS:

Dimensions: 9.5" h x 12" w x 9.5 depth (24.1 x 30.5 x 24.1 cm)

Interior Length of Phantom: 180 mm

Fillable Spheres (six) Inner Diameter: 10 mm, 13 mm, 17 mm, 22 mm, 28 mm and 37 mm

Distance From Sphere Plane to Inside Wall: 70 mm

Volume of Empty D Shaped Cylinder: 9.7 L

Cylindrical Insert Dimension: O.D.: 51 mm dia x 180 mm length

Shipping Weight: 11 lb (4.9 kg)


043-767 Phantom, PET, NEMA 2007/IEC 2008

Performance Measurements of Positron Emission Tomographs, NEMA Standards Publication No. NU2, National Electrical Manufacturers Association (NEMA), Washington, DC, 2007

International Standard: Radionuclide imaging devices- Characteristics and test conditions - Part I: Positron Emission Tomographs, International Electrotechnical Commission (IEC), 61676-1, Geneva, Switzerland, 1998 and IEC 61675-1 1.1 (2008)

BAR PHANTOMS

Determines resolution of scintillation cameras



Four-quadrant Bar Phantoms offer precise determination of camera intrinsic resolution, collimator spatial resolution, field size and linearity. We offer a range of sizes manufactured to the highest quality standards.

RECTANGULAR BAR PHANTOM

Dimensions: 22.25" l x 17" w x .5" h (56.5 x 43.2 x 1.27 cm)
Lead Bar Widths: .079", .098", .118" and .138" (2, 2.5, 3 and 3.5 mm)
Field Across Bar Configurations: 21" l x 15.9" w (53.3 x 40.5 cm)
Shipping Weight: 21 lb (10 kg)

243-935 Bar Phantom, Rectangular

STANDARD BAR PHANTOM

Dimensions: 16.875" l x 16.875" w x .5" h (43 x 43 x 1.27 cm)
Lead Bar Widths: .25", .187", .156", and .375" (6.4, 4.8, 4 and 9.5 mm)
Field Across Bar Configurations: 15.875" l x 15.875" w (40.3 x 40.3 cm)
Shipping Weight: 14 lb (6 kg)

243-850 Bar Phantom, Standard

STANDARD HIGH RESOLUTION BAR PHANTOM

Dimensions: 16.875" l x 16.875" w x .5" h (43 x 43 x 1.27 cm)
Lead Bar Widths: .25", .187", .156", and .125" (6.4, 4.8, 4 and 3.2 mm)
Field Across Bar Configurations: 15.875" l x 15.875" w (40.3 x 40.3 cm)
Shipping Weight: 14 lb (6 kg)

243-800 Bar Phantom, Standard,
High Resolution

SYMBIA AND E-CAM BAR PHANTOM

Includes two removable screw knobs for insertion/removal of phantom from camera head.

Dimensions: 16" l x 21.4" w x .5" h (40.6 x 54.4 x 1.27 cm)
Lead Bar Widths: .079", .098", .118" and .138" (2, 2.5, 3 and 3.5 mm)
Field Across Bar Configurations: 20.3" l x 14.875" w (51.6 x 37.8 cm)
Shipping Weight: 19 lb (8.7 kg)

243-986 Bar Phantom, Symbia and E-Cam

CARDIAC BAR PHANTOM

Dimensions: 15.5" l x 9.25" w x .5" h (39.4 x 23.5 x 1.27 cm)
Lead Bar Widths: .079", .098", .118" and .138" (2, 2.5, 3 and 3.5 mm)
Field Across Bar Configurations: 14.5" l x 8.25" w (36.8 x 21 cm)
Shipping Weight: 15 lb (6.80 kg)

243-955 Bar Phantom, Cardiac,
High Resolution


TRIPLE HEAD RECTANGULAR BAR PHANTOM

Dimensions: 11" l x 18" w x .5" h (28 x 45.7 x 1.27 cm)
Lead Bar Widths: .079", .098", .118" and .138" (2, 2.5, 3 and 3.5 mm)
Field Across Bar Configurations: 10" l x 17" w (25.4 x 43.2 cm)
Shipping Weight: 19 lb (8.7 kg)

243-975 Bar Phantom, Rectangular, Triple Head

THYROID UPTAKE NECK PHANTOM

A must for every department performing thyroid uptakes



The Neck Phantom is designed to simulate a patient's neck. The phantom is constructed of lucite. It has a two part insert that allows counting from a bottle, vial or capsule. A capsule holder is supplied to enable the user to count capsules directly, without having to dissolve them. The phantom's cylinder and carrier have scribelines for accurate alignment. A flat surface on the cylinder allows either vertical or horizontal positioning. Twelve 30 ml bottles are included with the phantom.

Proposed by the International Atomic Energy Agency (I.A.E.A.) and the American National Standards Institute

SPECIFICATIONS:

Dimensions: 5" h x 5" dia (127 x 127 cm)

I.D.: 4" h x 2" dia (10 x 5 cm)

043-365 Thyroid Uptake Neck Phantom
Includes bottle carrier, capsule holder and 12 polyethylene bottles

Replacement:

043-361 Polyethylene Bottles, 30 ml, 50/pkg

FLOOD PHANTOMS



Determines field uniformity of scintillation cameras

Flood Phantoms provide a means of lighting scintillation camera's crystal to determine response uniformity over the entire field.

Our Flood Phantoms feature extra strength side walls and clear lucite for easy positioning. Easy to fill and easy to drain, the phantoms are leak proof and are excellent for transmission imaging.

RECTANGULAR FLOOD PHANTOM

Dimensions: 20.5" x 28" x 1.25" thick (52 x 71.1 x 3.2 cm)

Cavity: 16.5" x 24" x .5" (41.9 x 60.9 x 1.3 cm)

043-054 Flood Phantom, Rectangular




Screws directly into the flood phantom for safe, easy, fast filling

131-010 Phantom Funnel, Flood

MODEL 14C SURVEY METER WITH PANCAKE GM PROBE

Satisfies NRC requirements for nuclear medicine departments



- Monitors alpha, beta, and gamma
- 0-2 R/hr range
- External pancake GM probe (39" cable included)
- Internal energy compensated GM detector
- Five counting scales (x0.1, x1, x10, x100, x1000)

The portable Model 14C Survey Meter with Pancake GM Probe meets the essential monitoring and surveying needs of most nuclear medicine facilities. The external pancake probe is used to check hands, clothing, floors, furniture, equipment, and package surfaces for contamination. For background readings, the unit's built-in energy-compensated GM detector handles counts up to 2 R/hr. For beta emitters, the pancake probe has approximately twice the counting efficiency as an end-window detector. For high-range gamma detection, the internal detector with the x1000 multiplier range is used. The meter features a built-in speaker with ON/OFF switch, front-access calibration controls, push-button reset to zero the meter, and a toggle switch to select FAST (4 sec.) or SLOW (22 sec.) response.

SPECIFICATIONS:

Survey Meter:

Compatible Detectors: GM
 Meter Face: 0-2 mR/hr, 0-2 mR/hr, 0-6.6k cpm
 Meter Face Dimensions: 2.43" l x 1.43" w (6.1 x 3.6 cm)
 Threshold: 30 mV \pm 10 mV
 Multiplier Ranges: x0.1; x1; x10; x100 for external detector; x1000 for internal detector
 Internal Detector: Energy-compensated GM, for high range gamma detection only; 2000 mR/hr
 High Voltage: 900V
 Response: Toggle switch for FAST (4 seconds) or SLOW (22 seconds) for 90% of final reading
 Reset: Push button to zero meter after over-range exposure
 Audio: Built-in unimorph speaker with ON/OFF switch
 Sensitivity: 2100 cpm/mR/hr for Cs-137
 Batteries: Two each, size "D", typical life 600 hours
 Construction: Painted aluminum housing
 Dimensions: 8.5" l x 3.5" w x 6.5" h (22 x 9 x 16.5 cm)
 Weight: 3.5 lb (1.6 kg) including batteries

Pancake GM Probe:

Indicated use: Alpha, beta, and gamma surveying; sample counting
 Detector: Pancake type halogen quenched GM
 Window: 1.7 \pm 0.3 mg/cm² mica
 Window Diameter: 1.75" (4 cm)
 Window Area: 15.5 cm² active, 12 cm² open
 Efficiency (4pi): 5% for C-14, 22% for Sr-90/Y-90, 19% for Tc-99, 32% for P-32, 15% for Pu-239
 Gamma Sensitivity: 3300 cpm/mR/hr for Cs-137
 Energy Response: Energy dependent
 Operating Voltage: 900 volts
 Protective Screen: 79% open, stainless steel
 Housing: Painted aluminum
 Cable: 39" l (99 cm), type "C" connectors
 Dimensions: 1.8" h x 2.7" w x 10.7" l (4.6 x 6.9 x 27.2 cm)
 Weight: 1 lb (0.5 kg)

051-013 Survey Meter, Model 14C

Includes: Internal Energy-Compensated GM Detector, External Pancake GM Probe, cable and mounted check source (Cs-137, .25 μ Ci)

MODEL 14C**SURVEY METER WITH END-WINDOW GM PROBE**

Satisfies NRC requirements for nuclear medicine departments



- Monitors alpha, beta, and gamma
- 0-2 R/hr range
- External end window GM probe (39" cable included)
- Internal energy-compensated GM detector
- Five counting scales (x0.1, x1, x10, x100, x1000)
- Alpha efficiency – 15%
- Beta efficiency – 5% for C-14, 20% for Sr-90
- Gamma efficiency – 14% for Tc-99
- Sensitivity – 2100 cpm/mR/hr for Cs-137

A standard survey instrument for nuclear medicine departments, the portable Model 14C Survey Meter with End-window Probe meets the essential surveying needs of most nuclear medicine facilities. The end-window probe is used to check hands, clothing, packages, and other surfaces for contamination, and for background readings. For high-range gamma detection, the built-in energy-compensated GM detector with the x1000 multiplier range handles counts up to 2R/hr. The meter features a built-in speaker with ON/OFF switch, front-access calibration controls, push-button reset to zero the meter, and a toggle switch to select FAST (4 sec.) or SLOW (22 sec.) response.

SPECIFICATIONS:*Survey Meter:*

Compatible Detectors: GM, scintillation

Meter Face: 0-2 mR/hr, 0-4.2k cpm, battery test

Threshold: 30 mV \pm 10 mV

Multiplier Ranges: x0.1; x1; x10; x100 for external detector; x1000 for internal detector

Internal Detector: Energy-compensated GM, for high range gamma detection only; 2000 mR/hr

High Voltage: 900V

Response: Toggle switch for FAST (4 seconds) or SLOW (22 seconds) for 90% of final reading

Reset: Push button to zero meter after over-range exposure

Audio: Built-in unimorph speaker with ON/OFF switch

Sensitivity: 2100 cpm/mR/hr for Cs-137

Batteries: Two each, size "D", typical life 600 hours

Construction: Painted aluminum housing

Dimensions: 3.5" w x 6.5" h x 8.5" l (9 x 16.5 x 22 cm)

Weight: 3.5 lb (1.6 kg) including batteries

End-Window GM Probe:

Indicated Use: Alpha, beta, and gamma surveying; sample counting

Detector: End-window halogen quenched GM

Window: 1.7 \pm 0.3 mg/cm² mica

Window Area: 6 cm² active, 5 cm² open

Efficiency (4pi): 2% for C-14, 10% for Sr-90 and Y-90, 7% for Tc-99, 7% for Pu-239

Gamma Sensitivity: 2100 cpm/mR/hr for Cs-137

Energy Response: Energy dependent

Operating Voltage: 900 volts

Protective Screen: 79% open, stainless steel

Housing: Anodized aluminum

Dimensions: 1.8" diameter x 5.8" l (4.6 x 14.7 cm)

Weight: 1 lb (0.5 kg)

Detector Cable: 39" (99 cm), type "C" connectors

051-014 Survey Meter, Model 14C

Includes: Internal Energy-Compensated GM Detector, External End-Window GM Probe, cable and mounted check source (Cs-137, .25 μ Ci)

ALARM RATEMETER WITH PANCAKE GM PROBE AND FOOT MONITOR

Portable multi-purpose unit for area monitoring and surveying



- *Monitors alpha, beta, and gamma radiation*
- *Audible alarm with red indicator light for continuous area monitoring*
- *Portable surveying with Pancake GM Probe (39" cable included)*
- *0-500K cpm range*
- *Operates on wall current or rechargeable battery*
- *Equipped for simultaneous hand and shoe contamination monitoring*

This multi-function instrument can be used as a continuous area monitor, a portable survey meter, or a hand and shoe contamination monitor. It can be powered by wall current, or by its built-in rechargeable battery. Audio features include a built-in speaker with volume control providing click-per-event or an audible full-volume alarm. The alarm can be set at any point on the meter scale and locked. A red lamp on the meter's front panel also indicates an alarm condition. A push-button resets the alarm and zeros the meter. The built-in battery is continuously trickle charged when the unit is operating on line power.

The external pancake probe monitors alpha, beta, and gamma radiation, and is the preferred detector for the majority of medical surveying applications. The Pancake GM Foot Monitor serves as a convenient shoe contamination monitor. Four multiplier settings (x1, x10, x100, x1000), threshold adjustment from 10 to 100 mV, and a toggle switch for FAST (2.2 sec.) or SLOW (22 sec.) response make this unit suitable for a wide range of medical applications.

SPECIFICATIONS:

Alarm Ratemeter:

Dimensions: 8" w x 5" h x 6" l (20 x 13 x 15 cm)
Range: 0-500 cpm with multipliers of x1, x10, x100, x1000
Meter: 2.5" (6 cm) panel meter reading 0-500 cpm and BAT TEST
Response: FAST (2.2 sec.) or SLOW (22 sec.) for 90% of full scale
Input Sensitivity: Factory set at 40 mV, ± 15 mV
Linearity: Typically $\pm 2\%$ of full scale
High Voltage: Adjustable from 400-1500 volts with HV readout
Reset: Zeros meter after over range or alarm reading
Speaker: Built-in unimorph speaker
Alarm: Indicated by a red lamp on the front panel (when in an alarm state, the speaker goes to full volume overriding the volume control); non-locking alarm available upon request
Alarm Set: Located on front panel and may be checked with Test Switch; alarm reading will be indicated by audio and visual alarms
Recorder: Correlated to meter movement and is adjustable to 1.5 volts at 1 mA
Data-out Connector: 9-pin type "D" type series plug with connections as follows: Pin 1: Battery terminal; Pin 2: Unregulated supply; Pin 3: Instrument common ground; Pin 4: Alarm sink; Pin 5: Pulse out; Pin 6: Unbuffered output; Pin 7: Recorder output; Pins 8 and 9: Spare pins.
Power: 115 VAC; 220 VAC or battery power
Battery: 6 volt Gel-Cell rechargeable, with 50 hours typical life (not included)
Battery Charger: Continuously trickle charged when the power switch is in the ON position and the instrument is connected to AC power; optional external charger required for fast charge
Battery Dependence: Less than 3% calibration change when batteries read within battery test limits on meter
Finish: Polyurethane paint, beige
Weight: 4.2 lb (2 kg) including batteries

Pancake GM Probe:

Indicated use: alpha, beta, gamma surveying
Detector: Pancake type halogen quenched GM
Window: 1.7 ± 0.3 mg/cm² mica
Window Diameter: 1.75" (4 cm)
Window Area: Active area approximately 15.5 cm²
Efficiency: 10% for C-14, 45% for Sr-90, 38% for Tc-99, 70% for P-32, Alpha 30%
Sensitivity: Typically 3300 cpm/mR/hr for Cs-137 gamma
Protective Screen: 79% open, stainless steel
Mounting: Aluminum holder, handle and window protector
Detector Cable: 39" (99 cm), type "C" connectors
Dimensions: 8.5" l x 2.75" diameter (21 x 7 cm), handle diameter: 1" (3 cm)
Weight: 1 lb (0.5 kg)

Pancake GM Foot Monitor:

Indicated use: alpha, beta-gamma foot monitor
Detector: 3 ea. pancake type halogen quenched GM
Window: 1.7 ± 0.3 mg/cm² mica
Window Area: Active 46 cm², open 27 cm²
Sensitivity: Typically 9900 cpm/mR/hr
Energy Response: Energy dependent
Efficiency (2pi geometry): 10% for C-14, 45% for Sr-90, 38% for Tc-99, 65% for P-32, Alpha 30%
Housing: Aluminum
Dimensions: 10.8" h x 7" w x 9" l (27 x 18 x 23 cm)
Weight: 4.4 lb (2 kg)
Cable: 60" (152 cm), type "C" connectors ("T" connector included for simultaneous use with Pancake Probe)

051-027 Ratemeter, Alarm, 115 VAC

Includes: Pancake GM Probe and Foot Monitor

051-028 Ratemeter, Alarm, 230 VAC

Includes: Pancake GM Probe and Foot Monitor

DIGITAL AREA MONITOR

With built-in GM Gamma Detector for continuous monitoring



Digital Area Monitor - available in mR/hr or μ Sv/h

- Detects gamma radiation
- LOW and HIGH alarms with yellow and red lights, adjustable audible tones
- 0.1 mR/hr – 1000 mR/hr or 1 μ Sv/h – 9999 μ Sv/h
- Detector failure light and tone to signal detector overload or instrument failure
- Threshold adjustment from 2 to 100 mV
- Battery continuously trickle charges when unit is operating and connected to wall current
- Low battery warning

The Digital Area Monitor with built-in GM Gamma Detector provides continuous gamma radiation monitoring of rooms where radionuclides are received, stored, or dispensed, and in waste management areas where there is the possibility of radioactive contamination. The monitor is wall-mountable.

SPECIFICATIONS:

Internal Detector: Halogen quenched GM gamma detector, sensitivity: 1000 cpm/mR/hr (Cs-137 gamma), energy response (60 keV – 3 MeV): within $\pm 25\%$ of true value

Display: Four-digit LED display with 0.8" (2 cm) character height

Display Range: 000.0-9999

Display Units: Available in mR/hr or μ Sv/h

Linearity: Reading within $\pm 10\%$ of true value with detector connected

Response: Typically three seconds from 10% - 90% of final reading

Status: (Green light) instrument functioning properly

Alarm: Low alarm - indicated by yellow light and slow beep (one per second) audible tone (can be set at any point from 0.0 - 9999); high alarm - indicated by red light and fast beep (4 per second) audible tone (can be set at any point from 0.0 - 9999); NOTE: Audible indicators can be configured as a single beep if desired

Det Fail: Detector overload, no count from detector, or instrument failure - indicated by red light and audible tone (greater than 68 dB at 2 feet)

Low Bat: (Yellow) indicates less than two hours of battery power remaining

High Voltage: Adjustable from 200 - 2500 volts

Threshold: Adjustable from 2 - 100 mV

Dead Time: Adjustable to compensate for dead time of detector and electronics (can be read on display)

Overload: Senses detector saturation (indicated by display reading "-OL-")

Overrange: Indicates radiation field being measured has exceeded counting range of instrument (indicated by display reading "- - -")

Data Output: 9 pin connector providing 5 decade logarithmic output, RS-232 output, signal ground connection, FAIL and ALARM signals (current sink), and direct connection to battery and ground

Calibration Controls: Accessible from front of instrument (protective cover provided)

Power: 90-260 VAC auto ranging or battery power

Battery Life: Typically 48 hours in non-alarm condition; 12 hours in alarm condition

Battery Charger: Battery is continuously trickle charged when instrument is connected to line power and turned on

Battery Dependence: <3% change in readings to battery endpoint

Construction: Wall mount aluminum housing with ivory polyurethane enamel paint

Dimensions: 7.4" h x 9.7" w x 2.5" depth (18.7 x 24.6 x 6.4 cm)

Weight: 6.7 lb (3 kg)

051-275 Monitor, Digital Area, 90-260 VAC
*Includes: Built-in GM Gamma Detector
display in mR/hr*

051-273 Monitor, Digital Area, 90-260 VAC
*Includes: Built-in GM Gamma Detector
display in μ Sv/h*

DIGITAL ENTRANCE/EXIT AREA MONITOR

For facility entrances, waste management areas, inter-department corridors



- *Detects gamma radiation*
- *0.1 Kcpm to 4,000,000 cpm*
- *Digital display*
- *Includes two wall-mount NaI (Tl) scintillation detectors with enclosures and cables*
- *Separate LOW and HIGH alert/alarm with yellow/red lights, adjustable audible tones*
- *Detector failure light and tone to signal detector overload or instrument failure*
- *Low battery warning*
- *Threshold adjusts from 2 to 100 mV*
- *Battery continuously trickle charges when unit is operating and connected to wall current*



This Area Monitor serves as a radiation sentry/alert for building entrances, halls, and doorways. Equipped with a pair of NaI (Tl) scintillation detectors that mount on both sides of a doorway or hall, the unit continuously monitors incoming and outgoing materials and people

for low-level gamma radiation. Typical applications are the monitoring of areas where waste material leaves departments or buildings, and monitoring emergency room entrances in the event of radiation-related emergencies. The Area Monitor can be set to give a visual and/or audible alert (yellow light, slow beep) when a selected background level has been exceeded, and give an alarm signal (red light, fast beep) when a selected radiation level is reached. Lights and audibles also indicate detector overload, no count from detectors, or instrument failure. A yellow LED illuminates when the battery has less than 2 hours of power remaining. The monitor is wall-mountable. Two 25-foot detector connection cables are included.

SPECIFICATIONS:

Detectors (2): 3" dia x 1" thick (7.6 x 2.5 cm) shielded NaI (Tl) scintillation detectors in NEMA4x enclosures
Display: Four-digit LED display with 0.8" (2 cm) character height
Display Units: 0.1 Kcpm to 4,000,000 cpm
Linearity: Reading within $\pm 10\%$ of true value with detector connected
Response: Typically three seconds from 10% - 90% of final reading
Status: (green light) instrument functioning properly
Alarm: Low Alarm - indicated by yellow light and slow beep (one per second) audible tone (can be set at any point from 0.0 - 9999); High Alarm - indicated by red light and fast beep (four per second) audible tone (can be set at any point from 0.0 - 9999); **NOTE:** Audible indicators can be configured as a single beep if desired
Det Fail: Detector overload, no count from detector, or instrument failure - indicated by red light and audible tone (greater than 68 dB at 2 feet)
Low Bat: (yellow) indicates less than two hours of battery power remaining
High Voltage: Adjustable from 200 - 2500 volts
Threshold: Adjustable from 2 - 100 mV
Dead Time: Adjustable to compensate for dead time of detector and electronics (can be read on display)
Overload: Senses detector saturation (indicated by display reading "-OL-")
Overrange: Indicates radiation field being measured has exceeded counting range of instrument (indicated by display reading "---")
Data Output: Nine-pin connector providing five decade logarithmic output, RS-232 output, signal ground connection, FAIL and ALARM signals (current sink), and direct connection to battery and ground
Calibration Controls: Accessible from front of instrument (protective cover provided)
Power: 90-260 VAC auto ranging; 220 VAC or battery power
Battery Life: Typically 48 hours in non-alarm condition, 12 hours in alarm condition
Battery Charger: Battery is continuously trickle charged when instrument is connected to line power and turned on
Battery Dependence: <3% change in readings to battery endpoint
Monitor Construction: Wall mount aluminum housing with ivory polyurethane enamel paint
Dimensions:
 Monitor electronics unit: 7.4" h x 9.7" w x 2.5" depth (18.7 x 24.6 x 6.4 cm)
 Detectors (each): 13" h x 17" w x 8.5" l (33 x 43 x 22 cm)
 Connecting cables (each): 25-foot, type "C" connector
Weight:
 Monitor electronics unit: 6.5 lb (2.3 kg)
 Detectors (each): 32 lb (14.5 kg)

- 051-022** Monitor, Area, Entrance/Exit,
90-260 VAC
Includes: Two NaI (Tl) Scintillation Detectors
- 051-024** Monitor, Area, Entrance/Exit,
90-260 VAC
Includes: Two NaI (Tl) Scintillation Detectors

MONITOR 4 SURVEY METER

Pocket-sized monitor with meter, LED count light, and beeper



- Monitors alpha, beta, gamma, and x-ray radiation
- Dual-scale analog meter (CPM and mR/hr)
- Operating range: 0-50 mR/hr, 0-50,000 cpm
- Multiplier selections: $\times 1$, $\times 10$, $\times 100$
- Red LED flashes with counts
- Audible beep indicates rising radiation level
- Operates up to 2,000 hours on 9V battery
- Weighs less than 9 ounces (including battery)

The Monitor 4 Survey Meter is lightweight, ergonomically designed and includes its own carrying case. Used to measure alpha, beta, gamma and x-ray radiation, it provides extremely accurate results with both audible and visual indication. The easy-to-read scale features three ranges for maximum sensitivity and runs up to 2,000 hours on a 9V battery.

SPECIFICATIONS:

Detector: Halogen-quenched uncompensated GM tube with thin mica window 1.5-2.0 mg/cm² thick
Energy Sensitivity: Detects alpha down to 2.5 MeV; typical detection efficiency at 3.6 MeV is greater than 80%
 Detects beta at 50 keV with typical 35% efficiency
 Detects beta at 150 keV with typical 75% efficiency
 Detects gamma and x-rays down to 10 keV typical through the end window, 40 keV minimum through the case
Accuracy: $\pm 15\%$ of full scale (referenced to Cs-137)
Visual Indicators: .875" x 1.75" (2 x 4 cm) dual scale analog meter, marked 0-500 cpm and 0-0.5 mR/hr; LED count light
Audio Indicator: Internally mounted beeper (can be switched off for silent operation)
Range Switch: $\times 1$, $\times 10$, $\times 100$, battery check
Ranges: 0-500, 0-5,000, 0-50,000 cpm and 0-0.5, 0-5, 0-50 mR/hr
Current Drain: Typically 190 μ A at background radiation levels
Operating Voltage Range: 7-11 Volts DC
Voltage Regulation: High and low voltage fully regulated
Power Requirement: One 9-volt alkaline battery, NEDA #1604A, or equivalent
Battery Life: Up to 2,000 hours at normal background radiation levels
Temperature Range: -4°F to +131°F (-20°C to 55°C)
Dimensions: 8.25" x 2.75" x 1.87" (210 x 70 x 48 cm)
Weight: 9.7 oz (.28 kg) including battery, 7.2 oz (.21 kg) without battery

078-400 Survey Meter, Monitor 4

Includes: Protective case (battery not included)

INSPECTOR SURVEY METER

Low-level radiation sensitivity in a compact hand-held unit



Excellent sensitivity to low levels of alpha, beta, gamma and x-ray radiation

- Monitors alpha, beta, gamma, and x-ray radiation
- Range: 0-300,000 cpm
- 4-digit LCD readout
- Red LED count light
- Audible beeper (can be switched off)
- External calibration controls
- Adjustable timer with:
 - 1-minute sampling periods from 1 to 10 minutes
 - 10-minute sampling periods from 10 to 60 minutes
 - 1-hour sampling periods from 1 to 24 hours
- Operates 200 hours on 9V battery
- Wipe Test Plate positions wipe directly in front of detector at 1 cm distance; protects GM from damage/contamination when meter is not in use
- User selectable display for mR/hr or sieverts

The Inspector is a microprocessor controlled radiation measuring instrument which offers excellent sensitivity to low levels of alpha, beta, gamma and x-rays. The digital readout is displayed with a red count light and audible beep, providing instant indications of the radiation level. Other benefits include an adjustable timer and external calibration controls.

Shielded Holder

Maintain full function of the 078-513 Inspector Survey Meter while shielding it laterally to reduce the influence of back-ground radiation. Constructed of steel with 0.5" thick lead lining built into left and right sides, the holder conveniently supports the Inspector while allowing full operation of the optional sliding wipe test plate -- without having to lift the meter. This accessory is especially suited for use in mobile units, or where the meter is frequently used for standardized procedures in a single location. The unit is supplied with strips of reclosable fastening material for affixing to countertops or other surfaces.

SPECIFICATIONS:

Detector: Halogen-quenched uncompensated GM tube with thin mica window, 1.4-2.0 mg/cm² area density.

Window: 1.75" (45 mm) effective diameter

Display: Four-digit liquid crystal display with indicators

Averaging Periods: Display will update every three seconds; at low background levels, the update is the moving average for the past 30-second time period; time period for moving average decreases as radiation level increases

Gamma Sensitivity: 3500 cpm/mR/hr (referenced to Cs-137); smallest detectable level for I-125 is .02 µCi at contact

Operating Range:

mR/hr - .001 (1µR) to 100 mR/hr

µSv/hr - .01 to 1000

CPM - 0 to 300,000

CPS - 0 to 5000

Total/Timer - 1 to 9,999,000 counts

Accuracy (Referenced to Cs-137):

mR/hr ± 15% from 0 to 30 mR/hr, ± 20% from 30 to 50 mR/hr

CPM ± 15% from 0 to 100,000 cpm, ± 20% from 100,000 to 300,000 cpm

Timer: Can set 1-minute sampling periods from 1 to 10 minutes, 10-minute sampling periods from 10 to 60 minutes, 1-hour sampling periods from 1 to 24 hours

Count Light: Red LED flashes with each radiation event

Audio Indicator: Internally mounted beeper (can be switched off for silent operation)

Outputs: Dual miniature jack drives CMOS or TTL devices for counting to computer or data logger; (Submini jack input allows for electronic calibration)

Utility Menu: Options accessed by holding down "+" button while turning instrument on: 3-second (fast response), switch from mR/hr and CPM to µSv/hr and CPS, reset cal factor to 100, adjust probe settings, adjust the CAL factor, reset all setting to original factory setting (default)

Anti-Saturation: Meter will hold at full scale in fields as high as 100 times the maximum reading

Power Requirements: (1) 9-volt alkaline battery

Temperature Range: 14°F to 122°F (-10°C to +50°C)

Battery Life: Typically 200 hours at normal background radiation levels

Dimensions: 5.9" h x 3.3" w x 1.2" d (15 x 8 x 3 cm)

Weight: 10.6 oz (300 grams) without battery

078-514 Inspector Survey Meter Shield:

Dimensions: 10.3" l x 4.8" w x 1.5" h (26 x 12 x 4 cm)

Construction: Steel, with lead lining in left and right side walls

Lead Shielding: 0.5" thick (1.3 cm) laterally

Weight: 7.8 lb (3.5 kg)

078-510 Survey Meter, Inspector

078-513 Survey Meter, Inspector

Includes: Wipe test plate and adapter (installed)

Related:

078-514 Survey Meter Shield, Inspector

For 078-513

SURFACE SURVEY METER

Easy contamination monitoring of bench-tops, clothing, and hands



- *Monitors alpha, beta, and gamma*
- *Built-in pancake detector*
- *Range: 0 to 50,000 cpm, 0 to 15 mR/hr*
- *Three multiplier ranges: $\times 1$, $\times 10$, $\times 100$*
- *Built-in speaker*
- *Anti-saturation circuit prevents false "zero" readings*
- *One-handed operation*

This compact 3-range surface rate meter is equipped with an internal 2-inch diameter pancake GM detector and built-in speaker. The detector's thin window is recessed and protected by an aluminum grill. Its small size, light weight, and one-hand operation make this unit an ideal tool for surveying bench tops and other surfaces, and for checking hands and clothing for radioactive contamination. The meter face reads in both cpm and mR/h. Anti-saturation circuitry keeps meter needle at full scale in high radiation fields.



SPECIFICATIONS:

Meter Dial: 2.5" (6 cm) rectangular
Ranges: Three linear: 0-500; 0-5,000; 0-50,000 cpm (0-0.15; 0-1.5; 0-15 mR/h)
Switch Positions: Off, Battery Test, X100, X10, X1
Audio: Internal speaker
Detector: Halogen-quenched pancake GM tube
Diameter: 2" (5 cm)
Window Diameter: 1.75" (4.5 cm)
Window Thickness: 1.5 mg/cm²
Background: Typical 50 cpm
Efficiency: 100% for all betas and alphas that have the energy to penetrate the thin window
Gamma Sensitivity: Nominal is 3000 cpm/mR/h (based on Cs 137)
Calibration: Single master calibration potentiometer, individual potentiometers for each range
Voltage: 900V nominal
Current Drain: 3 mA typical
Power: 9-volt battery (Eveready 1222 carbon, E146X mercury, or equivalent)
Battery Life: Typically 100 hours under normal operation
Feet: Neoprene feet for placement on surface without contaminating bottom surface of unit
Handle: Swivel type, polished anodized aluminum
Dimensions: 3" w x 5.25" l x 2.25" h (8 x 13 x 6 cm) excluding meter and handle
Weight: 22 oz (625 g)

069-310 Survey Meter, Surface

DIGITAL CUTIE PIE SURVEY METER

Ion chamber detector for fast and reliable measurements of exposure and dose



- Detects alpha, beta, gamma, x-ray
- Axial detection of gamma or x-ray below 5 keV
- Rate Range: 0.1 mR/hr – 9.999 R/hr
- Dose Range: 0.01 mR – 99.9 R
- Fast, flat response – free air ion chamber
- Digital LCD display – 8-digit rate, 8-digits integrate
- Lightweight - only 26 ounces (with batteries)

Based on stable and essentially drift-free electrometer technology, this sensitive ion chamber instrument has high sensitivity for alpha and low-to-high energy beta particles, and to gamma and x-ray radiation. The compact and lightweight Digital Cutie Pie is useful for measuring exposure and dose rates, determining shielding effectiveness, checking source containers, monitoring radiation areas, and checking results following decontamination procedures. The unit contains an air ion chamber coupled to a solid-state MOSFET input electrometer with built-in A-D converter. Readout is in mR/hr or mR. Rate range is 0.1 mR/hr – 9.999 R/hr in a single range. Dose range is 0.01 mR – 99.9 R in a single range. Accurate air equivalence is assured by 180 mg/cm² graphite-lined methacrylate walls. A thin (0.5 mg/cm²) Mylar window allows high-sensitivity readings of alpha and low-energy beta such as C-14. Ion chamber survey meters are highly stable over time, making this instrument useful to reliably measure calibration sources, and radiation therapy around machines and sources. Because of their energy independent response, ion chamber survey meters are recommended for any dose rate measurements made for regulatory compliance (i.e. licensing, state regulations).

SPECIFICATIONS:

Detector: Free air ion chamber 2.5" diameter x 3.5" l (6.4 x 8.9 cm), 260 cc internal volume

Wall, Cap: Graphite-lined 180 mg/cm² walls, 540 mg/cm² cap

Window: 2.0" diameter x 0.5 mg/cm² Mylar

Readout: LCD 8-digit

Indicator Lamp: Green LED, 10 pulses/sec per mR/h

Range:

Rate: 8-digit, 0.1 mR/hr to 9.999 R/hr

Integrate: 8-digit, 0.01mR-99.9 R in a single range

Electrometer: Solid State MOSFET input

Electronics: A-D converter, LCD drivers

Batteries: 10 ea. (button) NEDA CR-1220 (7-yr shelf life), 6 ea. (AA)

NEDA 15A (typically 1000 hr)

Dimensions: 5.5" h x 3.5" w x 8" l (14 x 8.9 x 20.3 cm) includes handle

Weight: 26 oz (.74 kg) with batteries

051-366 Survey Meter, Cutie Pie, Digital,
mR/h

ELECTRONIC PERSONAL DOSIMETER



- *Ergonomic rugged design*
- *Excellent response to gamma, beta and x-ray radiation*
- *Loud, audible alarm*
- *Easy to read display*

The Electronic Personal Dosimeter combines unequalled radiological performance with advanced software and hardware features. It is suitable for use as a single stand-alone dosimeter or as part of a comprehensive dosimetry management system.

SPECIFICATIONS:

Radiation: Gamma, beta and x-ray radiation

Dimensions: 3.3" l x 2.5" w x .8" depth (85 x 63 x 19 mm), excluding clip

Direct Readout: Dose equivalents of deep dose (Hp 10) and shallow dose (Hp 0.07)

Display Resolution: 0.1 mR up to 1 R

Dose Rate Display: 0 mR/h to >400 R/h auto ranging, resolution 2 most significant digits or 0.1 mR/h deep dose (Hp 10), 1mR/h shallow dose (Hp 0.07)

Energy Response:

Deep Dose (Photon Hp 10):

- ±50% 15 keV to 17 keV (ref. ¹³⁷Cs)
- ±20% 17 keV to 1.5 MeV (ref. ¹³⁷Cs)
- ±30% 1.5 MeV to 6 MeV (ref. ¹³⁷Cs)
- ±50% 6 MeV to 10 MeV (ref. ¹³⁷Cs)

Shallow Dose (Photon Hp 0.07):

- ±30% 20 keV to 6 MeV (ref. ¹³⁷Cs)
- ±50% 6 MeV to 10 MeV (ref. ¹³⁷Cs)

Shallow Dose (Beta Hp 0.07):

- ±30% 250 keV to 1.5 MeV E (ref. ⁹⁰Sr/⁹⁰Y)

Alarm Warnings: Audible and visual alarms for dose, dose rate, count down time, read time and failure modes; high/low tone, loud/quiet volume in 7 combinations of continuous or fast/slow intermittent

Audible Alarm: Fully sealed typically 97 dB (A) at 20 cm with multiple modes

Memory: Ten-year data retention without battery; short-term dose registers for deep dose (Hp 10) and shallow dose (Hp 0.07) can be reset after each Radiation Area entry, day, month, etc. by radiation protection staff or dose management system, according to local procedures

Battery: 1.5 V AA Alkaline (typically 8 weeks continuous operation)

Operating Conditions:

Temperature: 15° F to 105° F (-10° C to + 40° C)

Relative Humidity: 20% to 90%, non-condensing

Weight: 3.2 oz (95 g)

019-016 Dosimeter, Electronic Personal, mR
Not available for export.

DIRECT READING DOSIMETERS



For personnel working in radiation areas, Direct Reading Dosimeters provide an accurate and instantaneous indication of accumulated exposure.

Accurate to 10% of True Dose for Cs-137 or Co-60 gamma. Less than 0.5% leakage of full scale in 24 hours at 50° C. Hermetically sealed.

Resembling fountain pens in size and appearance, they contain a quartz fiber electrometer and compound microscope.

The fiber is set at zero with a battery operated charger. As X- or gamma rays strike the charged Dosimeter, the fiber moves up scale in proportion to the radiation exposure.

019-201 Dosimeter, Direct Reading, 200 mR

019-500 Dosimeter, Direct Reading, 500 mR

019-010 Dosimeter, Direct Reading, 2 R

DOSIMETER CHARGER



The Dosimeter Charger is used to "zero" all Direct-Reading Dosimeters. It requires one 1 1/2 volt D Cell battery, sufficient for thousands of chargings.

020-001 Dosimeter Charger

Note: Battery not included

Radiacwash™ has been used extensively in hospitals, universities, laboratories and reactor facilities since 1951. It is the first and most popular general purpose decontamination solution specifically created for the fast and safe removal of the entire spectrum of nuclidic radioactivity.

Radiacwash is a concentrated solution designed to rapidly control radioactive contamination and remove radioactive particles from surfaces by a two-way action. First, it will sequester metallic ions which contaminate surfaces. Second, it lifts up and firmly suspends the contaminating particles, allowing contamination to be rinsed away with hard, soft or salt water.

Radiacwash will remove general laboratory contaminants such as soil, grease, oil, blood, residues, resins, and tissue and can be used

safely on all surfaces, either straight or diluted, including skin, cloth, all metals, glass, floors, walls, leather, rubber, porcelain, plastic, laboratory instruments, utensils and equipment.

Radiacwash is a synergic liquid compound that optimally combines a number of different chemical and physical principles causing it to act as a surface-wetting sequestering agent, chelater, carrier, ion-exchanger, emulsifier, solvent, complexer, peptizer and detergent.

Radiacwash has a pH of 5, less than .008% Halides, contains no phosphates, chromates, silicates, enzymes, borates, aluminates, carbonates and inert fillers that can interfere with sensitive analytical procedures.

Radiacwash is non-alkaline, non-corrosive and biodegradable.



RADIACWASH™

This compound is formulated to decontaminate the widest possible range of contamination from radioisotopes and fission products, without affecting surfaces as acids and reagents do.

005-100 Radiacwash™, bottle, 1 gal
005-155 Radiacwash™, drum, 55 gal

SPRAY MIST



Radiacwash™ Spray Mist combines the power and economy of Radiacwash with the convenience of aerosols. Radiacwash Spray Mist is a Radiacwash solution packed into a special high pressure mist applicator. When the activator is depressed, over 200 lb of pressure is created at the nozzle. The pressure misting effect allows Radiacwash to penetrate around and under contaminants and lift them off the surface into solution to be wiped up and disposed.

005-400 Radiacwash™, Spray Mist, 1 L bottle

TOWELETTES



Radiacwash Towelettes are the most efficient and safest way to remove radiocontamination from hands and small objects. Radiacwash Towelettes are individually packaged paper towels saturated in a special Radiacwash solution.

Simply remove the towelette from the foil packet, scrub hands thoroughly and wash the contamination away with running water.

005-300 Radiacwash™, Towelettes, 100/box

WIPE TEST KITS



Wipe Test Kits may be used to perform wipe tests as indicated in NRC and Agreement State Regulations. They provide an efficient, convenient means of sampling contaminated areas with radioactivity on either

wet or dry surfaces.

Each box of wipes includes 500 record folders which may be used to: (a) identify each sample wipe, (b) prevent cross-contamination of smears, (c) transfer the wipe to the counting facility, and (d) store the wipe until discarded.

Wipes are 1.75" dia cloth discs packaged in a box of five hundred.

006-350 Wipe Test Kits, 1.75" dia, 500/pkg

DECONTAMINATION KIT



The essentials for decontamination

The Atomlab Decontamination Kit contains all the equipment needed to cope with a radioactive spill or routine decontamination problem in the laboratory. The drum serves as a container for the kit components and as a waste transfer/storage vessel.

Contents of Kit:

Quantity	Item
1	30 gallon fiber drum
2 pr	Coverall, Disposable
2 pr	Shoe Covers, Disposable
2	Respirators
4	Filters
2 pr	Gloves, Reusable
1 gal	Radiacwash
1 box	Radiacwash Towelettes
1 bottle	Radiacwash Spray Mist, 1 liter
10	Poly Bags
1	12" Niptong
1 ea	Sponge, mop, scrub brush, pail, rope, assorted signs

SPECIFICATIONS:

Drum Dimensions: 20" d x 29.5" h (50.8 x 74.9 cm)
Shipping Weight: 40 lb (18.5 kg)

121-180 Decontamination Kit

MINOR EMERGENCY SPILL KIT



The Emergency Spill Kit is based on suggested contents described in NRC Regulations Guide 10.8, Revision 2: Medical Use Programs, Appendix J: Medical Spill Procedures. Designed to handle minor radioactive spills and routine contamination problems in the lab, it also includes the emergency procedures and forms to document the spill and decontamination clean up.

Contents of Kit:

Quantity	Item
6 pr	Gloves, Disposable
1 pr	Gloves, Anti-C/housekeeping
2	Lab Coats, Disposable
2	Head Covers, Disposable
2	Shoe Covers, Disposable
2	Mess Up Mitts
6	Plastic Trash Bags with Twist Ties
1 roll	Tape ("Caution Radioactive Material")
3	Pre-strung Tags ("Caution Radioactive Material")
1 pkg	Contamination Wipe Sample Supplies
1	Clipboard with "Emergency Procedures" instructions and forms

SPECIFICATIONS:

Shipping Weight: 6 lb (2.72 kg)

121-190 Emergency Spill Kit, Minor

COLLIMATOR PROTECTORS



Shown with Wall Dispenser, 114-017

Disposable plastic Collimator Protectors are designed to keep the face of the scintillation camera collimators from becoming contaminated by infectious fluids.

Protector sheets are backed with an adhesive for quick and easy application. After use,

peel away the contaminated sheet and discard. Apply a new protector between each patient.

Protectors are packaged in perforated rolls. The optional dispenser is designed of heavy gauge polished stainless steel with an aluminum rod. Mount to any wall or door for easy access.

SPECIFICATIONS:

Dimensions: 17" x 17" (43 x 43 cm)


Weight: 9 lb (4 kg)

114-177 Collimator Protectors, 250/roll

Related:

114-017 Wall Dispenser, Collimator Protectors
Wall mounted

NIPTONGS



These low-cost tongs are used to handle small radioactive or otherwise dangerous objects up to 1" diameter. The tongs have a 45° v cut groove on each jaw. The compression spring maintains a strong grip

on the object until the tension is released by squeezing the finger bar. Niptongs are made of chrome-plated, high carbon steel with hardwood handles and are easily disassembled for decontamination and cleaning.

011-012 Niptongs, 12" (30 cm)

011-036 Niptongs, 36" (91 cm)

ABSORBENT PAPER



Highly absorbent with plastic lining for fast, safe clean up

Protect any work surface with plastic lined Absorbent Paper. Spills and splashes are immediately contained without messy clean up. Available in pre-cut sheets or by the roll, Absorbent Paper is efficient for any work surface.

Together with the optional paper dispenser, the 300-foot Absorbent Paper Roll can be conveniently located to any wall or bench. It works like a common paper towel dispenser, yet is extra-strong to support the oversized 11 lb roll.

SPECIFICATIONS:

033-300 Absorbent Paper Roll

Dimensions: 300' l x 20" w (91.5 m x 51 cm)

033-013 Absorbent Paper Sheets


Dimensions: 17" l x 13" w (43 x 33 cm)

033-300 Absorbent Paper, Roll, 300 ft

033-013 Absorbent Paper, Sheets, 50/pkg

Related:

046-275 Dispenser, Absorbent Paper Roll
Wall mounted



Radiacwash™ has been used extensively in hospitals, universities, laboratories and reactor facilities since 1951. It is the first and most popular general purpose decontamination solution specifically created for the fast and safe removal of the entire spectrum of nuclidic radioactivity.

See page 114 for details.

IMAGING CHAIR

Fully adjustable to accommodate all imaging procedures



Chair back and arms fold away from imaging area with the touch of a button




Armrests and seat back easily drop away from procedure area

For most imaging procedures, the patient lies on a table, which is then placed in the proper position. In some cases, the table restricts accessibility of the imaging equipment and patient comfort. The Imaging Chair puts the patient at ease, and provides alternative patient positioning for nuclear medicine procedures.

The Imaging Chair is fully adjustable. With the back support and armrests in place, it is a chair, which can swivel 360°. Once wheeled to the proper location, the chair can be anchored while the patient is seated by engaging the foot-operated wheel locks on the four casters. The height of the chair can be increased from 41" to 48" by stepping on the foot pump arm. Hold the arm down and the chair slowly lowers to its minimum height. If you raise the foot pump arm, it locks into that position, which restricts the swivel action of the chair. At the rear of the seat cushion is a locking handle which allows you to move the seat of the chair 4" to the left or right, while the base remains anchored. For lateral and posterior access, the back support and armrests drop out of the way with a simple tug of the corresponding pull pin.

The seat and back are upholstered with 2" foam and vinyl cover. The wide armrest allows the Imaging Chair to double as an injection chair.



SPECIFICATIONS:

Dimensions: 24" w x 41" h, (64 x 104 cm)
 Seat Dimensions: 17" w x 11" deep (43 x 28 cm)
 Height Adjustable: 19" h to 26" h (48 to 66 cm)
 Backrest: 19" above seat top (48 cm)
 Upholstery: 2" foam with gray vinyl cover
 Finish: Black baked enamel
 Wheels: Four locking swivel casters, 2" dia (5 cm)
 Patient Capacity: 350 lb (159 kg)
 Weight: 75 lb (34 kg)
 Shipping Weight: 85 lb (38.5 kg)

214-610 Chair, Imaging

EASY MOVER

One person can transfer patient from table to table



The Easy Mover is a high-density polyethylene patient transport board that offers high-impact strength, excellent abrasion resistance, tensile strength and resistance to stress cracks.

Semi-rigid and radiolucent, the Easy Mover allows the patient to remain on the board during imaging procedures. Often only one person is required to complete a move. Four strong straps, two on each side, are positioned for gripping and sliding the mover while the patient is on it. The straps are easily replaced for sanitation.

The Easy Mover is durable, light-weight and can be stored away when not in use.



SPECIFICATIONS:

Dimensions: 23" w x 73 1/2" x .19" thick (58 x 183 x .48 cm)
Weight: 2 lb (1 kg)

056-352 Easy Mover

Replacement:

056-350 Straps, Handle, 4/set

HIGH HANDLE FOOT STOOL



Non-skid pad for patient safety

This lightweight, yet durable stool is the safe and easy way for patients to get on or off examination tables. Steel construction assures long life. Rubber caps on the legs and a non-slip surface further assure stability and patient safety.

SPECIFICATIONS:

Shipping Weight: 12 lb (5.5 kg)
Patient Capacity: 250 lb (113.4 kg)

214-728 Foot Stool, High Handle

ADJUSTABLE HEIGHT STOOLS

Preferential seating



This Adjustable Height Stool is made of a 1" diameter chrome plated steel. For versatility within the department, the stool features precision machined screw height adjustment from 19" to 27".

The cushioned 15" diameter seat is a plush 4" thick poly-foam for maximum comfort. The vinyl upholstery is extremely durable and easy to maintain. Beneath the seat is a protective

ABS plastic shroud for asepsis.

The 4-legged stool sits atop a 19.5" diameter base with casters or rubber "tips" for added stability. Foot ring is standard.

SPECIFICATIONS:

Dimensions: 19.5" dia x 27" h (48 x 69 cm)
Height Adjustable: 19" to 27" (48 to 69 cm)
Seat Dimensions: 15" dia (38 cm)
Construction: Chrome plated steel
Weight: 18 lb (8.3 kg)

214-130 Stool, Adjustable Height with Casters

214-132 Stool, Adjustable Height with Rubber Tips

ARM SUPPORT STRAP



- *Reduces patient fatigue*
- *Holds patients still during long procedures*

Patient fatigue can often result in procedure interruption. The Arm Support Strap is designed to reduce patient fatigue by immobilizing and supporting the arms at patient's side. This is especially important for procedures performed on narrow tables, like SPECT, CT or MRI.

The Arm Support Strap is a comfortable cotton with fastening straps that wrap around both the table and the patient. Its construction offers no attenuation to the image.

Wide enough to comfortably support the weight, the strap measures 8.5" wide and comes in an easy to identify standard yellow and orange so that it is easy to find within the department.

SPECIFICATIONS:

Dimensions: 71" l x 8.5" w (180 x 21.6 cm)

056-786 Arm Support Strap

TABLE PAD



Virtually indestructible, the Table Pad is ideal for all patient handling applications. The Table Pad is a high density 1" polyurethane foam that ensures patient comfort by insulating from cold, hard table surfaces. The non-stick cover is made of "Gortex"-type material

which resists liquid penetration for easy stain removal. The six double-stitched strap handles and the slick "Gortex"-type cover permit easy one-person patient handling on the x-ray table and also make transfers from bed to gurney to table easier on patients as well as personnel.

Nylon hook eyes are positioned for quick and easy hanging storage to keep the pad clean and out of the way when not in use.

SPECIFICATIONS:

Dimensions: 25" w x 81" l (63.5 x 205.7 cm)

Construction: 1" thick polyurethane foam

056-358 Mattress, Table Pad, 1" thick

ARM SUPPORT




Patient comfort is essential, especially during long nuclear, CT or PET procedures that can be costly if interrupted. The Arm Support is a proven

design that provides anatomical comfort with correct arm and shoulder support. Accommodates single or dual head imaging systems.

056-791 Arm Support

For single and dual head imaging systems

SPECT ARM BOARDS



- *Enhanced patient comfort*
- *Stable arm support*

SPECT Arm Boards are designed to fit on all SPECT tables. They are positioned under the patient while lying on the table to effectively widen the surface. The patient's arms can then be held comfortably next to the body while imaging the head. The amount of overhang on the side of the table is adjusted by sliding the arm board further under the patient or pulling them further out for a wider surface.

The arm boards have been built with a convenient bend in them to aid in patient comfort and positioning.

SPECIFICATIONS:

Dimensions: 20" l x .25" thick (50.8 x .6 cm)

Width: Tapered 13.75" to 10.75" (34.9 to 27.3 cm)

Construction: Plexiglas


Weight: 4 lb (2 kg)

056-767 Arm Boards, SPECT, pair

POSITIONING WEDGES

Constructed of fine-cell, charcoal grey polyurethane foam, the positioners include a vinyl cover. Designed to position extremities, they can accommodate almost every conceivable position needed, including skulls, hips and spines. Vinyl covers have a smooth surface that is easy to clean without the crevices from stitched seams. Shipping weight for all positioners is 1 lb (.5 kg).


RECTANGULAR BLOCK



120-253 Positioner, Vinyl Covered,
Rectangular
8" w x 12" l x 3" h


BOLSTER BLOCK

Stabilizes patient for myelogram, needle insertion.




120-274 Positioner, Vinyl
Covered, Bolster
9" w x 16" l x 6.625" h

TRIANGLE



120-280 Positioner, Vinyl Covered,
Triangle
10" w x 10" l x 10" h


DECUBITUS BLOCK



120-251 Positioner, Vinyl Covered,
Decubitus
18" w x 24" l x 4" h


27° TORSO BLOCK

For dorsal and lumbar spine




120-250 Positioner, Vinyl Covered,
27° Torso
16" w x 12" l x 7" h

45° SPINAL WEDGE



120-265 Positioner, Vinyl Covered,
45° Spinal
7" w x 21.25" l x 7" h

45° SPINAL WEDGE



120-291 Positioner, Vinyl Covered,
45° Spinal Wedge
9" w x 30" l x 9" h



Sof-Skin Coat Apron
See page 32 for details.

PRESSURE SENSITIVE SHIPPING LABELS




- *Pressure Sensitive*
- *Peel-off backing*
- *Package of 20*

Comply with Federal Regulations, Title 49 - Transportation of Hazardous Materials Regulations, Dept. of Transportation, 49CFR, Part 173 as to wording, symbols, size and colors.

Labels, Shipping (Pressure Sensitive):

- 030-001** Label, Radioactive I, White, 20/pkg
030-002 Label, Radioactive II, Yellow, 20/pkg
030-003 Label, Radioactive III, Yellow, 20/pkg
 4.25" w x 4.25" h (10.8 x 10.8 cm)

PRESSURE SENSITIVE WARNING LABELS



028-002

Pre-cut labels are mounted on a paper backing. Simply peel off the pre-cut label and apply.

Labels, Warning (Pressure Sensitive):

- 028-002** Label, Caution, Radioactive Material, 20/pkg
 5" w x 6" h (12.7 x 15.2 cm)
026-015 Label, Radioactive Material, 500/roll
 3" w x 1" h (7.6 x 2.5 cm)
026-014 Label, Caution, Radioactive Material, 500/roll
 2" w x 1" h (5 x 2.5 cm)

PERFORATED WARNING TAPES



026-013

026-012



026-011



026-005


Pressure Sensitive adhesive plastic tapes stick to almost everything and are moisture resistant. All wording, symbols and colors conform to government regulations. Use ball-point pen on tape labels for proper container identification. Continuous 180 ft. rolls.

Perforated Warning Tape:

- 026-013** Tape, Radioactive Material, 180 ft/roll
 1" w x 3" l (2.5 x 8 cm)
026-012 Tape, Caution, Radioactive Materials, 180 ft/roll
 1" w x 2" l (2.5 x 5 cm)
026-011 Tape, Biohazard, Hazard Identity, 180 ft/roll
 1" w x 3" l (2.5 x 8 cm)
026-005 Tape, Caution, Radioactive Materials, 180 ft/roll
 .5" w 1.38" l (1.3 x 3.5 cm)
Suitable for test tubes, bottle necks, etc.

PRE-STRUNG TAG

Extra strength cardboard with reinforced string hole



029-001


- 029-001** Tags, Container, 100/pkg
 2.63" w x 5.25" h (7 x 13 cm)

RADIOACTIVE WARNING LABELS

Protects staff and visitors

CAUTION: RADIOACTIVE MATERIAL

PATIENT _____
 ISOTOPE _____
 ACTIVITY _____ mBq _____ mCi
 DATE _____




026-106

 **RADIOACTIVE WASTE
STORE FOR DECAY**

ISOTOPE _____ DATE IN _____ INITIALS _____

DATE	SURVEY METER	INITIALS
_____	_____	_____
_____	_____	_____
_____	_____	_____

026-109 For documenting periodic checks of radioactive waste

 **RADIOACTIVE
WASTE**

STORE FOR DECAY

DATE IN _____

INITIALS _____

END STORAGE DATE: _____

026-108

Labels, Radioactive Warning:

026-106 Label, Caution, Radioactive Material, 320/roll

3" w x .875" h (7.6 x 2.2 cm)

026-109 Label, Radioactive Waste, 240/roll

4" w x 2.63" h (10.2 x 6.7 cm)

026-108 Label, Radioactive Waste, 320/roll

2" w x 3" h (5 x 7.6 cm)

WARNING ROPE



- *Brightly colored - magenta and yellow*
- *Triple-strand strength*

Made of high quality polypropylene, this triple strand rope assures high visibility and remarkable strength.

121-073 Rope, Warning, 100'/pkg

.31" dia (.8 cm)

Conveniently packaged in 100' (31 meters) lengths

HAZARD TAPE



An inexpensive compliance with OSHA's code for nuclear hazard


Hazard Tape is a highly visible magenta and yellow warning tape that's perfect for identifying physical hazard areas and materials. Constructed of 6 mil. vinyl, the tape is durable to withstand long-term placement. The tape will easily affix to any clean, smooth, dry surface.

026-020 Tape, Hazard


2" w x 54' l (5 cm x 16.5 m)

NUCLEAR MEDICINE LABELS


These labels are used for patients who have received gamma-emitting radionuclides. Printed in magenta ink on yellow background, labels are backed with a special adhesive for temporary adherence, without leaving a sticky residue. Copy recommended by National Council on Radiation Protection and Measurements.




026-100




026-101



129-160



129-161



129-162

- 026-100** Label, Radionuclide Administered, 200/roll
4" w x 2.63" h (10.2 x 6.7 cm)
- 026-101** Label, Radionuclide Administered, 200/roll
4" w x 2.63" h (10.2 x 6.7 cm)
- 129-160** Label, Patient Injected With, 320/roll
3" w x 2" h (7.6 x 5 cm)
- 129-161** Label, Patient Injected For, 320/roll
3" w x 2" h (7.6 x 5 cm)
- 129-162** Label, Radionuclide Administered, 320/roll
3" w x 2" h (7.6 x 5 cm)

PLASTIC CAUTION SIGNS

Convey the warning clearly


These signs are constructed of heavy-duty polyethylene with convenient holes for indoor or outdoor mounting.



024-914



024-970



024-957



024-906



024-948



024-581



024-999



024-900




024-923



024-910



024-911



024-916

Signs, Plastic Caution:

- 024-914** Sign, Caution, Radiation Area
14" w x 10" h (35.6 x 25.4 cm)
- 024-970** Sign, Caution, Radiation Area
14" w x 10" h (35.6 x 25.4 cm)
- 024-957** Sign, Radiation Monitors Required
14" w x 10" h (35.6 x 25.4 cm)
- 024-906** Sign, Caution, High Radiation Area
14" w x 10" h (35.6 x 25.4 cm)
- 024-948** Sign, Caution, Radioactive Materials
14" w x 10" h (35.6 x 25.4 cm)
- 024-916** Sign, Caution, If you are Breastfeeding
14" w x 10" h (35.6 x 25.4 cm)
- 024-999** Sign, Caution, Radioactive Materials
7" w x 10" h (17.8 x 25.4 cm)
- 024-900** Sign, Caution, Radiation Area
7" w x 10" h (17.8 x 25.4 cm)
- 024-923** Sign, Caution, X-Ray Caution
7" w x 10" h (17.8 x 25.4 cm)
- 024-581** Sign, Danger, Radioactive Materials
10" w x 7" h (25.4 x 17.8 cm)
- 024-910** Sign, Caution, If you are Pregnant
14" w x 10" h (35.6 x 25.4 cm)
- 024-911** Sign, Precaucion, Si usted esta Encinta
Spanish version of 024-910
14" w x 10" h (35.6 x 25.4 cm)

001-001	30	001-881	49	007-957	10	033-300	116	043-767	101	066-536	30	127-999	4	150-780	60	212-026	44
001-075	42	002-246	23	007-973	8	039-106	27	043-768	100	069-310	111	129-160	123	150-781	61	212-027	44
001-100	24	002-248	23	007-974	8-9	039-110	27	043-769	100	075-594	69, 73	129-161	123	150-782	61	212-028	44
001-130	24	005-100	114	007-975	8	039-244	59	043-771	101	075-596	71, 73	129-162	123	150-783	61	212-029	44
001-132	24	005-155	114	007-978	58	039-284	59	043-772	97	078-400	109	130-019	90	150-785	60	212-030	44
001-135	24	005-300	114	007-980	8	039-286	59	043-777	99	078-510	110	130-020	90	150-951	61	212-031	44
001-140	24	005-400	114	007-983	8	039-325	17	043-790	100	078-513	110	130-100	89	150-952	61	212-032	44
001-142	24	006-350	114	007-984	57	039-326	18	043-795	98	078-514	110	130-551	89	150-960	61	214-129	63
001-145	24	007-007	16	007-985	8	039-330	16	043-830	94	086-241	75	130-639	90	150-961	61	214-130	118
001-147	24	007-040	13	007-986	57	039-335	17	043-840	94	086-242	74	130-900	91	150-971	61	214-132	118
001-148	24	007-241	7	007-988	57	039-338	16-18	043-845	94	086-243	75	130-901	91	150-991	61	214-200	14
001-149	24	007-242	7	007-989	57	039-341	17	043-855	94	086-301	75	130-902	91	154-090	30	214-210	15
001-179	12	007-243	7	007-990	8	039-350	17	043-860	94	086-330	66-69	130-903	91	177-075	81	214-220	14
001-180	12	007-244	7	007-995	8	039-387	16-18	046-275	116	086-331	70-71	130-904	91	177-090	80-81	214-610	117
001-181	12	007-245	7	007-997	13	039-412	18	050-200	16	086-332	72-73	131-010	103	177-091	81	214-728	118
001-182	12	007-246	7	007-999	11	039-413	18	050-205	16	086-333	69	132-319	85	177-092	81	243-800	102
001-200	24	007-275	15	008-015	75	039-514	59	050-250	16	086-334	65	132-503	84-85	177-095	82	243-850	102
001-236	43	007-276	158	008-300	11	039-599	59	051-013	104	086-335	64-65	132-555	85	177-124	82	243-935	102
001-280	46	007-612	6	008-400	11	040-108	27	051-014	105	086-336	69	132-680	86	177-324	82	243-955	102
001-281	46	007-620	6	009-205	11	040-315	16	051-022	108	086-338	75	132-681	86	177-325	82	243-975	102
001-282	46	007-635	6	011-012	30	042-224	22	051-024	108	086-339	69	132-684	86	185-302	90	243-986	102
001-283	47	007-652	6	011-036	30	042-226	22	051-027	106	086-341	69	132-685	88	187-130	76-79	244-003	39
001-284	47	007-661	42	019-010	113	042-227	22	051-028	106	086-342	71	132-690	86	187-135	76-79	244-006	26
001-300	24	007-663	42	019-016	113	042-228	22	051-273	107	086-343	69	132-691	86	187-140	76-79	244-007	35
001-400	24	007-665	42	019-201	113	042-407	25	051-275	107	086-344	69	132-692	87	187-145	76-79	244-009	39
001-500	24	007-670	6	019-500	113	042-413	19	051-366	112	086-423	74	132-695	88	187-246	77	244-110	34
001-600	24	007-671	6	020-001	113	042-417	25	053-607	43	086-435	74	132-698	88	187-256	77	244-111	34
001-706	43, 54	007-675	6	024-581	123	042-419	21	053-805	43	086-509	74	132-699	87	190-210	26	244-112	34
001-707	54	007-680	6-7	024-900	123	042-425	25	055-903	23	087-112	26, 39	132-750	89	190-215	62-63	244-120	38
001-708	53	007-685	6	024-906	123	042-426	25	056-350	118	099-289	93	132-770	87	190-217	63	244-121	38
001-721	48-55	007-691	42	024-910	123	042-428	21	056-352	118	099-291	93	132-771	87	190-218	63	244-122	38
001-723	53	007-693	42	024-911	123	042-433	20	056-358	119	101-103	93	132-772	90	190-219	63	244-130	37
001-724	54	007-695	42	024-914	123	042-434	25	056-767	119	101-356	95	132-773	90	190-220	63	244-131	37
001-726	45	007-700	4	024-916	123	042-435	25	056-786	119	103-701	32	132-774	87	208-065	31	244-132	37
001-730	50	007-701	4	024-923	123	042-448	41	056-791	119	104-030	23	132-781	86	208-070	31	244-140	34
001-739	52	007-711	9	024-948	123	042-449	21	060-133	92	112-036	26	132-784	86	208-075	31	244-141	34
001-754	47	007-712	9	024-957	123	042-455	56-57	060-137	92	114-017	116	132-793	87	208-080	31	244-142	34
001-756	47	007-713	9	024-970	123	042-456	57	060-139	91	114-177	116	134-772	87	208-085	31	244-150	36
001-757	47	007-716	9	024-999	123	042-460	58	063-100	93	116-010	32	135-022	14	208-090	31	244-151	36
001-758	47	007-717	9	026-005	121	042-466	56-57	063-137	93	117-425	32	136-755	85	212-001	44	244-152	36
001-759	47	007-718	9	026-011	121	042-467	56-57	063-138	93	120-250	120	139-101	90	212-002	44	244-160	36
001-771	43	007-723	5	026-012	121	042-519	28	063-139	93	120-251	120	139-102	90	212-003	44	244-161	36
001-779	47	007-734	5	026-013	121	042-522	28	063-140	93	120-253	120	139-104	90	212-004	44	244-162	36
001-785	50-52	007-735	5	026-014	121	042-580	29	063-261	95	120-265	120	139-677	88	212-005	44	244-170	37
001-786	50	007-736	5	026-015	121	042-582	29	063-350	95	120-274	120	139-680	90	212-006	44	244-171	37
001-787	51	007-738	5	026-020	122	042-710	23	063-360	95	120-280	120	150-000	60	212-007	44	244-172	37
001-788	45	007-747	5	026-100	123	043-054	103	063-361	95	120-291	120	150-001	61	212-008	44	244-181	38
001-789	45	007-755	5, 7	026-101	123	043-274	93	063-362	95	121-073	122	150-005	61	212-009	44	244-182	38
001-790	45	007-800	4	026-106	122	043-361	103	063-363	95	121-180	115	150-025	61	212-011	44	244-190	35
001-793	55	007-801	4	026-108	122	043-365	103	063-364	95	121-190	115	150-122	61	212-012	44	244-191	35
001-794	55	007-900	4	026-109	122	043-730	99	063-562	95	123-500	93	150-125	61	212-013	44	244-192	35
001-797	55	007-901	4	028-002	121	043-740	98	063-586	95	124-230	94	150-130	61	212-014	44	244-200	40
001-850	58	007-909	4	029-001	121	043-750	96	063-700	94	124-231	94	150-160	61	212-015	44	244-205	40
001-855	57	007-911	4	030-001	121	043-759	97	063-701	94	127-691	42	150-315	85	212-016	44		
001-860	47-48	007-913	4	030-002	121	043-762	96	063-720	95	127-693	42	150-771	61	212-017	44		
001-861	48	007-945	8	030-003	121	043-763	96	066-533	30	127-700	4	150-772	61	212-019	44		
001-880	48	007-956	10	033-013	116	043-765	97	066-535	30	127-789	4	150-773	61	212-023	44		

A bsorbent Paper	116	D ecay and Storage Cabinet	36	Fume Hood, Radioiodine	26
Aerotech I Radioaerosol System	82	Decay Cabinet	34	G enerator and Storage Cabinet	38
Aluminum Breakthru Kit	60	Decay Drums, Shielded	59	Gloves, Lead Neoprene	32
Apron, Radiation Protection	32	D econtamination and Sampling	114-116	Gloves, Radiation Attenuating	31
Arm Support, Imaging Tables	119	Decontamination Kit	115	Goggles, Lead Glass	32
Atomlab Dose Calibrators	64-73	Dippers	74-75	GRX Rapid Exchange Glove System	63
Atomlab Thyroid Uptake Systems	76-79	Disposables, Lung Ventilation	80-92		
Atomlab Wipe Test Counter	70-71	Document Protector	45, 48-55		
B acteria Filter, Disposable	89	D ose Calibration	64-75	H azard Tape	122
Bar Phantoms	102	Atomlab 400	64-65	Hoffman 3-D Brain Phantom	100
Barrier, Clear-Lead, Mobile	29	Atomlab 500	66-69	Hypalon Sleeves	63
Bench Top Shield	23	Atomlab 500 Plus	72-73		
B EXXAR Accessories:		Atomlab Wipe Test Counter	70-71	I -131 Vial Pig Shipping System	49
Dose Preparation Stand	58	Chamber Dose Calibrator	69, 71, 73	Imaging Accessories	117-120
Dose Preparation Syringe Shield	58	Chamber, Wipe Test	69, 73	Imaging Chair	117
Dose Preparation System	58	Lead Shield, Wipe Test Chamber	71, 73	Injection Chairs	14-15
Dose Preparation Vial Shield	58	Shielding Rings	75	Injection Pump Shield, Graseby	13
Graseby Injection Pump Shield	13	Well Insert	75	Injection Stand	14
Syringe Dipper	75	Dose Drawing System, PET	56-57	Injector, Dose, Manual	13
Unit Dose Pig	47	Dose Drawing Vial Shield, PET	57	Intego PET Vial Shipping System	53
Unit Dose Shipping System	48	Dose Injector, Manual	13	Interlocking Lead Bricks	24
Bricks, Lead	23-25	Dose Preparation Stand	56-57	Isolator, Shielded	62-63
C abinet, Storage	30	Dose Preparation System, BEXXAR	58		
Cart, Transport, PET Shipping System	50	Dosimeters	113	L abels:	
Caves, Brick Lead	25	Drierite and Soda-Lime	90	Dose Calibrator	67-69
Chair, Imaging	117	Drums, Decay	59	Nuclear Medicine	123
Chair, Injection	14-15	E asy Mover	118	Shipping	121
Chromatography	60-61	Emergency Spill Kit	115	Syringe, Color-Coded	44
Clear-Lead Mobile Barriers	29	F ace Mask Harness	88	Warning	121-122
Collimator Protectors	116	Face Masks	88	L-Block Shields, PET	19-21
Container, Rectangular, Shielded	30	Filter, Bacteria, Disposable	89	Lead Bricks	23
Container, Shipping, Zevalin	45	Flood Phantoms	103	Lead Foil	23
Containers, Sharps	16-18	Flood Sources	94	Lead Vinyl Sheets	23
Containers, Storage	16	Foil, Lead	23	Leaded Neoprene Gloves	32
Cyclotron Workbench	57	Foot Stool, High Handle	118	L ead-Lined Laboratory Furniture	33-41
		Forceps	30	Decay and Storage Cabinet	36
				Decay Cabinet	34
				Generator and Storage Cabinet	38

CONDITIONS OF PURCHASE

CREDIT TERMS: Net 30 days, available to qualified accounts only. Accounts not paid within 30 day period are subject to a 1.5% charge for each additional month.

PRICE: All prices in this catalog are F.O.B. Shipping Point and are subject to change without notice.

INTERNATIONAL CUSTOMERS: Prices in this catalog are for U.S. customers only. Prices do not include duty, shipping, VAT, international packing, international labeling, language options, warranty, installation, service or support. Biodef products are sold and supported by a worldwide network of qualified distributors. Visit our web site for the name of your local distributor, www.biodef.com.

MINIMUM ORDER: \$50.00

METHOD OF PAYMENT: Visa, MasterCard, American Express and personal, prepaid checks accepted. 

SHIPPING POLICY

SHIPPING: Shipment of orders is handled by a computerized system to assure the lowest cost possible. Under 70 lb - UPS. Over 70 lb - Motor Freight. (All motor freight charges are collect unless otherwise stated). All orders will be charged actual shipping fees. When ordering by prepaid check, call to confirm shipping charges and applicable tax.

RUSH ORDERS: Emergency orders can be shipped next-day air or second-day air, incurring additional freight charges.

DAMAGE IN SHIPMENT: At time of shipment, all goods are inspected and warranted to be in first-rate condition. Our packaging of all products complies with, and in most cases, exceeds requirements of Motor Freight Classifications. Despite the care used in packaging, damages occasionally happen in transportation. Upon receipt of shipment, open immediately and inspect for concealed damages. If damage is sighted, stop unpacking. Hold shipment container and ask delivering carrier for an inspection for concealed damages. This notification to the carrier must be made within 10 days of receipt of shipment for claim purposes.

CLAIMS: Any claim for damages should be made against the transportation company. Our responsibility ceases upon receipt of shipment by the carrier, but if we can be of any assistance in processing a claim, please advise us by calling 1-800-224-6339.

SATISFACTION GUARANTEED

EXCHANGES AND RETURNS: Thank you for choosing Biodef. If you are not satisfied with your order, it can be returned for credit or exchanged for a more suitable product. The customer is responsible for return freight charges and for certain transactions, subject to a restocking charge. In the event of an incorrect product shipment, due to Biodef's error, Biodef is responsible for return freight charges. Products returned must be unused, in resalable condition and in original packaging. For complete terms of Biodef Exchange and Return Policy, www.biodef.com/conditions.

LIMITED WARRANTY: Biodef Medical Systems, Inc., (Biodef), hereby warrants its products sold to buyer to be free of defects in workmanship and material for a period of one year from date of shipment, unless stated otherwise in the specifications of the specific product. THE FOREGOING WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES EXPRESSED, STATUTORY OR IMPLIED. INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. Final determination of the suitability of the products for the use contemplated by Buyer is the sole responsibility of Buyer, and Biodef shall in no way be responsible for the suitability of the Products for any particular use. Claims for defects in workmanship or material must be presented to Biodef in writing and received by Biodef prior to the expiration of the warranty period. BUYER'S SOLE AND EXCLUSIVE REMEDY, AND BIODEx'S SOLE OBLIGATION, UNDER THIS WARRANTY SHALL BE REPAIR OR REPLACEMENT OF THE DEFECTIVE PRODUCT AS DETERMINED AT BIODEx'S SOLE OPTION. IT IS UNDERSTOOD THAT UPON REPAIR OR REPLACEMENT, BIODEx SHALL HAVE NO FURTHER OBLIGATION TO BUYER WITH RESPECT TO THE DEFECTIVE PRODUCT. UNDER NO CIRCUMSTANCES SHALL BIODEx BE LIABLE TO BUYER OR ANY THIRD PARTY FOR SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, LOSSES OR EXPENSES IN CONNECTION WITH, OR BY REASON OF THE INABILITY TO USE THE PRODUCTS FOR ANY PURPOSE, OR FOR DAMAGES AND/OR INJURIES INCURRED BY BUYER, ITS AGENTS ~ 1/4 AND EMPLOYEES, BY VIRTUE OF DEFECTS IN THE PRODUCTS, EVEN IF BIODEx HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES. Use or misuse by the buyer and/or agents, employees or others will, at the option of Biodef, void all product warranties.

This catalog information herein is the property of Biodef Medical Systems. Copying in whole or in part, transmittal to other individuals or companies, and use for any purpose other than that for which it is intended, are expressly prohibited.

- Phantom Cabinet39
 Preparation Enclosure Base Cabinet35
 Preparation Enclosure35
 Radioisotope Storage Cabinet34
 Refrigerator39
 Sink and Waste Cabinet37
 Unit Dose Cabinet38
 Unit Dose Cabinet, PET40
 Unit Dose Table, PET41
 Waste and Storage Cabinet37
 Waste Cabinet36
 Lineator74
 Liver Marker/Ruler93
 Lung Ventilation Systems and Accessories80-92
- M**attress, Table Pad119
- Meters:**
 Ratemeter, Alarm106
 Survey104-105, 109-112
 Mobile Barriers, Clear-Lead29
 Mobile Radiation Shields28
 Moly Assay Shields74
- Monitors:**
 Digital Area107
 Entrance/Exit Area108
 Xenon Trap85
 Mouthpieces89
- N**iptongs30, 116
 Nose Clips89
- P**aper, Absorbent116
 Patient Transfer Board, Easy Mover118
 Personal Radiation Protection31-32
- PET Accessories:**
 Atomlab Dose Calibrators64-73
 Dose Calibrator Shielding Rings75
 Dose Drawing Syringe Shields57
 Dose Drawing System56-57
 Dose Drawing Vial Shields57
 Dose Preparation Stand56-57
 L-Block Shields19-21
 Manual Dose Injector13
 Mobile Radiation Shields28
 Pig Rack47
 Sharps Container Shield18
 Shipping Systems50-55
 Syringe Shields8-9
 Unit Dose Cabinet40
 Unit Dose Pigs47, 50-52
 Unit Dose Table, PET41
 Vial Pig43, 54
 Vial Shield53
- Phantoms96-103**
 Anthropomorphic Torso98
 Bar102
 Flood, Rectangular103
 Hoffman 3-D Brain100
 Lung-Spine98
 PET96-97, 100-101
 SPECT, Jaszczak96-97
 Thyroid, Neck103
- Phantom Cabinet39
 Phantom Inserts96, 99
- Pigs and Accessories:**
 BEXXAR Unit Dose Pig48
 PET Double-Ended Pig55
 PET Shipping Systems50-55
 PET Unit Dose Pigs50-52
 PET Vial Pig43, 54
 Pro-Tec Pig Rack47
 Pro-Tec Pig Shipping Bag47
 Pro-Tec Unit Dose Pigs46
 Unit Dose Pig Wall Rack11
- Zevalin Unit Dose Pig45
 Positioning Wedges120
 Preparation Enclosure35
 Preparation Enclosure Base Cabinet35
 Pulmonex Hose Kit91
 Pulmonex II Xenon System84-85
- R**adiacwash114
 Radiation Attenuating Gloves31
 Radiation Detection104-113
 Radiation Shields, Mobile28
 Radioaerosol Delivery Systems80-82
 Radioiodine Fume Hood26
 Radioisotope Storage Cabinet34
 Radiopharmaceutical QC Manual60
 Radiopharmacy42-63
 Ratemeter, Alarm106
 Refrigerator, Lead-Lined39
- S**afe, Storage, Lead-Lined26
 Sharps Shields and Containers16-18
 Shield, Bench Top23
 Shield, Injection Pump, Graseby13
 Shielded Decay Drums59
 Shielded Isolator62-63
 Shielded Rectangular Container30
 Shielded Waste Containers27
 Shielding and Storage19-30
 Shields, L-Block19-21
 Shields, Moly Assay74
 Shields, Tabletop22
 Shipping Bag, Pig, Pro-Tec47
 Shipping Container, Intego53
 Shipping Container, Zevalin45
 Shipping System, BEXXAR48
 Shipping Systems, PET50-55
 Shipping System, Vial Pig, I-13149
 Signs, Labels and Tags121-123
 Sink and Waste Cabinet37
 Soda-Lime and Drierite90
- Sources93-95**
 Check93
 Dose Calibrator Reference95
 Fillable Point Marker93
 Flood, Rectangular94
 Liver Marker/Ruler93
 Pen Point Marker94
 Radioactive Ruler94
 Rod93
 Spot Marker93
 SPECT Phantoms96-97
 Stand, Injection14
 Stools118
 Storage Cabinet, Shielded30
 Storage Containers, Shielded16
 Storage Safe26
- Survey Meters:**
 Cutie Pie Digital112
 Inspector110
 Model 14C104-105
 Monitor 4109
 Surface111
 Syringe and Vial Shields4-18
 Syringe Carriers, Shielded12
 Syringe Dippers74-75
 Syringe Holder, Shielded11
 Syringe Labels, Color-Coded44
 Syringe Recapper11
- Syringe Shields:**
 Beta10
 Color-Coded42
 Dose Drawing47
 Dose Drawing, PET57
 Dose Preparation, BEXXAR58
- High Density Lead Glass6
 Myoview7
 PET8-9
 Pro-Tec II4
 Pro-Tec III5
 Pro-Tec IV6
 Quadramet7
 Zevalin10
- T**able Pad119
 Table, Steel25
 Table, Unit Dose, PET41
 Tabletop Shields22
 Tags, Pre-Strung121
 Tape, Hazard and Warning121-122
 Tec-Control Chromatography Systems60-61
 Thyroid Phantom103
 Thyroid Uptake System76-79
 Tubing Splice90
 Tubing, Corrugated90
- U**nit Dose Cabinet38
 Unit Dose Pig, BEXXAR47
 Unit Dose Pig, PET47
 Unit Dose Pig, Pro-Tec46
 Unit Dose Pig, Zevalin45
 Unit Dose Pig Wall Rack11
 Unit Dose Table, PET41
- V**einlite LED15
 Venti-Scan IV80-81
 Vial Pig, PET43
- Vial Shields:**
 .25" Lead43
 Dose Drawing, PET57
 Dose Preparation, BEXXAR58
 High Density Lead Glass42
 PET53
 Tungsten43
 Zevalin43
- Vinyl Sheets, Lead23
- W**all Rack, Unit Dose Pig11
 Warning Tape & Rope121-122
 Waste Containers, Shielded27
 Wedges, Positioning120
 Well Counters77
 Well Insert75
 Wipe Test Counter, Atomlab70-71
 Wipe Test Kits114
 Workbench, Cyclotron57
- X**enon Convenience Kits86-87
 Xenon Dispenser, Automatic85
 Xenon Rebreathing System92
 Xenon System, Pulmonex II84-85
 Xenon Trap Monitor85
 Xenovap Administration System91
- "Y"** Connectors90
- Z**evalin Accessories:
 Manual Dose Injector13
 Shipping Container45
 Syringe Shields10
 Unit Dose Pig45
 Vial Shield43

"The Clinical Advantage"™

ATOMLAB

Dose Calibrators & Wipe Test Counters



Proven performance for fast, accurate measurements.

See pages 64-75 for details.

Visit www.biomedex.com/dosecalibrators

BIODEX

www.biomedex.com

1-800-224-6339

Int'l 631-924-9000

BIODEX

Biodex Medical Systems, Inc.
20 Ramsey Road, Shirley, New York, 11967-4704

CHANGE SERVICE REQUESTED

**SUBSCRIBE to
BIODEX**

- NEWSLETTERS
- BLOGS
- EMAIL

www.biodex.com/eregister



Follow Biodex-Nuclear Medicine

Let's get a discussion going...

FN: 12-187 7/12

"The Clinical Advantage"™

SERVING THE
NUCLEAR MEDICINE COMMUNITY
FOR MORE THAN 60 YEARS

www.biodex.com/nuclearmedicine

BIODEX

www.biodex.com

1-800-224-6339

Int'l 631-924-9000