



Lung kV, 2 VISICOIL™ 0.75 mm x 1 cm



Prostate kV, 2 VISICOIL™ 0.5 mm x 1 cm








Pancreas kV, 3 VISICOIL™ 0.35 mm x 1 cm




**The new 0.5 mm VISICOIL™ in pre-loaded 19 Gauge Needle:**


- ✓ Reduced patient discomfort during implantation
- ✓ Designed to minimize edema
- ✓ Increased implantation safety especially for lung cases
- ✓ Recommended for: KV, CBCT, OBI, Ultrasound, Proton Therapy



Breast



Lung



Prostate

**Safety:** Ribbed flexible VISICOIL™ designed NOT to migrate

**Efficiency:** Linear marker, less implantation (1 VISICOIL™ = 2 reference points)

**Visibility:** Hollow core reduces imaging artifacts

**Flexibility:** Multiple sizes to match your imaging, treatment indication and workflow